



An ISO 9001 : 2015 Certified Company



LABORATORY GLASSWARE

Made from Borosilicate 3.3 Glass

PRICE LIST
1st April, 2022

Manufactured by :

J-SIL SCIENTIFIC INDUSTRIES

46, 48A & 48B, INDUSTRIAL ESTATE, NUNHAI, AGRA-282006 (INDIA)

Phone : 0562-2281967, 2281887, 2280186

jsilindia@gmail.com; sales.agra@j-sil.in

www.j-sil.com, www.j-sil.in

Authorised Stockist

J-SIL TRADING CORPORATION

21/29, Popat Wadi, Kalbadevi Road,

Mumbai - 400 002, INDIA

Phone : 022 - 2205 7568 / 2200 4293

E-mail : sales@j-sil.in



J-SIL
www.j-sil.in

Certificate



E-Certificate are available
on our website www.j-sil.in

Certificate for Volumetric Flask

J-SIL SCIENTIFIC INDUSTRIES
48-B INDUSTRIAL ESTATE, NUNHAI, AGRA-282006 (INDIA)
E-mail : sales@j-sil.in URL : www.j-sil.in

CLASS "A" TEST CERTIFICATE OF ACCURACY

Item Name : • Volumetric Flask
Test Method : • Gravimetric Method
Item Serial No. : • 183/ND
Normal Volume : • 500 ml.
Tolerance± : • 0.15 ml.
Certified Volume : • 500.03 ml. at 27°C
Delivery Time : • N/A
Environment Conditions : • Temperature : • 24.5°C Relative Humidity : • 50%

The product has been calibrated with following devices. All the below devices have valid Calibration traceable to NPL, India, Traceable to International Standards.

Instruments	Make	Serial No.	Validity Dt.	Range	Least Count
Balance	Mettler Toledo	4749506904	05/03/2021	1600 g	0.01 g
Thermometer	J.R.M.	R-18122	05/03/2021	(0-50) °C	0.1 °C

The above product has been tested with distilled water according to the IP / INDIAN STANDARD Standard Specification & found to the satisfied Tolerance of CLASS 'A'

The above calibration results is valid at the time of calibration under lab. condition. This Report shall not be reproduced , except in full, unless with written permission from M/s. J-SIL SCIENTIFIC INDUSTRIES.

Tested By : ROHIT
Date : 27/01/2021


ISO 9001 :2015

Jaipal Agarwal
Technical Officer
E-Certificate

Certificate for Measuring Cylinder

J-SIL SCIENTIFIC INDUSTRIES
48-B INDUSTRIAL ESTATE, NUNHAI, AGRA-282006 (INDIA)
E-mail : sales@j-sil.in URL : www.j-sil.in

CLASS "A" TEST CERTIFICATE OF ACCURACY

Item Name : • Measuring Cylinder (Stoppered) Hexagonal Base
Test Method : • Gravimetric Method
Item Serial No. : • 280/ND
Normal Volume : • 500 ml.
Tolerance± : • 2.50 ml.
Certified Volume : • 499.65 ml. at 27°C
Delivery Time : • N/A
Environment Conditions : • Temperature : • 24.5°C Relative Humidity : • 52%

The product has been calibrated with following devices. All the below devices have valid Calibration traceable to NPL, India, Traceable to International Standards.

Instruments	Make	Serial No.	Validity Dt.	Range	Least Count
Balance	Mettler Toledo	4749506904	05/03/2021	1600 g	0.01 g
Thermometer	J.R.M.	T-131021	20/09/2021	(0-50) °C	0.1 °C

The above product has been tested with distilled water according to the IP / INDIAN STANDARD Standard Specification & found to the satisfied Tolerance of CLASS 'A'

The above calibration results is valid at the time of calibration under lab. condition. This Report shall not be reproduced , except in full, unless with written permission from M/s. J-SIL SCIENTIFIC INDUSTRIES.

Tested By : TRILOK
Date : 09/02/2021


ISO 9001 :2015

Jaipal Agarwal
Technical Officer
E-Certificate

Certificate for Graduated Pipette

J-SIL SCIENTIFIC INDUSTRIES
48-B INDUSTRIAL ESTATE, NUNHAI, AGRA-282006 (INDIA)
E-mail : jsil@sancharnet.in URL : www.j-sil.com

CLASS "A" TEST CERTIFICATE OF ACCURACY

Item Name : • Graduated Pipette (Serological)
Test Method : • Gravimetric Method
Item Serial No. : • 375/CD
Normal Volume : • 10 ml.
Tolerance± : • 0.05 ml.
Certified Volume : • 9.991 ml. at 27.00°C
Delivery Time : • 18 sec
Environment Conditions : • Temperature 24.5 °C, Relative Humidity, 53 %.

The product has been calibrated with following devices. All the below devices have valid Calibration traceable to NPL, India, Traceable to International Standards.

Instruments	Make	Serial No.	Validity Dt.	Range	Least Count
Balance	Mettler Toledo	1226200463	05/03/2021	210 g	0.001 g
Thermometer	J.R.M.	M-13222	20/09/2021	(0-50) °C	0.1 °C

The above product has been tested with distilled water according to the IP / INDIAN STANDARD Standard Specification & found to the satisfied Tolerance of CLASS 'A'

The above calibration results is valid at the time of calibration under lab. condition. This Report shall not be reproduced , except in full, unless with written permission from M/s. J-SIL SCIENTIFIC INDUSTRIES.

Tested By : DASHRATH
Date : 01/02/2021


ISO 9001 :2008

Jaipal Agarwal
Technical Officer
E-Certificate

Certificate for Burette

J-SIL SCIENTIFIC INDUSTRIES
48-B INDUSTRIAL ESTATE, NUNHAI, AGRA-282006 (INDIA)
E-mail : jsil@sancharnet.in URL : www.j-sil.com

CLASS "A" TEST CERTIFICATE OF ACCURACY

Item Name : • Burette With PTFE Key Stop Cock (USP)
Test Method : • Gravimetric Method
Item Serial No. : • 843/CG
Normal Volume : • 50.00 ml.
Subdivision : • 0.1 ml.
Tolerance± : • 0.05 ml.

	From	To	Certified Volume	Temperature
Upto the gradation	0.00	• 10.00 ml.	• 9.99 ml.	• 20.00 at °C
Upto the gradation	0.00	• 25.00 ml.	• 24.97 ml.	• 20.00 at °C
Upto the gradation	0.00	• 40.00 ml.	• 39.96 ml.	• 20.00 at °C
Upto the gradation	0.00	• 50.00 ml.	• 49.96 ml.	• 20.00 at °C

Delivery Time : • 76 Sec.
Environment Conditions : • Temperature 25.2 °C, Relative Humidity, 50 %.

The product has been calibrated with following devices. All the below devices have valid Calibration traceable to NPL, India, Traceable to International Standards.

Instruments	Make	Serial No.	Validity Dt.	Range	Least Count
Balance	Mettler Toledo	1226200463	05/03/2021	210 g	0.001 g
Thermometer	J.R.M.	M-13222	20/09/2021	(0-50) °C	0.1 °C

The above product has been tested with distilled water according to the USP Standard Specification & found to the satisfied Tolerance of CLASS 'A'

The above calibration results is valid at the time of calibration under lab. condition. This Report shall not be reproduced , except in full, unless with written permission from M/s. J-SIL SCIENTIFIC INDUSTRIES.

Tested By : DASHRATH
Date : 05/02/2021


ISO 9001 :2008

Jaipal Agarwal
Technical Officer
E-Certificate






Vials & Septa



51 Chromatography Graduated, Screw Vials 9 mm wide Opening with **White Patch**

- 9mm screw thread vials are designed especially for use with any rotating or robotic are samplers 2ml, 12 x 32mm vials are manufactured of clear & amber
- Made from USP glass 1 and include a write - in patch for sample identification

	Type	Cap	Pkt	Price/ Pkt
A	Clear Glass	2 ml	100 Pc	700.00
B	Amber Glass	2 ml	100 Pc	800.00



56 Screw Cap with Septa for 9 mm, Screw Vials (Closure)

- Caps are made of High quality polypropylene to exact manufacturing tolerances and lined in a controlled manufacturing environment
- Septa use only the highest quality material to ensure proper function and can be Pre-slit to ease needle Penetration
- Absence of elastomer layer reduces Potential for Contamination

	Cap	Septa	Pkt	Price/ Pkt
A	Blue	Non Slit/Non Bonded	100 Pc	700.00
B	Blue	Pre-Slit/Non Bonded	100 Pc	800.00
C	Blue, Pre-Slit/ Bonded Red PTFE/White Silicon		100 Pc	1150.00
D	Blue, Pre-Slit/ Bonded White PTFE/Beige Silicon		100 Pc	1200.00



61 Clear Glass Screw Vials 9mm with Screw Blue Cap & PTFE-Silicon Septa (**Combo**)

	Cap	Septa	Pkt	Price/ Pkt
A	2 ml	Non Slit/Non Bonded	100 Pc	1400.00
B	2 ml	Pre-Slit/Non Bonded	100 Pc	1500.00
C	2 ml, Pre-Slit/ Bonded Red PTFE/White Silicon		100 Pc	1850.00
D	2 ml, Pre-Slit/ Bonded White PTFE/Beige Silicon		100 Pc	1900.00



66 Amber Clear Glass Screw Vials 9mm with Screw Blue Cap & PTFE-Silicon Septa (**Combo**)

	Cap	Septa	Pkt	Price/ Pkt
A	2 ml	Non Slit/Non Bonded	100 Pc	1500.00
B	2 ml	Pre-Slit/Non Bonded	100 Pc	1600.00
C	2 ml, Pre-Slit/ Bonded Red PTFE/White Silicon		100 Pc	1950.00
D	2 ml, Pre-Slit/ Bonded White PTFE/Beige Silicon		100 Pc	2000.00



J-SIL
www.j-sil.in



71 Nylon Syringe Filter (Non Sterile) Packed in Transparent Plastic Jar

High quality and cost effective, Nylon syringe filters are designed for sterilization or clarification of aqueous and solvent samples in a wide array of analytical and research applications. It is primarily used to filter small aqueous and organic samples with low hold-up volume in order to inject them directly into the HPLC sample loop system.

Disposable

- Luer-Lok and pressure fitting
- Ultrasonically bonded for reliability
- Clearly inscribed with the filter media type on barrel
- Excellent unit-to-unit and lot-to-lot reproducibility and consistency
- Made of high purity virgin polypropylene housing, no dyes or pigments added
- Unique with a concentric funnel shaped barrel to increase volume and decrease back pressure



	DIA	MICRON	Pkt	Price/ Pkt
A	13 mm	0.22 µm	100 Pc	2000.00
B	13 mm	0.45 µm	100 Pc	2000.00
C	25 mm	0.22 µm	100 Pc	2150.00
D	25 mm	0.45 µm	100 Pc	2050.00
E	30 mm	0.22 µm	50 Pc	1500.00
F	30 mm	0.45 µm	50 Pc	1500.00



Dessicators



Refer Cat. No. 1906





Syringe Filter

76 PVDF Syringe Filter (Non Sterile) Hydrophilic, Packed in Transparent Plastic Jar

Syringe filters are designed for aqueous and organic solution samples. The sample prep syringe filters are used for clarification and particle removal from samples before instrumental analysis such as HPLC, IC, and GC. These syringe filters are available in a variety of filter membrane materials, offered in 0.2 μm or 0.45 μm pore sizes.

Disposable

- Luer-Lok and pressure fitting
- Ultrasonically bonded for reliability
- Clearly inscribed with the filter media type on barrel
- Excellent unit-to-unit and lot-to-lot reproducibility and consistency
- Made of high purity virgin polypropylene housing, no dyes or pigments added
- Unique with a concentric funnel shaped barrel to increase volume and decrease back pressure



	DIA	MICRON	Pkt	Price/ Pkt
A	13 mm	0.22 μm	100 Pc	3000.00
B	13 mm	0.45 μm	100 Pc	3000.00
C	25 mm	0.22 μm	100 Pc	3300.00
D	25 mm	0.45 μm	100 Pc	3200.00
E	30 mm	0.22 μm	50 Pc	3100.00
F	30 mm	0.45 μm	50 Pc	3100.00



Dessicators



Refer Cat. No. 1911





81 PTFE Syringe Filter (Non Sterile) Packed in Transparent Plastic Jar

High quality and cost effective, these Polytetrafluoroethylene (PTFE) syringe filters are HPLC certified and available in hydrophobic version.

- PTFE membranes are inherently hydrophobic; they are ideal for the filtration of nonpolar liquids and gases.
- Very resistant to solvents as well as acids or bases.
- Can be made more hydrophilic by flushing with alcohol, followed by water
- Used for filtering and de-gassing chromatography solvents
- Permits optimal utilization of small sample volumes
- Excellent chemical compatibility for filtering harsh chemicals that destroy other membrane materials
- Moisture barrier for venting applications and vacuum pump line protection
- Ideal for aerosol sampling
- Filter viscous organic based HPLC samples
- Glass microfiber prefilter raises absorption capacity and loading
- After pre-wetting with alcohol, may be used for fine particle removal of aqueous and mild organic solutions in Ion Chromatography



	DIA	MICRON	Pkt	Price/ Pkt
A	13 mm	0.22 μ m	100 Pc	3100.00
B	13 mm	0.45 μ m	100 Pc	3100.00
C	25 mm	0.22 μ m	100 Pc	3400.00
D	25 mm	0.45 μ m	100 Pc	3300.00
E	30 mm	0.22 μ m	50 Pc	3000.00
F	30 mm	0.45 μ m	50 Pc	3000.00



S.S. Tong



Refer Cat. No. 1370





Micropipettes




JFE Micropipette, **Fixed** Volume Fully Autoclavable, As Per NABL Certified

- The single channel Fixed Volume Micropipette is recommended **for institutes & university**
- A calibration certificate is included inside the product package.
- A calibration tool is also included for quick in-lab recalibration.
- The instrument is fully autoclavable at 121 °C and 15 psi for a duration of 10-15 minutes.

Cat. No.	Cap. μ l	Accuracy $\pm\mu$ l	CV $\pm\mu$ l	Price each pc
1	1 μ l	0.050	0.050	1500.00
2	2 μ l	0.080	0.040	1500.00
5	5 μ l	0.100	0.050	1500.00
10	10 μ l	0.100	0.050	1500.00
20	20 μ l	0.160	0.080	1500.00
25	25 μ l	0.200	0.100	1500.00
50	50 μ l	0.400	0.200	1500.00
100	100 μ l	0.600	0.200	1500.00
200	200 μ l	1.200	0.400	1500.00
250	250 μ l	1.500	0.500	1500.00
500	500 μ l	3.000	1.000	1500.00
1000	1000 μ l	6.000	2.000	1500.00
2000	2000 μ l	12.000	4.000	1500.00
5000	5000 μ l	30.000	10.000	1500.00
10000	10000 μ l	60.000	20.000	1500.00

 Recommended for Schools, Colleges & Education Labs

Forcep



Refer Cat. No. 1497 & 1499



Micropipettes



JVE Micropipette, **Variable** Volume Fully Autoclavable, As Per NABL Certified

- The single channel Variable Volume Micropipette is recommended **for institutes & university**
- A calibration certificate is included inside the product package.
- A calibration tool is also included for quick in-lab recalibration.
- The instrument is fully autoclavable at 121 °C and 15 psi for a duration of 10-15 minutes.

Cat. No.	Volume range µl	Increment µl	Accuracy ±µl	CV ±µl	Price each pc
2	0.2-2µl	0.01	0.040	0.024	2800.00
10	0.5-10µl	0.1	0.100	0.050	2800.00
20	2-20µl	0.1	0.160	0.080	2800.00
50	5-50µl	1	0.400	0.200	2800.00
100	10-100µl	1	0.600	0.200	2800.00
200	20-200µl	1	1.200	0.400	2800.00
1000	100-1000µl	10	6.000	2.000	2800.00
5000	500-5000µl	100	30.000	10.000	2800.00
10000	1000-10000µl	100	60.000	20.000	2800.00

👉 Recommended for Schools, Colleges & Education Labs

Spatulla



Refer Cat. No. 1492, 1494 & 1495





Micropipettes




JFR Micropipette, **Fixed** Volume Fully Auto-clavable, As Per ISO 8655 (**Research series**) NABL Certified

- A range of High Precision Micropipettes that offers a sophisticated blend of features and functionality. Built by our in-house team of product design engineers, it is categorised into Single Channel (Fixed and Variable Volume) Micropipettes. These instruments are designed with ergonomics and intuitive handling in mind facilitate remarkable user experiences in practical laboratory environments.
- The Single Channel Fixed Volume Micropipette is recommended for Clinical Diagnostics, Control Analysis, etc.
- A calibration certificate is included inside the product package.
- A calibration tool is also included for quick in-lab recalibration.
- The instrument is fully autoclavable at 121 °C and 15 psi for a duration of 10-15 minutes.

Cat. No.	Volume range µl	Accuracy ±µl	CV ±µl	Price each pc
1	1µl	0.050	0.050	2950.00
2	2µl	0.080	0.040	2950.00
5	5µl	0.100	0.050	2950.00
10	10µl	0.100	0.050	2950.00
20	20µl	0.160	0.080	2950.00
25	25µl	0.200	0.100	2950.00
50	50µl	0.400	0.200	2950.00
100	100µl	0.600	0.200	2950.00
200	200µl	1.200	0.400	2950.00
250	250µl	1.500	0.500	2950.00
500	500µl	3.000	1.000	2950.00
1000	1000µl	6.000	2.000	2950.00
2000	2000µl	12.000	4.000	2950.00
5000	5000µl	30.000	10.000	2950.00
10000	10000µl	60.000	20.000	2950.00

 Recommended for QC & RND Labs

Boss Head



Price on Request



J-SIL
www.j-sil.in

Micropipettes




JVR Micropipette, **Variable** Volume Fully Auto-clavable, As per ISO 8655 (**Research series**) NABL certified

- A, range of High Precision Micropipettes that offers a sophisticated blend of features and functionality
- it is categorised into Single Channel Variable Volume
- The Single Channel Variable Volume Micropipette is prominently used for Molecular Biology, Microbiology, Immunology, Cell Culture, Analytical Chemistry, Biochemistry, Genetics, etc.
- The instrument includes a Holder that enables easy, efficient and safe storage.
- The instrument is fully autoclavable at 121 °C and 15 psi for a duration of 10-15 minutes.

Cat. No.	Volume range µl	Increment µl	Accuracy ±µl	CV ±µl	Price each pc
2	0.2-2 µl	0.002	0.040	0.024	5500.00
10	0.5-10µl	0.02	0.100	0.050	5500.00
20	2-20µl	0.02	0.160	0.080	5500.00
50	5-50µl	0.1	0.400	0.200	5500.00
100	10-100µl	0.2	0.600	0.200	5500.00
200	20-200µl	0.2	1.200	0.400	5500.00
1000	100-1000µl	1	6.000	2.000	5500.00
5000	500-5000µl	10	30.000	10.000	5500.00
10000	1000-10000µl	20	60.000	20.000	5500.00

 Recommended for QC & RND Labs

3 Finger Clamp



Price on Request





Micropipette & Safety Rings



JNP Nano, Fixed Volume Micropipette, As per ISO 8655 NABL Certified

- Nano miniature fixed volume micropipette, is a low cost, high precision liquid handling instrument. It offers Two-step Plunger Operation and is specially for use with diagnostic kits.
- Fully autoclavable at 121 °C and 15 psi for a duration of 10 – 15 minutes.
- Comes with one tip


Cat. No.	Volume range μ l	Accuracy $\pm\mu$ l	CV $\pm\mu$ l	Price each pc
5	5 μ l	0.100	0.050	450.00
10	10 μ l	0.100	0.050	450.00
20	20 μ l	0.160	0.080	450.00
25	25 μ l	0.200	0.100	450.00
50	50 μ l	0.400	0.200	450.00
100	100 μ l	0.600	0.200	450.00
200	200 μ l	1.200	0.400	450.00
250	250 μ l	1.500	0.500	450.00
500	500 μ l	3.000	1.000	450.00
1000	1000 μ l	6.000	2.000	450.00

 Recommended for Diagnostic Kits


RVF Safety Rings For Volumetric Flask

- It is fully autoclavable at 121°C and 15 psi for a duration of 10-15 minutes.
- Use high quality silicon material for long lasting

Cat. No.	Cap. ml	Joint size	Minimum Quantity	Price each pc
01	5ml	10/15	12 Pc	80.00
02	10ml	10/15	12 Pc	85.00
03	20ml	10/15	12 Pc	85.00
04	25ml	10/15	12 Pc	85.00
05	50ml	10/15	12 Pc	85.00
06	100ml	14/23	12 Pc	110.00
07	200ml	14/23	12 Pc	110.00
08	250ml	14/23	12 Pc	110.00
09	500ml	19/26	12 Pc	140.00
10	1000ml	19/26	12 Pc	140.00



Burettes



101 Burettes, with Straight Bore Glass Stop-cock, Accuracy as per Class 'B' of IS 1997 :1982


- Generally used in Titration process
- In clear glass & Straight bore stop cock with screw cap
- Printed in permanent amber colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	1	0.01	0.01	2 Pc	700.00
B	2	0.02	0.02	2 Pc	700.00
C	5	0.05	0.02	2 Pc	700.00
F	10	0.05	0.05	2 Pc	600.00
H	25	0.10	0.10	2 Pc	500.00
I	25	0.05	0.06	2 Pc	1100.00
J	50	0.10	0.10	2 Pc	525.00
K	100	0.20	0.20	1 Pc	900.00
L	25	White Graduation		2 Pc	560.00
M	50	White Graduation		2 Pc	585.00

106 Burettes, Single Bore Teflon Key (Rota Flow), Accuracy as per Class 'B' of IS 1997 :1982

- Generally used in Titration process
- Comes with Rotaflow PTFE key, lubricate free stop cock
- Printed in permanent amber colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	1	0.01	0.01	2 Pc	650.00
B	2	0.02	0.02	2 Pc	650.00
C	5	0.05	0.02	2 Pc	650.00
F	10	0.05	0.05	2 Pc	550.00
G	25	0.10	0.10	2 Pc	450.00
H	25	0.05	0.06	2 Pc	1100.00
I	50	0.10	0.10	2 Pc	460.00
J	100	0.20	0.20	1 Pc	860.00
K	25	White Graduation		2 Pc	500.00
L	50	White Graduation		2 Pc	530.00





Burettes



107 Burettes **Amber Colour** Single Bore Teflon Key (**Rota Flow**) Stop Cock Accuracy as per **Class 'B'**

- Generally used in Titration process
- Comes with Rotaflow PTFE key, lubricate free stop cock
- Printed in white colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
D	10	0.05	0.05	1 Pc	2500.00
E	25	0.10	0.10	1 Pc	2350.00
F	50	0.10	0.10	1 Pc	2400.00
G	100	0.20	0.20	1 Pc	2600.00

111 Burettes Straight Bore with **Teflon Key** Stop Cock **Class 'B'**

- Generally used in Titration process
- Comes with straight bore PTFE key, lubricate free stop cock
- Printed in permanent Amber colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	10	0.05	0.05	1 Pc	750.00
B	25	0.10	0.10	1 Pc	710.00
C	50	0.10	0.10	1 Pc	800.00
D	100	0.20	0.20	1 Pc	1000.00



112 Burette clamp, Mild steel (**Single**)

- Mild steel, Powder Coated, Fisher Type with Boss Head

Cat No.	Type	Max. Temp.	Pack per box	Price each pc
A	Silicon Sleeve	280°C	12 pc	190.00

113 Burette clamp, Mild steel (**Double**)

- Mild steel, Powder Coated, Fisher Type with Boss Head

Cat No.	Type	Max. Temp.	Pack per box	Price each pc
A	Silicon Sleeve	280°C	6 pc	335.00



J-SIL
www.j-sil.in



116 Burettes Straight Bore with **Teflon Key Stop-cock Class 'A'** with work's Certificate & Calibrated @ 27°C.

- Generally used in Titration process
- Comes with straight bore PTFE key, lubricate free stop cock
- Comes with Individual Work Certificate in clear Glass
- Printed in permanent Amber colour


	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	10	0.05	0.02	1 Pc	1500.00
B	25	0.10	0.05	1 Pc	1500.00
C	50	0.10	0.05	1 Pc	1550.00
D	100	0.20	0.10	1 Pc	1700.00

117 Burettes Straight Bore with **Teflon Key Stop-cock Accuracy as per USP (US Pharmacopoeia) Class A** with work's Certificate & Calibrated @ 20°C

- Generally used in Titration process
- Comes with straight bore PTFE key, lubricate free stop cock
- Comes with Individual Work Certificate in clear Glass
- Printed in permanent Amber colour


	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	10	0.05	0.02	1 Pc	1600.00
B	25	0.10	0.03	1 Pc	1600.00
C	50	0.10	0.05	1 Pc	1650.00
D	100	0.20	0.10	1 Pc	1800.00

ASTM



Refer Cat. No. 112





121 Burettes, with Straight Bore Glass Stop-cock, Accuracy as per of I.S. 1997 : 1982 Class A with work's Certificate & Calibrated @ 27°C

- Generally used in Titration process
- In clear Glass & Straight bore stop cock with screw cap
- Comes with Individual Work Certificate
- Printed in permanent Amber colour



	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	1	0.01	0.006	1 Pc	1500.00
B	2	0.02	0.01	1 Pc	1500.00
C	5	0.05	0.02	1 Pc	1500.00
F	10	0.05	0.02	1 Pc	1400.00
G	25	0.10	0.05	1 Pc	1400.00
H	25	0.05	0.03	1 Pc	2100.00
I	50	0.10	0.05	1 Pc	1550.00
J	100	0.20	0.10	1 Pc	1600.00

125 Burettes Single Bore with Teflon Key (Rota Flow) Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C

- Generally used in Titration process
- Comes with straight bore PTFE key, lubricate free stop cock
- Comes with Individual Work Certificate in clear Glass
- Printed in permanent Amber colour

	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	10	0.05	0.02	1 Pc	1500.00
B	25	0.10	0.03	1 Pc	1500.00
C	50	0.10	0.05	1 Pc	1550.00
D	100	0.20	0.10	1 Pc	1850.00

ASTM




Refer Cat. No. 113



J-SIL
www.j-sil.in

Burettes



126 Burettes, Single Bore Teflon Key (Rota Flow) Class A with work's Certificate & Calibrated @ 27°C.

- Generally used in Titration process
- Comes with Rotaflow PTFE key, lubricate free stop cock
- Comes with Individual Work Certificate in clear Glass
- Printed in permanent Amber colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	1	0.01	0.006	1 Pc	1500.00
B	2	0.02	0.01	1 Pc	1500.00
C	5	0.05	0.02	1 Pc	1500.00
F	10	0.05	0.02	1 Pc	1400.00
H	25	0.10	0.05	1 Pc	1400.00
I	25	0.05	0.03	1 Pc	2000.00
J	50	0.10	0.05	1 Pc	1450.00
K	100	0.20	0.10	1 Pc	1950.00

127 Burettes Amber Colour Single Bore Teflon Key (Rota Flow) Accuracy as per Class A with work's Certificate & Calibrated @ 27°C

- Generally used in Titration process
- Comes with Rotaflow PTFE key, lubricate free stop cock
- Comes with Individual Work Certificate in Amber Glass
- Printed in white colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	1	0.01	0.006	1 Pc	3000.00
B	2	0.02	0.01	1 Pc	3000.00
C	5	0.05	0.01	1 Pc	3000.00
D	10	0.05	0.02	1 Pc	2500.00
E	25	0.10	0.05	1 Pc	2400.00
F	50	0.10	0.05	1 Pc	2500.00
G	100	0.20	0.10	1 Pc	3000.00





Burettes



ASTM



128 Burettes Amber Colour Single Bore with Teflon Key (Rota Flow) Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C


- Generally used in Titration process
- Comes with Rotaflow PTFE key, lubricate free stop cock
- Comes with Individual Work Certificate in Amber Glass
- Printed in white colour

	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	10	0.05	0.02	1 Pc	3000.00
B	25	0.10	0.03	1 Pc	2600.00
C	50	0.10	0.05	1 Pc	2800.00
D	100	0.20	0.10	1 Pc	3500.00

131 Burettes, for Pinch Cock with Jet, Rubber Tube, Accuracy as per Class 'B' of I.S.1997 :1982

- Generally used in Titration process
- Comes with Pinch Clip, rubber tube & glass jet
- Printed in permanent Amber colour


	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	10	0.05	0.05	6 Pc	200.00
B	25	0.10	0.10	5 Pc	230.00
C	50	0.10	0.10	5 Pc	260.00
D	100	0.20	0.20	4 Pc	500.00
E	Spare Glass Jet with Rubber Tube				70.00
F	Spare Pinch Clip for Burettes				30.00




132 Fisher clamp

- These clamps are moulded in polypropylene.
- The spring loaded jack in the middle provides a vice like grip to the burette and prevents it from slipping.

Cat No.	Cap. ml.	Pack per box	Price per pc
A	Fisher clamp- single	12 pcs.	144.00
B	Fisher clamp- double	6 pcs	288.00





176 Burettes, Automatic Zero Mounted on Borosilicate glass Reservoir with Glass Stop Cock Accuracy as per Class 'B' of I.S. 1997 : 1982 complete with rubber pump


- With auto zero adjustment
- In clear Glass & Straight bore Glass stop cock with screw cap
- Printed in permanent Amber colour
- Comes with Rubber Bellow

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Cap of Reservoir	Price each pc
A	10	0.05	0.05	500 ml	5200.00
B	25	0.10	0.10	2000 ml	5500.00
C	50	0.10	0.10	2000 ml	5950.00
D	100	0.20	0.20	2000 ml	6700.00
Spare Rubber Bellow					700.00

181 Burettes, Automatic Zero Mounted on Borosilicate glass Reservoir with Teflon Stop cock (Rota Flow). Accuracy as per Class 'B'

- With auto zero adjustment
- In clear Glass with Rotaflow PTFE key, lubricate free stop cock
- Printed in permanent Amber colour
- Comes with Rubber Bellow

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Cap of Reservoir	Price each pc
A	10	0.05	0.05	500 ml	5200.00
B	25	0.10	0.10	2000 ml	5500.00
C	50	0.10	0.10	2000 ml	5950.00
D	100	0.20	0.20	2000 ml	6500.00
Spare Rubber Bellow					700.00





182 Burettes Amber Colour Automatic Zero Mounted on Reservoir Accuracy as per Class 'B' with Teflon Stop cock (Rota Flow) complete with rubber pump


- With auto zero adjustment
- In Amber Glass with Rotaflow PTFE key, lubricate free stop cock
- Printed in permanent Amber colour
- Comes with Rubber Bellow


	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Cap of Reservoir	Price each pc
A	10	0.05	0.05	500 ml	17500.00
B	25	0.10	0.10	2000 ml	18000.00
C	50	0.10	0.10	2000 ml	19000.00
Spare Rubber Bellow					700.00

186 Burettes, Automatic Zero Mounted on Borosilicate glass Reservoir with Glass Stop Cock Accuracy as per Class 'A' of I.S. 1997 : 1982 complete with rubber pump & work's Certificate & Calibrated @ 27°C

- With auto zero adjustment
- In clear Glass & Straight bore Glass stop cock with screw cap
- Comes with Individual Work Certificate
- Printed in permanent Amber colour
- Comes with Rubber Bellow

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Cap of Reservoir	Price each pc
A	10	0.05	0.02	500 ml	7500.00
B	25	0.10	0.05	2000 ml	7800.00
C	50	0.10	0.05	2000 ml	8000.00
D	100	0.20	0.10	2000 ml	8500.00
Spare Rubber Bellow					700.00





191 Burettes, Automatic Zero Mounted on Reservoir with Teflon Key (Rota Flow) Stop-cock. Accuracy as per Class 'A' of I.S. specification No. 1997 : 1982 with work's Certificate & Calibrated @ 27°C

- With auto zero adjustment
- In clear Glass & Straight bore Glass stop cock with screw cap
- Comes with Individual Work Certificate
- Printed in permanent Amber colour
- Comes with Rubber Bellow


	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Cap of Reservoir	Price each pc
A	10	0.05	0.02	500 ml	7000.00
B	25	0.10	0.05	2000 ml	7300.00
C	50	0.10	0.05	2000 ml	7500.00
D	100	0.20	0.10	2000 ml	8000.00
Spare Rubber Bellow					700.00

192 Burettes, Automatic Zero Amber Colour Mounted on Reservoir with Teflon Key (Rota Flow) Stop-cock. Accuracy as per Class 'A' of I.S. specification No. 1997 : 1982 with work's Certificate & Calibrated @ 27°C

- With auto zero adjustment
- In Amber Glass with Rotaflow PTFE key, lubricate free stop cock
- Comes with Individual Work Certificate
- Printed in White colour
- Comes with Rubber Bellow

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Cap of Reservoir	Price each pc
A	10	0.05	0.02	500 ml	18000.00
B	25	0.10	0.05	2000 ml	18500.00
C	50	0.10	0.05	2000 ml	19500.00
Spare Rubber Bellow					700.00





196 Cylinders Graduated Single Metric Scale with Pour out Accuracy as per **Class 'B'** of I.S. 878 :1975


- Comes with **Round base**
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	5	0.1	0.1	2 Pc	155.00
B	10	0.2	0.2	2 Pc	170.00
C	25	0.5	0.5	2 Pc	230.00
D	50	1.0	1.0	2 Pc	250.00
E	100	1.0	1.0	1 Pc	310.00
F	250	2.0	2.0	1 Pc	650.00
G	500	5.0	5.0	1 Pc	950.00
H	1000	10.0	10.0	1 Pc	1470.00


197 Cylinders Graduated Single Metric Scale with Pour out Accuracy as per **Class 'B'** of I. S. 878 :1975 Hexagonal Base

- Comes with **Hexagonal base**
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	5	0.1	0.1	2 Pc	230.00
B	10	0.2	0.2	2 Pc	255.00
C	25	0.5	0.5	2 Pc	300.00
D	50	1.0	1.0	2 Pc	335.00
E	100	1.0	1.0	1 Pc	400.00
F	250	2.0	2.0	1 Pc	710.00
G	500	5.0	5.0	1 Pc	950.00
H	1000	10.0	10.0	1 Pc	1470.00
I	2000	20.0	20.0	1 Pc	3200.00



Cylinders



202 Cylinders Graduated Single Metric Scale with Pour out Accuracy as per **Class 'A'** of I. S. 878 : 1975 with work's Certificate & Calibrated @ 27°C with Hexagonal Base

- Comes with **Hexagonal base**
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration


	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	5	0.1	0.05	2 Pc	390.00
B	10	0.2	0.10	2 Pc	445.00
C	25	0.5	0.25	2 Pc	500.00
D	50	1.0	0.50	2 Pc	550.00
E	100	1.0	0.50	1 Pc	635.00
F	250	2.0	1.00	1 Pc	890.00
G	500	5.0	2.50	1 Pc	1400.00
H	1000	10.0	5.00	1 Pc	1800.00
I	2000	20.0	10.00	1 Pc	3750.00
any special size & sub. div.				price on request	

203 **Tapped Density** Measuring Cylinders with pour out Accuracy as per **USP** (US pharmacopoeia) **Class 'A'** Amber Stain, with work's Certificate & Calibrated @ 20°C

- Comes with Hexagonal base **for tapped density** apparatus
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration




	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	100	1.0	0.5	1 Pc	900.00
B	250	2.0	1.0	1 Pc	1300.00





Cylinders



204 Cylinders Graduated Single Metric Scale with Pour out Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C.

- Comes with Hexagonal base
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration

	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	5	0.1	0.05	2 Pc	425.00
B	10	0.2	0.10	2 Pc	520.00
C	25	0.5	0.15	2 Pc	550.00
D	50	1.0	0.25	2 Pc	650.00
E	100	1.0	0.50	1 Pc	700.00
F	250	2.0	1.00	1 Pc	990.00
G	500	5.0	1.50	1 Pc	1500.00
H	1000	10.0	3.00	1 Pc	1900.00
I	2000	20.0	6.00	1 Pc	3850.00

205 Cylinders Measuring Amber Colour with pour out Class 'A' with work's Certificate & Calibrated @ 27°C with Hexagonal Base


- Comes with Hexagonal base
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in White colour
- With Unique serial number & Certificate of Calibration

	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	5	0.1	0.05	2 Pc	1200.00
B	10	0.2	0.10	2 Pc	1200.00
C	25	0.5	0.25	2 Pc	1500.00
D	50	1.0	0.50	2 Pc	1750.00
E	100	1.0	0.50	1 Pc	1980.00
F	250	2.0	1.00	1 Pc	2600.00
G	500	5.0	2.50	1 Pc	3500.00
H	1000	10.0	5.00	1 Pc	5200.00



J-SIL
www.j-sil.in

Cylinders



206 Cylinders Graduated Single Metric Scale with Interchangeable Stopper Accuracy as per **Class 'B'** of I.S. 878 :1975


- Comes with **Round base & Penny head solid stopper**
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour

	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Joint Size	Pack per box	Price each pc
A	5	0.1	0.1	10/19	2 Pc	380.00
B	10	0.2	0.2	10/19	2 Pc	400.00
C	25	0.5	0.5	14/23	2 Pc	440.00
D	50	1.0	1.0	14/23	2 Pc	500.00
E	100	1.0	1.0	19/26	1 Pc	600.00
F	250	2.0	2.0	24/29	1 Pc	1000.00
G	500	5.0	5.0	24/29	1 Pc	1400.00
H	1000	10.0	10.0	34/35	1 Pc	2100.00

207 Cylinders Graduated Single Metric Scale with interchangeable Stopper Accuracy as per **Class 'B'** of I.S. 878 : 1975 Hexagonal Base


- Comes with **Hexagonal base & Penny head solid stopper**
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour

	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Joint Size	Pack per box	Price each pc
A	5	0.1	0.1	10/19	2 Pc	380.00
B	10	0.2	0.2	10/19	2 Pc	400.00
C	25	0.5	0.5	14/23	2 Pc	440.00
D	50	1.0	1.0	14/23	2 Pc	500.00
E	100	1.0	1.0	19/26	1 Pc	600.00
F	250	2.0	2.0	24/29	1 Pc	1000.00
G	500	5.0	5.0	24/29	1 Pc	1400.00
H	1000	10.0	10.0	34/35	1 Pc	2100.00





Cylinders



212 Cylinders Graduated Single Metric Scale with interchangeable Stopper Accuracy as per Class 'A' of I. S. 878 : 1975 with work's Certificate & Calibrated @ 27°C with Hexagonal Base

- Comes with Hexagonal base & Penny head solid stopper
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration


	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Joint Size	Pack per box	Price each pc
A	5	0.1	0.05	10/19	2 Pc	635.00
B	10	0.2	0.10	10/19	2 Pc	665.00
C	25	0.5	0.25	14/23	2 Pc	800.00
D	50	1.0	0.50	14/23	2 Pc	910.00
E	100	1.0	0.50	19/26	1 Pc	1050.00
F	250	2.0	1.00	24/29	1 Pc	1650.00
G	500	5.0	2.50	24/29	1 Pc	2200.00
H	1000	10.0	5.00	34/35	1 Pc	2900.00

213 Cylinders Graduated Single Metric Scale with interchangeable Stopper Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C.

- Comes with Hexagonal base & Penny head solid stopper
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration



	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Joint Size	Pack per box	Price each pc
A	5	0.1	0.05	10/19	2 Pc	635.00
B	10	0.2	0.10	10/19	2 Pc	665.00
C	25	0.5	0.15	14/23	2 Pc	800.00
D	50	1.0	0.25	14/23	2 Pc	910.00
E	100	1.0	0.50	19/26	1 Pc	1050.00
F	250	2.0	1.00	24/29	1 Pc	1650.00
G	500	5.0	1.50	24/29	1 Pc	2200.00
H	1000	10.0	3.00	34/35	1 Pc	2900.00



J-SIL
www.j-sil.in

Cylinders



216 Cylinders Measuring for **Suspensibility Test**

- Comes with Round base & Penny head solid stopper
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour

	Description	Price/ pc
A	250 ml Capacity with Stopper distance between 0 and 250ml mark shall be 210 ± 5 mm. There shall be an extra marks at 25 ml. to easy withdrawal of nine-tenths of the suspension	2500.00

217 Cylinders Measuring for **Suspensibility Test** with **Class 'A'** work's Certificate & Calibrated @ 27°C.

- Comes with Round base & Penny head solid stopper
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration

	Description	Price/ pc
A	250 ml Capacity with Stopper distance between 0 and 250ml mark shall be 210 ± 5 mm. There shall be an extra marks at 25 ml. to easy withdrawal of nine-tenths of the suspension	3000.00



218 Measuring Cylinder Hexagonal (**Polypropylene**)


- Made from polypropylene.
- These autoclavable cylinder have raised prominent graduations for easy reading.

	Cap ml	Least Count	Pack per box	Price per pc
A	10	0.2	24 Pc	36.00
B	25	0.5	12 Pc	54.00
C	50	1	12 Pc	66.00
D	100	1	12 Pc	75.00
E	250	2	6 Pc	144.00
F	500	5	4 Pc	225.00
G	1000	10	3 Pc	440.00





Cylinders



221 Cylinders Colour Comparison Nessler flat Bottom, Graduated as per I.S.4161 :1967



- For colour comparison of water samples
- Manufactured from 3.3 Low Expansion Borosilicate glass

	Cap ml	Mark	Pack per box	Price each pc
A	50	25 & 50 ml.	12 Pc	170.00
B	100	50 & 100 ml.	12 Pc	240.00
C	50	Match Pair	1 Set	350.00
D	100	Match Pair	1 Set	500.00

222 Cylinders Colour Comparison Class 'A' with work's Certificate & Calibrated @ 27°C

- For colour comparison of water samples
- With Unique serial number & Certificate of Calibration

	Cap ml	Mark	Pack per box	Price each pc
A	50	25 & 50 ml.	2 Pc	350.00
B	100	50 & 100 ml.	2 Pc	500.00




226 Crow Receivers Class 'B'

- Manufactured from 3.3 Low Expansion Borosilicate glass
- Graduation in permanent Amber colour
- Crow Receivers provide a sturdy and reliable choice for the collection of distillate in oil and tar tests.

	Cap ml	Pack per box	Price each pc
A	100	1 Pc	650.00



J-SIL
www.j-sil.in



231 Crow Receivers with I/C stopper **Class 'B'**


- Comes with **Round base & Penny head solid stopper**
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- Crow Receivers provide a sturdy and reliable choice for the collection of distillate in oil and tar tests.

	Cap ml	Pack per box	Price each pc
A	100	1 Pc	750.00

236 Crow Receivers with **Class 'A'** work's Certificate & Calibrated @ 27°C

- Comes with Round base
- Manufactured from 3.3 Low Expansion Borosilicate
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration
- Crow Receivers provide a sturdy and reliable choice for the collection of distillate in oil and tar tests.

	Cap ml	Tol. $\pm$ ml	Pack per box	Price each pc
A	100	0.40	1 Pc	850.00



241 Crow Receivers with I/C stopper **Class 'A'** with work's Certificate & Calibrated @ 27°C

- Comes with **Round base & Penny head solid stopper**
- Manufactured from 3.3 Low Expansion Borosilicate glass
- Graduation in permanent Amber colour
- With Unique serial number & Certificate of Calibration
- Crow Receivers provide a sturdy and reliable choice for the collection of distillate in oil and tar tests.

	Cap ml	Tol. $\pm$ ml	Pack per box	Price each pc
A	100	0.40	1 Pc	1050.00



Volumetric Flask



246 Volumetric Flask I/C Penny Head Solid stopper Accuracy as per Class 'B' of I.S. 915 : 2006

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Comes with Penny Head Solid glass stopper
- 3 to 20ml Supplied with TRAPEZOIDAL shape for better stability

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	1..	10/15	10/19	0.02	6 pc	200.00
B	2	10/15	10/19	0.02	6 pc	200.00
C	3	10/15	10/19	0.03	2 pc	200.00
D	4	10/15	10/19	0.03	2 pc	200.00
E	5	10/15	10/19	0.04	2 pc	160.00
F	10	10/15	10/19	0.04	2 pc	160.00
G	20	10/15	10/19	0.06	2 pc	170.00
H	25	10/15	10/19	0.06	2 pc	160.00
I	50	10/15	12/21	0.08	2 pc	170.00
J	100	14/15	14/23	0.15	2 pc	180.00
K	200	14/15	14/23	0.20	2 pc	260.00
L	250	14/15	14/23	0.20	2 pc	270.00
M	500	19/20	19/26	0.30	2 pc	440.00
N	1000	19/20	24/29	0.80	2 pc	725.00
O	2000	24/29	29/32	1.20	2 pc	1900.00

- Test Tube Shape

247 Volumetric Flask Amber Colour Class 'B'

- Manufactured from 3.3 Low Expansion Borosilicate Amber glass
- Comes with Penny Head Solid glass stopper
- 5 to 20ml Supplied with TRAPEZOIDAL shape for better stability


	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	1..	10/15	10/19	0.02	6 pc	350.00
B	2	10/15	10/19	0.02	6 pc	350.00
C	5	10/15	10/19	0.04	2 pc	350.00
D	10	10/15	10/19	0.04	2 pc	350.00
E	20	10/15	10/19	0.06	2 pc	360.00
F	25	10/15	10/19	0.06	2 pc	360.00
G	50	10/15	12/21	0.08	2 pc	380.00
H	100	14/15	14/23	0.15	2 pc	400.00
I	200	14/15	14/23	0.20	2 pc	450.00
J	250	14/15	14/23	0.20	2 pc	470.00
K	500	19/20	19/26	0.30	2 pc	730.00
L	1000	19/20	24/29	0.80	2 pc	1200.00
M	2000	24/29	29/32	1.20	2 pc	2600.00

- Test Tube Shape



J-SIL
www.j-sil.in

Volumetric Flask




251 Volumetric Flask with LDPE Stopper Accuracy as per Class 'B' of I.S.915 :2006

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Comes with PP Stopper
- 3 to 20ml Supplied with TRAPEZOIDAL shape for better stability

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	1..	10/15	10/19	0.02	6 pc	200.00
B	2	10/15	10/19	0.02	6 pc	200.00
C	3	10/15	10/19	0.03	2 pc	200.00
D	4	10/15	10/19	0.03	2 pc	200.00
E	5	10/15	10/19	0.04	2 pc	160.00
F	10	10/15	10/19	0.04	2 pc	160.00
G	20	10/15	10/19	0.06	2 pc	170.00
H	25	10/15	10/19	0.06	2 pc	160.00
I	50	10/15	12/21	0.08	2 pc	170.00
J	100	14/15	14/23	0.15	2 pc	180.00
K	200	14/15	14/23	0.20	2 pc	260.00
L	250	14/15	14/23	0.20	2 pc	270.00
M	500	19/20	19/26	0.30	2 pc	440.00
N	1000	19/20	24/29	0.80	2 pc	725.00
O	2000	24/29	29/32	1.20	2 pc	1900.00

- • Test Tube Shape

Weight Ring for Volumetric Flask



Price on Request





Volumetric Flask



256 Volumetric Flask I/C Penny Head Solid stopper Accuracy as per of I.S. 915 : 2006, with **Class 'A'** work's Certificate & Calibrated @27°C


- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Comes with Penny Head Solid glass stopper
- 2 to 20 and 200ml Supplied with TRAPEZOIDAL shape for better stability
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	1..	10/15	10/19	0.01	2 pc	500.00
B	2	10/15	10/19	0.01	2 pc	500.00
C	3	10/15	10/19	0.02	2 pc	500.00
D	4	10/15	10/19	0.02	2 pc	500.00
E	5	10/15	10/19	0.02	2 pc	450.00
F	10	10/15	10/19	0.02	2 pc	450.00
G	20	10/15	10/19	0.03	2 pc	500.00
H	25	10/15	10/19	0.03	2 pc	460.00
I	50	10/15	12/21	0.04	2 pc	465.00
J	100	14/15	14/23	0.06	2 pc	475.00
K	200	14/15	14/23	0.10	2 pc	620.00
L	250	14/15	14/23	0.10	2 pc	625.00
M	500	19/20	19/26	0.15	2 pc	900.00
N	1000	19/20	24/29	0.20	2 pc	1325.00
O	2000	24/29	29/32	0.60	2 pc	2700.00
P	900	19/20	19/26	0.20	2 pc	3000.00
Q	5000	34/35	34/35	2.00	1 pc	10000.00
Any other size						price on request

- • Test Tube Shape

👉 For Wide Mouth Volumetric Flask, Kindly refer Cat. No. 267

Safety Ring



Refer Cat. No. RVF



J-SIL
www.j-sil.in

Volumetric Flask



257 Volumetric Flask Amber Colour Class 'A' with work's Certificate & Calibrated @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate Amber glass
- Comes with Penny Head Solid glass stopper
- 5 to 20 and 200ml Supplied with TRAPEZOIDAL shape for better stability
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	1..	10/15	10/19	0.01	2 pc	750.00
B	2	10/15	10/19	0.01	2 pc	750.00
C	5	10/15	10/19	0.02	2 pc	700.00
D	10	10/15	10/19	0.02	2 pc	660.00
E	20	10/15	10/19	0.03	2 pc	650.00
F	25	10/15	10/19	0.03	2 pc	610.00
G	50	10/15	12/21	0.04	2 pc	625.00
H	100	14/15	14/23	0.06	2 pc	630.00
I	200	14/15	14/23	0.10	2 pc	850.00
J	250	14/15	14/23	0.10	2 pc	850.00
K	500	19/20	19/26	0.15	2 pc	1120.00
L	1000	19/20	24/29	0.20	2 pc	1710.00
M	2000	24/29	29/32	0.60	2 pc	2990.00
Any other size						price on request


- • Test Tube Shape

For Wide Mouth Volumetric Flask, Kindly refer Cat. No. 266

261 Volumetric Flask with Glass stopper, Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Comes with Penny Head Solid glass stopper
- 5 to 20 and 200ml Supplied with TRAPEZOIDAL shape for better stability
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	5	10/15	10/19	0.02	2 pc	510.00
B	10	10/15	10/19	0.02	2 pc	510.00
C	20	10/15	10/19	0.03	2 pc	550.00
D	25	10/15	10/19	0.03	2 pc	520.00
E	50	10/15	12/21	0.05	2 pc	540.00
F	100	14/15	14/23	0.08	2 pc	550.00
G	200	14/15	14/23	0.10	2 pc	750.00
H	250	14/15	14/23	0.12	2 pc	760.00
I	500	19/20	19/26	0.15	2 pc	960.00
J	1000	19/20	24/29	0.30	2 pc	1410.00
K	2000	24/29	29/32	0.50	2 pc	2800.00



For Wide Mouth Volumetric Flask, Kindly refer Cat. No. 268



Volumetric Flask



ASTM

262 Volumetric Flask Amber Colour with Glass stopper, Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C

- Manufactured from 3.3 Low Expansion Borosilicate Amber glass
- Comes with Penny Head Solid glass stopper
- 5 to 20ml Supplied with TRAPEZOIDAL shape for better stability
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	5	10/15	10/19	0.02	2 pc	770.00
B	10	10/15	10/19	0.02	2 pc	770.00
C	20	10/15	10/19	0.03	2 pc	750.00
D	25	10/15	10/19	0.03	2 pc	700.00
E	50	10/15	12/21	0.05	2 pc	740.00
F	100	14/15	14/23	0.08	2 pc	750.00
G	200	14/15	14/23	0.10	2 pc	910.00
H	250	14/15	14/23	0.12	2 pc	935.00
I	500	19/20	19/26	0.15	2 pc	1310.00
J	1000	19/20	24/29	0.30	2 pc	1790.00
K	2000	24/29	29/32	0.50	2 pc	3100.00

☞ For Wide Mouth Volumetric Flask, Kindly refer Cat. No. 269

266 Volumetric Flask Wide Mouth, Amber Colour Class 'A' with work's Certificate & Calibrates @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate Amber glass
- Comes with Penny Head Solid glass stopper
- Wide mouth to put bigger size of Tablet without break
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	20	14/15	14/23	0.03	2 pc	700.00
B	25	14/15	14/23	0.03	2 pc	725.00
C	50	14/15	14/23	0.04	2 pc	680.00
D	100	19/20	14/23	0.06	2 pc	750.00
E	200	19/20	19/26	0.10	2 pc	1000.00
F	250	19/20	19/26	0.10	2 pc	1020.00
G	10	14/15	14/23	0.03	2 pc	800.00



J-SIL
www.j-sil.in

Volumetric Flask



267 Volumetric Flask **Wide Mouth** I/C Stopper Accuracy as per **Class 'A'** with work's Certificate & Calibrated @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Comes with Penny Head Solid glass stopper
- Wide mouth to put bigger size of Tablet without break
- Each flask With Unique serial number & Certificate of Calibration


	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	25	14/15	14/23	0.03	2 pc	460.00
B	50	14/15	14/23	0.04	2 pc	480.00
C	100	19/20	19/26	0.06	2 pc	600.00
D	200	19/20	19/26	0.10	2 pc	710.00
E	250	19/20	19/26	0.10	2 pc	720.00
AA	20	14/15	14/23	0.03	2 pc	560.00
BB	10	14/15	14/23	0.03	2 pc	600.00

268 Volumetric Flask **Wide Mouth** I/C Stopper Accuracy as per **USP (US Pharmacopoeia) Class 'A'** with work's Certificate & Calibrated @ 20°C.

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Comes with Penny Head Solid glass stopper
- Wide mouth to put bigger size of Tablet without break
- Each flask With Unique serial number & Certificate of Calibration



	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	25	14/15	14/23	0.03	2 pc	550.00
B	50	14/15	14/23	0.05	2 pc	580.00
C	100	19/20	19/26	0.08	2 pc	600.00
D	250	19/20	19/26	0.12	2 pc	710.00
AA	20	14/15	14/23	0.03	2 pc	550.00
E	10	14/15	14/23	0.03	2 pc	600.00



Volumetric Flask



ASTM

269 Volumetric Flask Wide Mouth Amber Colour I/C Stopper Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C



- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Comes with Penny Head Solid glass stopper
- Wide mouth to put bigger size of Tablet without break
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Stopper Size (Old)	Stopper Size (New)	Tol ± ml.	Pack per box	Price each pc
A	20	14/15	14/23	0.03	2 pc	650.00
B	25	14/15	14/23	0.03	2 pc	650.00
C	50	14/15	14/23	0.05	2 pc	680.00
D	100	19/20	19/26	0.08	2 pc	780.00
E	200	19/20	19/26	0.10	2 pc	850.00
F	250	19/20	19/26	0.12	2 pc	850.00
G	10	14/15	14/23	0.03	2 pc	900.00

271 Sugar Flask, with two Marks, Class 'B'

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- With two Graduation mark

	Cap. ml.	Marked at	Tol ± ml.	Pack per box	Price each pc
A	50	50 & 55	0.10	2Pc	500.00
B	100	100 & 110	0.15	2Pc	600.00
C	200	200 & 220	0.25	2Pc	700.00



276 Sugar Flask, with two Marks, Class 'A' with work's Certificate & Calibrated @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- With two Graduation mark
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Marked at	Tol ± ml.	Pack per box	Price each pc
A	50	50 & 55	0.06	2 Pc	700.00
B	100	100 & 110	0.10	2 Pc	750.00
C	200	200 & 220	0.15	2 Pc	820.00



J-SIL
www.j-sil.in



281 Kohlrausch Flask (Mud) Class 'B'

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- It is for Sugar Examination

	Cap ml	Tol ± ml	Pack per box	Price each pc
A	50	0.08	2 Pc	450.00
B	100	0.15	2 Pc	550.00
C	200	0.20	2 Pc	650.00

286 Kohlrausch Flask (Mud) Class 'A' with work's Certificate & Calibrated @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Each flask With Unique serial number & Certificate of Calibration
- It is for Sugar Examination

	Cap ml	Tol ± ml	Pack per box	Price each pc
A	50	0.04	2 Pc	700.00
B	100	0.06	2 Pc	750.00
C	200	0.10	2 Pc	900.00




291 Viscosity Flask

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- With large hemispherical tip, heavy walls and reinforced rim
- Each flask With Unique serial number & Certificate of Calibration

	Cap. ml.	Class	Tol ± ml.	Pack per box	Price each pc
A	60	Class 'B'	0.10	2Pc	600.00
B	60	Class 'A'	0.06	2Pc	850.00





296 Cassia Flask Graduated long neck


- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Determine amount of unsaturated constituents in straw oil or similar wash oils, by sulfonation method

	Cap ml	Cap. x Sub. ml	Tol. $\pm$ ml	Pack per box	Price each pc
A	100	10 x 0.1	0.08	2 Pc	750.00
B	200	20 x 0.2	0.15	2 Pc	850.00

297 Cassia Flask Graduated long neck Class 'A' with work's Certificate & Calibrated @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- Determine amount of unsaturated constituents in straw oil or similar wash oils, by sulfonation method
- Each flask With Unique serial number & Certificate of Calibration

	Cap ml	Cap. x Sub. ml	Tol. $\pm$ ml	Pack per box	Price each pc
A	100	10 x 0.1	0.06	2 Pc	950.00
B	200	20 x 0.2	0.10	2 Pc	990.00




301 Chatelier Flask

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- The flask determines density of powdered materials, such as limes, slag and hydraulic cement.

	Description	Price each pc
	For specific Gravity test-Bulb capacity approx. 250 ml. Neck graduated from 0 to 1 ml. and 18 to 24 ml. in 0.1 ml divisions	3500.00





306 Chatelier Flask with Class 'A' work's Certificate & Calibrated @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- The flask determines density of powdered materials, such as limes, slag and hydraulic cement.
- Each flask With Unique serial number & Certificate of Calibration

	Description	Price each pc
A	For specific Gravity test-bulb capacity approx. 250 ml. Neck graduated from 0 to 1 ml. and 18 to 24 ml. in 0.1 ml. divisions	4000.00

311 Phenol Flask


- Manufactured from 3.3 Low Expansion Borosilicate clear glass
- These flasks are designed for the determination of phenols in essential oils,
- Graduation in permanent Amber colour

	Cap. ml.	Neck Cap.	Sub. Div. ml.	Pack per box	Price each pc
A	45	5	0.10	2 Pc	1350.00
B	150	10	0.10	2 Pc	1600.00
C	200	25	0.20	2 Pc	1800.00

Spatulla



Refer Cat. No.1492, 1494 & 1495





316 Volumetric Pipettes, Transfer Accuracy as per Class 'B' of I.S. 1117 :1975



- Manufactured from 3.3 Low Expansion Borosilicate glass
- They have precision forced jets and bevel ground ends to resist chipping.
- All of the pipettes are individually measured and then marked
- Transfer (Volumetric) pipettes are manufactured from specially selected, heavy wall machine drawn borosilicate accurate bore glass tubing

	Cap ml	Tol ± ml.	Pack per box	Price each pc
A	1	0.02	10 Pc	115.00
B	2	0.02	10 Pc	115.00
C	3	0.03	10 Pc	140.00
D	4	0.03	10 Pc	140.00
E	5	0.03	10 Pc	115.00
F	10	0.04	10 Pc	117.00
G	15	0.05	10 Pc	160.00
H	20	0.05	10 Pc	170.00
I	25	0.06	10 Pc	165.00
J	50	0.08	6 Pc	250.00
K	100	0.12	6 Pc	400.00

321 Volumetric Pipettes, Transfer Accuracy as per Class 'A' of I.S. 1117:1975 with work's Certificate & Calibrated @ 27°C

- Manufactured from 3.3 Low Expansion Borosilicate glass
- They have precision forced jets and bevel ground ends to resist chipping.
- All of the pipettes are individually measured and then marked
- Transfer (Volumetric) pipettes are manufactured from specially selected, heavy wall machine drawn borosilicate accurate bore glass tubing
- Each pipette With Unique serial number & Certificate of Calibration

	Cap ml	Tol ± ml.	Pack per box	Price each pc
A	1	0.01	1 Pc	240.00
B	2	0.01	1 Pc	250.00
C	3	0.015	1 Pc	335.00
D	4	0.015	1 Pc	335.00
E	5	0.02	1 Pc	280.00
F	10	0.02	1 Pc	320.00
G	15	0.02	1 Pc	375.00
H	20	0.03	1 Pc	450.00
I	25	0.03	1 Pc	410.00
J	50	0.04	1 Pc	500.00
K	100	0.06	1 Pc	800.00
N	6	0.02	1 Pc	400.00
Q	7	0.02	1 Pc	400.00
R	8	0.02	1 Pcs	400.00
S	9	0.02	1 Pcs	400.00



Pipettes



322 Volumetric Pipettes, Amber Colour Class 'A' with work's Certificate & Calibrated @ 27°C


- Manufactured from 3.3 Low Expansion Borosilicate Amber glass
- They have precision forced jets and bevel ground ends to resist chipping.
- All of the pipettes are individually measured and then marked
- Transfer (Volumetric) pipettes are manufactured from specially selected, heavy wall machine drawn borosilicate accurate bore glass tubing
- Each pipette With Unique serial number & Certificate of Calibration

	Cap ml	Tol ± ml.	Pack per box	Price each pc
A	1	0.01	1 Pc	800.00
B	2	0.01	1 Pc	800.00
C	3	0.015	1 Pc	950.00
D	4	0.015	1 Pc	950.00
E	5	0.02	1 Pc	800.00
F	10	0.02	1 Pc	825.00
G	15	0.02	1 Pc	850.00
H	20	0.03	1 Pc	880.00
I	25	0.03	1 Pc	900.00
J	50	0.04	1 Pc	1000.00
K	100	0.06	1 Pc	1150.00

323 Volumetric Pipettes, Transfer Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C.

- Manufactured from specially selected, heavy wall machine drawn borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping.
- All of the pipettes are individually measured and then marked
- Each pipette With Unique serial number & Certificate of Calibration

	Cap ml	Tol ± ml.	Pack per box	Price each pc
A	1	0.006	1 Pc	290.00
B	2	0.006	1 Pc	300.00
C	3	0.01	1 Pc	350.00
D	4	0.01	1 Pc	350.00
E	5	0.01	1 Pc	340.00
F	6	0.01	1 Pc	390.00
G	7	0.01	1 Pc	390.00
H	8	0.02	1 Pc	390.00
I	9	0.02	1 Pc	390.00
J	10	0.02	1 Pc	390.00
K	15	0.03	1 Pc	425.00
L	20	0.03	1 Pc	460.00
M	25	0.03	1 Pc	480.00
N	50	0.05	1 Pc	610.00
O	100	0.08	1 Pc	920.00





326 Pipettes, Measuring Serological, Accuracy as per Class 'B' of I. S. 4162 :1985 Colour Coated

- Manufactured from specially selected, medium wall borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping
- All of the pipettes are individually measured and then marked
- The graduation marks are in durable amber stain

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	0.1	0.01	0.01	10 Pc	140.00
B	0.2	0.01	0.01	10 Pc	140.00
C	0.5	0.20	0.01	10 Pc	140.00
D	1.0	0.10	0.01	10 Pc	140.00
E	1.0	0.01	0.02	10 Pc	130.00
F	2.0	0.10	0.02	10 Pc	150.00
G	2.0	0.02	0.06	10 Pc	135.00
H	5.0	0.10	0.10	10 Pc	155.00
I	5.0	0.05	0.10	10 Pc	145.00
J	10.0	0.10	0.15	10 Pc	160.00
K	20.0	0.20	0.20	10 Pc	230.00
L	25.0	0.20	0.20	10 Pc	240.00
M	1	White Graduation		10 Pc	180.00
N	2	White Graduation		10 Pc	180.00
O	5	White Graduation		10 Pc	190.00
P	10	White Graduation		10 Pc	210.00


331 Pipettes, Measuring Serological, Accuracy as per Class 'A' of I. S. 4162 :1985 with work's Certificate & Calibrated @27°C, Colour Coated

- Manufactured from specially selected, medium wall borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping
- All of the pipettes are individually measured and then marked
- The graduation marks are in durable amber stain
- Each pipette With Unique serial number & Certificate of Calibration

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	0.1	0.01	0.006	1 Pc	330.00
B	0.2	0.01	0.006	1 Pc	330.00
C	0.5	0.02	0.006	1 Pc	400.00
D	1.0	0.10	0.006	1 Pc	340.00
E	1.0	0.01	0.006	1 Pc	320.00
F	2.0	0.10	0.01	1 Pc	340.00
G	2.0	0.02	0.01	1 Pc	320.00
H	5.0	0.10	0.03	1 Pc	380.00
I	5.0	0.05	0.03	1 Pc	360.00
J	10.0	0.10	0.05	1 Pc	395.00
K	20.0	0.20	0.1	1 Pc	600.00
L	25.0	0.20	0.1	1 Pc	550.00



Pipettes



ASTM

332 Pipettes, Measuring (Serological) Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C.


- Manufactured from specially selected, heavy wall machine drawn borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping
- All of the pipettes are individually measured and then marked
- The graduation marks are in durable amber stain
- Each pipette With Unique serial number & Certificate of Calibration


	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	0.1	0.01	0.006	1 Pc	380.00
B	0.2	0.01	0.006	1 Pc	380.00
C	0.5	0.02	0.006	1 Pc	380.00
D	1.0	0.01	0.006	1 Pc	360.00
E	2.0	0.02	0.01	1 Pc	370.00
F	5.0	0.05	0.02	1 Pc	400.00
G	10.0	0.10	0.03	1 Pc	450.00
H	20.0	0.20	0.05	1 Pc	580.00
I	25.0	0.20	0.05	1 Pc	550.00

336 Pipettes, Graduated Serological Amber Colour Class 'A' with work's Certificate & Calibrated @ 27°C

- Manufactured from specially selected, heavy wall machine drawn borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping
- All of the pipettes are individually measured and then marked
- The graduation marks are in durable amber stain
- Each pipette With Unique serial number & Certificate of Calibration
- Serological pipette delivery from zero mark to the tip.

	Cap. ml.	Sub. Div. ml.	Tol ± ml.	Pack per box	Price each pc
A	0.1	0.01	0.006	1 Pc	1200.00
B	0.2	0.01	0.006	1 Pc	1200.00
C	0.5	0.02	0.006	1 Pc	1200.00
D	1.0	0.01	0.006	1 Pc	1200.00
E	2.0	0.02	0.01	1 Pc	1200.00
F	5.0	0.05	0.03	1 Pc	1250.00
G	10.0	0.10	0.05	1 Pc	1300.00
H	25.0	0.20	0.1	1 Pc	1500.00





341 Pipettes, Measuring (Mohr Type) Accuracy as per Class 'B' of I.S. 4162 : 1985 (These Pipettes are graduated for delivery from zero mark to the graduated mark) Colour Coated


- Manufactured from specially selected, heavy wall machine drawn borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping
- All of the pipettes are individually measured and then marked
- The graduation marks are in durable amber stain
- Mohr type pipette delivery from zero mark to the last graduation mark


	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	0.1	0.01	0.01	10 Pc	140.00
B	0.2	0.01	0.01	10 Pc	140.00
C	0.5	0.20	0.01	10 Pc	140.00
D	1.0	0.10	0.01	10 Pc	140.00
E	1.0	0.01	0.02	10 Pc	130.00
F	2.0	0.10	0.02	10 Pc	150.00
G	2.0	0.02	0.06	10 Pc	135.00
H	5.0	0.10	0.10	10 Pc	155.00
I	5.0	0.05	0.10	10 Pc	145.00
J	10.0	0.10	0.15	10 Pc	160.00
K	20.0	0.20	0.20	10 Pc	230.00
L	25.0	0.20	0.20	10 Pc	240.00

346 Pipettes, Measuring (Mohr Type) Accuracy as per Class 'A' of I. S. 4162 : 1985 with work's Certificate & Calibrated @ 27°C

- Manufactured from specially selected, heavy wall machine drawn borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping
- All of the pipettes are individually measured and then marked
- The graduation marks are in durable amber stain
- Mohr type pipette delivery from zero mark to the last graduation mark
- Each pipette With Unique serial number & Certificate of Calibration

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	0.1	0.01	0.006	1 Pc	330.00
B	0.2	0.01	0.006	1 Pc	330.00
C	0.5	0.02	0.006	1 Pc	400.00
D	1.0	0.10	0.006	1 Pc	340.00
E	1.0	0.01	0.006	1 Pc	320.00
F	2.0	0.10	0.01	1 Pc	340.00
G	2.0	0.02	0.01	1 Pc	320.00
H	5.0	0.10	0.03	1 Pc	380.00
I	5.0	0.05	0.03	1 Pc	360.00
J	10.0	0.10	0.05	1 Pc	395.00
K	20.0	0.20	0.10	1 Pc	600.00
L	25.0	0.20	0.10	1 Pc	550.00





ASTM



347 Pipettes, Measuring (Mohr Type) Accuracy as per USP (US Pharmacopoeia) Class 'A' with work's Certificate & Calibrated @ 20°C.

- Manufactured from specially selected, heavy wall machine drawn borosilicate 3.3 accurate bore glass tubing
- They have precision forced jets and bevel ground ends to resist chipping
- All of the pipettes are individually measured and then marked
- The graduation marks are in durable amber stain
- Mohr type pipette delivery from zero mark to the last graduation mark
- Each pipette With Unique serial number & Certificate of Calibration

	Cap. ml.	Sub. Div. ml.	Tol $\pm$ ml.	Pack per box	Price each pc
A	0.1	0.01	0.006	1 Pc	380.00
B	0.2	0.01	0.006	1 Pc	380.00
C	0.5	0.02	0.006	1 Pc	380.00
D	1.0	0.01	0.006	1 Pc	360.00
E	2.0	0.02	0.01	1 Pc	370.00
F	5.0	0.05	0.02	1 Pc	400.00
G	10.0	0.10	0.03	1 Pc	450.00
H	20.0	0.20	0.05	1 Pc	580.00
I	25.0	0.20	0.05	1 Pc	550.00

356 Microliter Pipettes Lambda, Class 'A' with work's Certificate & Calibrated @ 27°C

- Each pipette With Unique serial number & Certificate of Calibration

	Cap. UL lambda	Tol $\pm$ UL %	Pack per box	Price each pc
A	1	0.50	1 Pc	650.00
B	2	0.50	1 Pc	650.00
C	3	0.50	1 Pc	650.00
D	4	0.50	1 Pc	650.00
E	5	0.25	1 Pc	650.00
F	10	0.25	1 Pc	650.00
G	25	0.25	1 Pc	650.00
H	50	0.15	1 Pc	650.00
I	100	0.15	1 Pc	650.00
J	200	0.10	1 Pc	650.00
K	250	0.10	1 Pc	650.00
L	500	0.10	1 Pc	650.00






366 Milk Test Pipettes, Bacteriological

	Cap. ml.	Graduate at	Tol $\pm$ ml.	Pack per box	Price each pc
A	1.1	0.5, 1 & 1.1	0.025	10 Pc	135.00
B	2.2	1, 2, 2.1 & 2.2	0.04	10 Pc	140.00
C	11.0	10 & 11	0.20	10 Pc	200.00

371 Milk Test Pipettes, Gerber method

	Cap. ml.	To deliver	Tol $\pm$ ml.	Pack per box	Price each pc
A	1.0	Amyl Alcohol	0.05	10 Pc	175.00
B	10.0	Sulphuric	0.05	10 Pc	180.00
C	10.75	Milk	0.03	10 Pc	200.00



372 Pipette Stand (Horizontal), Polypropylene

- Detachable stand can hold 12 pipettes, 6 on each side, in a horizontal position
- Easy disassembling for space saving storage


Cat No.	Capacity	Pack per box	Price per pc
A	For 12 pipettes	6 pc	360.00

373 Pipette Stand (Vertical), Polypropylene

- Holds maximum 28 pipettes and rotates on a central vertical axis for convenient selection of any pipette
- Liquids draining off the pipettes are retained on the top of the stationary base by the outer ring
- Easy disassembling & reassembling

Cat No.	Capacity	Pack per box	Price per pc
A	For 28 pipettes	1 pc	1200.00







376 Centrifuge Tube, Plain with Rim Conical Bottom

	Cap. ml.	OD x Height mm	Pack per box	Price each pc
A	10	16 x 110	25 Pc	45.00
B	15	17 x 120	20 Pc	50.00
C	25	24 x 120	20 Pc	140.00
D	50	28 x 135	12 Pc	230.00

381 Centrifuge Tube Plain I/C Stopper

	Cap. ml.	OD x Height mm	Pack per box	Price each pc
A	10	16 x 125	25 Pc	190.00
B	15	17 x 135	20 Pc	250.00
C	25	24 x 135	20 Pc	350.00
D	50	28 x 150	12 Pc	450.00




386 Centrifuge Tube, **Graduated** with Rim Conical Bottom

	Cap. ml.	OD x Height mm	Pack per box	Price each pc
A	10	16 x 110	25 Pc	135.00
B	15	17 x 120	20 Pc	140.00
C	25	24 x 120	20 Pc	280.00
D	50	28 x 135	12 Pc	290.00

391 Centrifuge Tube, **Graduated &** I/C Stopper

	Cap. ml.	OD x Height mm	Pack per box	Price each pc
A	10	16 x 125	25 Pc	300.00
B	15	17 x 135	20 Pc	350.00
C	25	24 x 135	20 Pc	450.00
D	50	28 x 150	12 Pc	470.00
E	15	17 x 130	Screw Cap	450.00





396 Centrifuge Tube, Round bottom graduated

	Cap ml	OD x Height mm	Pack per box	Price each pc
A	12	16 x 110	12 Pc	180.00
B	30	24 x 110	12 Pc	200.00
C	40	28 x 110	12 Pc	250.00
D	50	30 x 110	12 Pc	315.00
E	100	45 x 110	6 Pc	450.00
F	200	55 x 116	6 Pc	700.00

397 Centrifuge Tube Conical Bottom, Polypropylene & Autoclavable

- Clear strong & safe
- Leakproof autoclavable and both printed graduation
- Can withstand forces upto 7500 x G 15ml & 1200 x G for 50 ml

Cat No.	Capacity	Pack per box	Price per pack
A	15 ml	100 pc	3000.00
B	50 ml	100 pc	3300.00



398 Centrifuge Tube Round Bottom

- Made of Polypropylene & Polycarbonate
- Are strong enough to be used upto 50,000xG
- Autoclavable

Cat No.	Capacity	Pack per box	Price per pack
A	16 ml. Tube PP.	50 pc	1800.00
B	50 ml. Tube PP.	50 pc	3000.00
C	16 ml. Tube PC.	50 pc	6600.00
D	50 ml. Tube PC.	50 pc	11000.00




399 OAK Ridge Centrifuge Tubes, Polypropylene

- Can be used upto 50,000xG
- Screw cap is fitted with silicone O-ring making it totally leakproof
- Should not be used below 4°C

Cat No.	Capacity	Pack per box	Price per pack
A	30	12 pc	1150.00
B	50	12 pc	1585.00
C	100	12 pc	3025.00





400 Rack for Micro Centrifuge Tube (Autoclavable)

- Made of **Polypropylene**.
- Made for holding twenty four 1.5 ml. Micro centrifuge tubes.
- Holes are numbered & lettered for easy identification.

Cat No.	Cap. ml.	Pack per box	Price each pc
A	100 places for 1.5 ml (folding)	6 pc	420.00

401 Centrifuge Vessel (Spherical)

- For the determination of water and sediment in petroleum products

	Cap. ml.	Height mm.	Price each pc
	100 Graduated 0 to 1.5 ml. in 0.1, 1.5 to 5 ml in 0.5, 5 to 10 ml, in 1.0 divisions and at 15, 20, 25, 50 and 100 ml.	157	800.00



402 Micro Centrifuge Tube (PP)


- Made of **Polypropylene**.
- Built in plugs make these tubes leakproof.


	Cap. ml.	Pack per box	Price per pack
A	0.5	1000	1700.00
B	1.5	500	1250.00

406 Centrifuge Vessel (Cylindrical)

- Tube used in the determination of water and sediment in crude mineral oils, fuel oils and petroleum products; of precipitate formed by centrifuging definite quantities of steam cylinder stocks, black oils and lubricating oils; and in testing for acidity of distillation residues or hydrocarbon liquids of gasoline/petroleum solvents.

	Cap. ml.	Height mm.	Price each pc
	100 Graduated 0 to 0.5 ml. in 0.05, 0.5 to 2 ml. in 0.1 to 3 ml. in 0.2, 3 to 5 ml. in 0.5, 5 to 10 ml. in 1.0, 10 to 25 ml. in 5.0 and at 50, 75 and 100 ml.	200	700.00






407 Centrifuge Tube Vessel, Class 'A' with work's Certificate & Calibrated @ 27°C

- tube used in the determination of water and sediment in petroleum products; of precipitate formed by centrifuging definite quantities of steam cylinder stocks, black oils and lubricating oils; and in testing for acidity of distillation residues or hydrocarbon liquids of gasoline/petroleum solvents.

		Price each pc
	Graduated 0 to 0.5 ml. in 0.05, 0.5 to 2 ml. in 0.1 to 3 ml. in 0.2, 3 to 5 ml. in 0.5, 5 to 10 ml. in 1.0, 10 to 25 ml. in 5.0 and at 50, 75 and 100 ml.	800.00

408 Test tube with Rim Borosilicate Glass

	Size mm.	Cap. ml	Pack per box	Price each pc
A	10 x 75	3 ml	100 Pc	8.50
B	12 x 75	5 ml	100 Pc	8.50
C	12 x 100	7 ml	100 Pc	12.00
D	15 x 125	13 ml	100 Pc	15.00
E	15 x 150	15 ml	100 Pc	17.00
F	18 x 150	27 ml	100 Pc	18.00
G	25 x 150	55 ml	50 Pc	27.00
H	25 x 200	70 ml	50 Pc	38.00



408A Test Tube W/o Rim Borosilicate Glass

	Size mm.	Cap. ml	Pack per box	Price each pc
A	10 x 75	3 ml	100 Pc	8.50
B	12 x 75	5 ml	100 Pc	8.50
C	12 x 100	7 ml	100 Pc	12.00
D	15 x 125	13 ml	100 Pc	15.00
E	15 x 150	15 ml	100 Pc	17.00
F	18 x 150	27 ml	100 Pc	18.00
G	25 x 150	55 ml	50 Pc	27.00
H	25 x 200	70 ml	50 Pc	38.00
I	16 x 100	12 ml	100 Pc	18.00
J	16 x 150	20 ml	100 Pc	23.00





409 E




409 Tubes-Culture, Media, Round Bottom, with Screw Cap and Teflon Liner.

	Cap. ml.	O.D. x Height mm.	Pack per box	Price each pc
A	5	16 x 75	100 Pc	42.00
B	10	16 x 125	100 Pc	53.00
C	30	25 x 100	100 Pc	65.00
D	60	25 x 200	50 Pc	80.00
E	150	38 x 200	25 Pc	250.00
F	150	Graduated	25 Pc	350.00


410 Tubes-Culture, Media, Flat Bottom, with Screw Cap and Teflon Liner

	Cap. ml.	O.D. x Height mm.	Pack per box	Price each pc
A	5	16 x 50	100 Pc	42.00
B	15	25 x 57	100 Pc	50.00
C	30	25 x 95	100 Pc	65.00
D	10	25 x 50	100 Pc	53.00
E	20	25 x 72	100 Pc	60.00
F	50	25 x 145	50 Pc	75.00



411 Test tube with Rim Graduated

	O.D. mm.	Length mm.	Cap. ml.	Pack per box	Price each pc
A	12	100	5	25 Pc	80.00
B	15	125	10	25 Pc	90.00
C	15	150	15	25 Pc	90.00
D	18	150	20	25 Pc	100.00
E	25	150	30	25 Pc	120.00
F	25	200	50	25 Pc	140.00




Refer Cat. No. 427



Refer Cat. No. 428





412 Glass Rod (For Stirring) Both End Fire Polish

- Made of Borosilicate 3.3. glass

	O.D. x Height mm.	Pack per box	Price each pc
A	7 x 150	50 Pc	25.00
B	8 x 200	50 Pc	34.00
C	8 x 250	50 Pc	35.00
D	9 x 300	50 Pc	60.00

For PTFE Blade Stirrers, Kindly refer Cat. No. 1066 & 1067

413 Tubes-Culture, Media, Flat Bottom, with Screw Cap and Teflon Liner Amber Colour

- Made of Borosilicate 3.3 glass Amber tube


	Cap. ml.	O.D. x Height mm.	Pack per box	Price each pc
A	5	16 x 50	100 Pc	85.00
B	15	25 x 57	100 Pc	110.00
C	30	25 x 95	100 Pc	140.00
D	10	25 x 50	100 Pc	108.00
E	20	25 x 72	100 Pc	135.00



414 Tubes - Culture, Media, Round Bottom, With Screw Cap and Teflon Liner, Amber Colour

	Cap. ml.	O.D. x Height mm.	Pack per box	Price each pc
A	5	16 x 75	100 Pc	85.00
B	10	16 x 125	100 Pc	95.00
C	30	25 x 100	100 Pc	125.00
D	60	25 x 200	50 Pc	165.00
E	150	38 x 200	25 Pc	750.00





416 Test Tube Stoppered with Hexagonal Hollow Stopper



- Made of Borosilicate 3.3 glass

	Cap. ml.	Dia. x Height mm.	Joint Size (Old)	Joint Size (New)	Pack per box	Price each pc
A	5	12 x 100	10/13	10/19	25 pc	110.00
B	10	15 x 125	10/13	10/19	25 pc	125.00
C	15	16 x 160	14/15	14/23	25 pc	140.00
D	20	18 x 160	14/15	14/23	25 pc	145.00
E	30	25 x 150	19/26	19/26	25 pc	155.00
F	40	25 x 175	19/26	19/26	25 pc	215.00
G	50	25 x 200	19/26	19/26	25 pc	220.00
H	75	32 x 200	24/29	24/29	12 pc	300.00
I	100	34 x 200	24/29	24/29	12 pc	340.00

417 Test Tube Stoppered with Penny Head Solid Stopper

- Made of Borosilicate 3.3 glass

	Cap. ml.	Dia. x Height mm.	Joint Size (Old)	Joint Size (New)	Pack per box	Price each pc
A	5	12 x 100	10/13	10/19	25 pc	110.00
B	10	15 x 125	10/13	10/19	25 pc	125.00
C	15	16 x 160	14/15	14/23	25 pc	155.00
D	20	18 x 160	14/15	14/23	25 pc	165.00
E	30	25 x 150	19/26	19/26	25 pc	180.00
F	40	25 x 175	19/26	19/26	25 pc	230.00
G	50	25 x 200	19/26	19/26	25 pc	250.00
H	75	32 x 200	24/29	24/29	12 pc	350.00
I	100	34 x 200	24/29	24/29	12 pc	390.00



Refer Cat. No. 427



Refer Cat. No. 428





421 Test Tubes Stoppered, Graduated with Hexagonal Hollow Stopper

- Made of Borosilicate 3.3 glass
- Printing with permanent Amber colour graduation

	Cap. ml.	Sub. Div. mm	Dia. x Height mm.	Joint Size (Old)	Joint Size (New)	Pack per box	Price each pc
A	10	0.2	15 x 125	10/13	10/19	25 pc	145.00
B	15	0.2	16 x 160	14/15	10/19	25 pc	150.00
C	20	0.5	18 x 160	14/15	14/23	25 pc	165.00
D	30	0.5	25 x 150	19/26	19/26	25 pc	200.00
E	40	1.0	25 x 175	19/26	19/26	25 pc	235.00
F	50	1.0	25 x 200	19/26	19/26	25 pc	240.00
G	100	2.0	34 x 200	24/29	24/29	12 pc	380.00

422 Test Tube Graduated, Stoppered with Penny Head Solid Stopper

- Made of Borosilicate 3.3 glass
- Printing with permanent Amber colour graduation

	Cap. ml.	Sub. Div. mm	Dia. x Height mm.	Joint Size (Old)	Joint Size (New)	Pack per box	Price each pc
A	10	0.2	15 x 125	10/13	10/19	25 pc	160.00
B	15	0.2	16 x 160	14/15	14/23	25 pc	165.00
C	20	0.5	18 x 160	14/15	14/23	25 pc	170.00
D	25	0.5	22 x 160	19/26	19/26	25 pc	220.00
E	30	1.0	25 x 150	19/26	19/26	25 pc	230.00
F	50	1.0	25 x 200	19/26	19/26	25 pc	260.00
G	100	2.0	34 x 200	24/29	24/29	12 pc	410.00



423 Storage Vial (Polypropylene)

- Cap is provided with a built in seal ring making the vial leakproof

Cat No.	Cap. ml.	Pack per box	Price per pack
A	2.0	500 pc	3300.00
B	5.0	500 pc	6000.00





424 Sv5 Vial, Polypropylene


- Have polypropylene cap with built in seal ring making it a leak proof vial

Cat No.	Capacity	Pack per box	Price per pack
A	SV - 5ml	500 pc	5100.00

425 Test Tube With Screw Cap, Polypropylene

- Tubes provide excellent option for RIA, coagulation & bacteriology etc.
- Available with screw cap.

Cat No.	Capacity	Pack per box	Price per pack
A	16 x 100 mm	500	9000.00



426 Blood Sugar Tubes, as per I. S. 3740 : 1966

	Cap. ml.	Dia mm.	Height mm.	Graduated at ml.	Pack per box	Price each pc
A	25	20	215	4.00, 12.5 & 25	25 Pc	175.00





Tubes Stand & Basket

427 Test Tube Stand (3 TIER), Polypropylene & Autoclavable

- Ingenious rack with 3-tier design provide clear view of the tube contents
- Convenient inserting & removal of tubes & Coloured & Opaque
- Withstand sub-freezing temperatures




Cat No.	Size/Capacity	Pack per box	Price per pc
A	13 mm x 31 Tubes	2 pc	330.00
B	16 mm x 31 Tubes	2 pc	330.00
C	20 mm x 20 Tubes	2 pc	480.00
D	25 mm x 12 Tubes	2 pc	330.00
E	25 mm x 18 Tubes	2 pc	480.00
F	32 mm x 12 Tubes	2 pc	540.00
G	25 mm x 36 Tubes	2 pc	960.00
H	13 mm x 18 Tubes	4 pc	225.00
I	15 mm x 18 Tubes	4 pc	225.00
J	18 mm x 12 Tubes	4 pc	225.00
K	13 mm x 62 Tubes	2 pc	660.00
L	16 mm x 62 Tubes	2 pc	660.00
M	20 mm x 40 Tubes	2 pc	960.00
N	13 mm x 48 Tubes	2 pc	510.00
O	13 mm x 72 Tubes	2 pc	660.00


428 Test Tube Baskets (Polypropylene)

- Autoclavable & non-corrosive.
Ideal for rinsing & draining glassware.
- Useful for storage of test tubes or culture tubes during sterilization or incubation.

Cat No.	Size mm	Pack per box	Price per pc
A	140 x 120 x 110 mm	6	225.00
B	160 x 160 x 160 mm	6	300.00



Refer Cat. No. 471



Refer Cat. No. 476



J-SIL
www.j-sil.in

Specific Gravity Bottles



431 Specific Gravity bottle with capillary bore I/C stopper

- These borosilicate glass specific gravity bottles are accurately adjusted at 27°C to determine the density of liquid samples. While the joints have Standard Taper dimensions, the capillary stoppers are not interchangeable since the final adjustment for volume is made on the stopper.

	Cap. ml.	Joint size	Pack per box	Price each pc
A	5	10/19	25 Pc	220.00
B	10	10/19	25 Pc	230.00
C	25	10/19	25 Pc	240.00
D	50	10/19	10 Pc	300.00
E	100	10/19	10 Pc	450.00

436 Specific Gravity bottle with capillary bore I/C Teflon Stopper

- These borosilicate glass specific gravity bottles are accurately adjusted at 27°C to determine the density of liquid samples. While the joints have Standard Taper dimensions, the capillary stoppers are not interchangeable since the final adjustment for volume is made on the stopper.

	Cap. ml.	Joint size	Pack per box	Price each pc
A	5	10/19	25 Pc	360.00
B	10	10/19	25 Pc	370.00
C	25	10/19	25 Pc	400.00
D	50	10/19	10 Pc	450.00
E	100	10/19	10 Pc	525.00



441 Specific Gravity bottle Class 'A' with capillary bore I/C Teflon stopper, with work's Certificate & Calibrated @ 27°C

- These borosilicate glass specific gravity bottles are accurately adjusted at 27°C to determine the density of liquid samples. While the joints have Standard Taper dimensions, the capillary stoppers are not interchangeable since the final adjustment for volume is made on the stopper.

	Cap. ml.	Joint size	Pack per box	Price each pc
A	5	10/19	2 Pc	825.00
B	10	10/19	2 Pc	850.00
C	25	10/19	2 Pc	950.00
D	50	10/19	2 Pc	990.00
E	100	10/19	2 Pc	1100.00





Pyknometer



446 Specific Gravity bottle with Thermometer (Pyknometer)

- Thermometer with Standard Taper No. 10/18 ground joint serves as stopper Side arm is fitted with cap Thermometer range 0° to 50°C in 0.2° divisions; tolerance $\pm 0.2^\circ\text{C}$
- Capacity is within 10% of nominal shown on the flask
- Each flask body is numbered for ready identification

	Cap. ml.	Joint size	Pack per box	Price each pc
A	10	10/19	1 Pc	900.00
B	25	10/19	1 Pc	950.00
C	50	10/19	1 Pc	995.00
D	Spare Thermometer		1 Pc	750.00

447 Specific Gravity bottle with Thermometer (Pyknometer) with Certified Thermometer

- Thermometer with Standard Taper No. 10/18 ground joint serves as stopper Side arm is fitted with cap Thermometer range 0° to 50°C in 0.2° divisions; tolerance $\pm 0.2^\circ\text{C}$
- Capacities are within 10% of nominal capacities
- Each flask body is numbered for ready identification

	Cap. ml.	Pack per box	Price each pc
A	10	1 Pc	2500.00
B	25	1 Pc	2600.00
C	50	1 Pc	2700.00



Refer Cat. No. 61



Refer Cat. No. 66



Beakers



451 Beaker, Euro Design. PP

- Made from Polypropylene Graduated and Autoclavable

Cat No.	Capacity ml.	Pack per box	Price per pc
A	25	12 pc	30.00
B	50	12 pc	30.00
C	100	12 pc	33.00
D	250	12 pc	60.00
E	500	12 pc	90.00
F	1000	6 pc	150.00
G	2000	6 pc	270.00
J	5000	2 pc	750.00

456 Measuring Jugs with handle (Euro Design). PP

- It is Autoclavable and Graduated

Cat No.	Capacity ml.	Pack per box	Price per pc
A	500	12 pc	90.00
B	1000	6 pc	138.00
C	2000	6 pc	234.00
D	3000	6 pc	450.00



466 Beakers, Low Form, with Spout


- Manufactured from 3.3 expansion, low extractable borosilicate glass
- Thick, beaded top, with a spout designed to have excellent pouring characteristics

	Cap. ml.	OD x Height mm	Pack per box	Price each pc
A	1	12x20	25 Pc	75.00
B	5	20x25	25 Pc	80.00
C	10	25x32	25 Pc	85.00
D	25	34x45	10 Pc	90.00





Beakers



471 Beakers low form with spout Double graduated

- Improved mechanical and thermal properties result from the uniform sidewall and bottom thickness design
- Thick, beaded top, with a spout designed to have excellent pouring characteristics
- Suitable for use on a hotplate
- Easy-to-read amber permanent graduated scale for measuring and/or mixing liquids is provided on all sizes from 50 to 2000mL
- Manufactured from 3.3 expansion, low extractable borosilicate glass

	Cap. ml.	OD x Height mm	Pack per box	Price each pc
A	50	42x60	10 Pc	72.00
B	100	50x73	10 Pc	78.00
C	150	56x83	10 Pc	85.00
D	400	80x110	6 Pc	150.00
E	250	68x95	10 Pc	95.00
F	500	83x115	6 Pc	145.00
G	600	90x125	6 Pc	198.00
H	1000	105x150	3 Pc	320.00
I	2000	131x185	1 Pc	650.00
K	1000	For tablet disintegration machine		485.00




472 Beakers, Griffin Low Form Amber Colour with Spout

- Improved mechanical and thermal properties result from the uniform sidewall and bottom thickness design
- Thick, beaded top, with a spout designed to have excellent pouring characteristics
- Suitable for use on a hotplate
- Easy-to-read white graduated scale for measuring and/or mixing liquids is provided on all sizes from 50 to 2000mL
- Manufactured from 3.3 expansion, low extractable borosilicate Amber Glass

	Cap. ml.	OD x Height mm	Pack per box	Price each pc
A	25	34x45	10 Pc	250.00
B	50	42x60	10 Pc	395.00
C	100	50x73	10 Pc	425.00
D	150	56x83	10 Pc	500.00
E	250	68x95	10 Pc	600.00
F	500	83x115	6 Pc	825.00
G	1000	105x150	3 Pc	1900.00
H	2000	131x185	1 Pc	3600.00



J-SIL
www.j-sil.in



474 Beakers Tall Form with Graduation.

- Manufactured from 3.3 expansion, low Extractable borosilicate glass.

	Cap. ml.	Spout	OD x Height mm	Price each pc
A	50	with	38x70	140.00
B	100	with	48x80	150.00
D	250	with	60x120	155.00
E	500	with	80x140	215.00
I	50	without	38x70	140.00
J	100	without	48x80	150.00
L	250	without	60x120	155.00
M	500	without	80x140	215.00




Refer Cat. No.1492



Refer Cat. No.1494



Refer Cat. No.889



Refer Cat. No.1495



Refer Cat. No.1497




Refer Cat. No.1499





Conical Flasks



476 Conical Flasks Erlenmeyer, N.M.

- Tops are reinforced and tooled with a rounded finish, containing more glass to give them maximum mechanical strength
- Body is thick-walled, with a long tapered outside contour to minimize chipping when struck or rubbed together
- All flasks have permanent amber scales to indicate approximate volumes at various levels, useful in measuring and mixing solutions where a high degree of accuracy is not necessary
- Manufactured from 3.3 expansion, low extractable borosilicate glass


	Cap. ml.	Neck x Dia x Height mm	Pack per box	Price each pc
A	5	14x28x45	4 Pc	65.00
B	10	15x37x60	4 Pc	70.00
C	25	22x40x70	4 Pc	75.00
D	50	22x51x82	4 Pc	80.00
E	100	26x64x105	4 Pc	90.00
F	150	28x72x124	4 Pc	110.00
G	250	34x85x140	4 Pc	135.00
H	250 W.M.	50x85x140	2 Pc	225.00
I	500	34x105x180	2 Pc	200.00
J	1000	42x131x220	2 Pc	380.00
K	2000	50x168x280	2 Pc	775.00

👉 W.M. = Wide Mouth

481 Conical Flasks with Screw Cap.


- Used to mix and store culture media, as well as with other culturing techniques.
- Manufactured from 3.3 expansion, low extractable borosilicate glass

	Cap. ml.	Size	Pack per box	Price each pc
A	100	GL25	4 Pc	250.00
B	150	GL32	4 Pc	280.00
C	250	GL45	4 Pc	350.00
D	500	GL45	2 Pc	410.00
E	1000	GL45	2 Pc	750.00



J-SIL
www.j-sil.in


Baffled Flasks



486 Baffled Erlenmeyer Flask

- Used to mix & store culture media, as well as with other Culturing techniques.
- Manufactured from 3.3 expansion, low extractable Borosilicate glass.

	Cap. ml.	Neck	Pack per box	Price each pc
A	250	plain	4	800.00
B	500	plain	2	950.00
C	1000	Plain	2	1300.00
D	250	GL 45	4	990.00
E	500	GL 45	2	1200.00
F	1000	GL 45	2	1800.00



Refer Cat. No. 256



Refer Cat. No. 257





Adapters



496 Adapter Reduction

- Manufactured from 3.3 expansion, low extractable borosilicate glass
- Top size socket in small & bottom cone in size big

	Socket size	Cone size	Price each pc
A	10/19	14/23	220.00
B	14/23	19/26	240.00
C	14/23	24/29	270.00
E	19/26	24/29	280.00
F	19/26	29/32	350.00
G	19/26	34/35	430.00
L	24/29	29/32	390.00
M	24/29	34/35	450.00
N	24/29	40/38	540.00
R	29/32	34/35	560.00
U	34/35	50/58	900.00

501 Adapter Expansion

- Manufactured from 3.3 expansion, low extractable borosilicate glass
- Top size socket is big & bottom size cone in small

	Socket size	Cone size	Price each pc
A	14/23	10/13	220.00
B	19/26	14/23	240.00
C	24/29	19/26	280.00
D	24/29	14/23	270.00
E	29/32	24/29	390.00
H	34/35	29/32	560.00
I	34/35	24/29	450.00



506 Adapter Multiple, Two Necks Parallel

	Socket size	Cone size	Price each pc
A	10/19	10/13	370.00
B	14/23	14/23	380.00
C	19/26	19/26	400.00
D	19/26	24/29	425.00
F	24/29	24/29	610.00



J-SIL
www.j-sil.in

Adapters



511 Adapter Multiple, Two Necks one Vertical and One at Angle

	Socket size	Cone size	Price each pc
A	19/26	19/26	450.00
B	19/26	24/29	525.00

516 Adapter Multiple, Three Necks, Two Parallel and One at Angle

	Parallel Neck Socket	Angle Neck Socket	Cone size	Price each pc
A	14/23	14/23	19/26	550.00
B	19/26	14/23	24/29	700.00
C	19/26	19/23	24/29	750.00
E	19/26	19/26	34/35	850.00



521 Adapter Swan Neck for use with air leak tube

	Socket size	Cone size	Price each pc
A	14/23	14/23	400.00
B	19/26	19/26	460.00
C	19/26	24/29	510.00
D	24/29	24/29	560.00


526 Adapter with 'T' connection

	Socket size	Cone size	Price each pc
A	14/23	14/23	290.00
B	14/23	19/26	300.00
C	14/23	24/29	330.00
D	19/26	19/26	340.00
E	19/26	24/29	350.00





Adapters



531 Cone Adapters, Straight Connection plain

	Cone size	Price each pc
A	10/19	165.00
B	14/23	168.00
C	19/26	170.00
D	24/29	180.00
E	29/32	235.00
F	34/35	330.00

536 Cone Adapters, Straight Connection with Stop-cock

	Socket size	Stop Cock	Price each pc
A	14/23	Glass Stopcock	440.00
B	19/26	Glass Stopcock	440.00
C	24/29	Glass Stopcock	450.00
D	29/32	Glass Stopcock	525.00
E	34/35	Glass Stopcock	620.00
G	19/26	Teflon Stopcock	700.00
I	24/29	Teflon Stopcock	800.00



541 Cone Adapters Angle Connection


	Cone size	Price each pc
A	10/19	165.00
B	14/23	168.00
C	19/26	170.00
D	24/29	180.00
E	29/32	235.00
F	34/35	330.00

546 Cone Adapters, Right Angle Connection with Stop-cock

	Cone size	Stop Cock	Price each pc
A	14/23	Glass Stopcock	440.00
B	19/26	Glass Stopcock	440.00
C	24/29	Glass Stopcock	450.00
D	29/32	Glass Stopcock	525.00
E	34/35	Glass Stopcock	620.00
G	19/26	Teflon Stopcock	700.00
I	24/29	Teflon Stopcock	800.00



Adapters



551 Socket Adapters, Straight Connection with Stop-cock

	Socket size	Price each pc
A	14/23	470.00
B	19/26	480.00
C	24/29	500.00
D	29/32	560.00
E	34/35	610.00

556 Adapters Receiver Plain Bend Short stem

	Socket size	Price each pc
A	10/19	175.00
B	14/23	180.00
C	19/26	190.00
D	24/29	220.00
E	29/32	290.00



561 Adapter Receiver Plain Bend Long stem

	Socket size	Price each pc
A	10/19	185.00
B	14/23	190.00
C	19/26	190.00
D	24/29	230.00

566 Adapter Receiver Straight Delivery

	Socket size	Price each pc
A	14/23	190.00
B	19/26	195.00
C	24/29	222.00





Adapters



571 Adapter Receiver, Plain Bend, Socket to Cone

	Socket size	Cone size	Price each pc
A	14/23	14/23	290.00
B	19/26	19/26	300.00
C	19/26	24/29	300.00
D	24/29	24/29	325.00
E	29/32	29/32	450.00

576 Adapter Receiver, Bent with Vent

	Socket size	Cone size	Price each pc
A	10/19	10/19	290.00
B	14/23	14/23	300.00
C	19/26	19/26	310.00
D	19/26	24/29	320.00
E	24/29	24/29	340.00



581 Adapter Receiver, Bend with Vacuum Connection


	Socket size	Cone size	Price each pc
A	14/23	14/23	440.00
B	14/23	19/26	440.00
C	19/26	19/26	450.00
D	19/26	24/29	450.00
E	24/29	24/29	520.00
F	29/32	29/32	770.00
G	34/35	34/35	1100.00

586 Adapter Receiver, Straight with Vacuum Connection

	Socket size	Cone size	Price each pc
A	19/26	19/26	450.00
B	19/26	24/29	450.00
C	24/29	24/29	520.00



Adapters



591 Adapter Receiver, Bend with Socket Connection

	Socket size	Cone size	Price each pc
A	14/23	14/23	480.00
B	19/26	19/26	485.00
C	19/26	24/29	500.00
D	24/29	24/29	560.00

621 Recovery Bend, Sloping End

	Flask Cone	Cond. Cone	Price each pc
A	14/23	14/23	250.00
B	19/26	19/26	260.00
C	24/29	14/23	260.00
D	24/29	19/26	260.00
E	24/29	24/29	265.00
F	29/32	19/26	350.00
G	29/32	29/32	480.00
H	34/35	19/26	480.00
I	34/35	24/29	500.00
J	34/35	34/45	650.00



626 Recovery Bend, Vertical

	Flask Cone	Cond. Cone	Price each pc
B	19/26	19/26	450.00
C	24/29	19/26	490.00
D	24/29	24/29	500.00





Adapters



631 Plain Still head, Sloping

	Flask Cone	Cond. Cone	Socket size	Price each pc
A	10/13	10/13	10/13	345.00
B	14/23	14/23	14/23	350.00
C	19/26	19/26	14/23	390.00
D	19/26	19/26	19/26	400.00
E	24/29	19/26	14/23	470.00
F	24/29	24/29	14/23	480.00
G	24/29	24/29	19/26	485.00
H	29/32	19/26	14/23	700.00
I	29/32	29/32	14/23	800.00
J	34/35	19/26	14/23	850.00
K	34/35	24/29	14/23	1050.00
L	34/35	34/35	19/26	1400.00

636 Claisen Heads, Sloping


	Flask Cone	Cond. Cone	Socket size	Price each pc
A	14/23	14/23	14/23	680.00
B	14/26	19/26	14/23	690.00
C	24/29	19/26	14/23	800.00
D	24/29	24/29	14/23	850.00



641 Splash Head Adapter

	Socket size	Cone size	Price each pc
A	19/26	24/29	600.00
B	24/29	24/29	650.00





646 Splash Head, Sloping

	Flask Cone	Cond. Cone	Price each pc
A	19/26	19/26	550.00
B	24/29	19/26	580.00
C	24/29	24/29	600.00
D	29/32	19/26	630.00
E	29/32	29/32	750.00
F	34/35	19/26	725.00
G	34/35	24/29	800.00
H	34/35	34/35	950.00

651 Splash Head, Vertical

	Flask Cone	Cond. Cone	Price each pc
A	19/26	19/26	630.00
B	24/29	19/26	670.00
C	24/29	24/29	710.00
D	34/35	24/29	1100.00



656 Steam Distillation Head, Sloping

	Flask Cone	Cond. Cone	Price each pc
A	24/29	19/26	980.00
B	24/29	24/29	990.00
C	34/35	24/29	1100.00


671 Thermometers with Stem and Teflon Cone

	Cone Size	Range	Price each pc
D	14/23	-10 to 110°C	850.00
E	14/23	0 to 250°C	850.00
F	14/23	0 to 360°C	850.00
G	19/26	-10° to 110°C	950.00
H	19/26	0 to 250°C	950.00
I	19/26	0 to 360°C	950.00





Adapters



676 Thermometers Pocket, Closed End

	Cone size	Stem Length mm	Price each pc
A	14/23	20	130.00
B	14/23	40	130.00
C	19/26	50	160.00
D	19/26	60	175.00
E	19/26	80	185.00
F	19/26	100	190.00
G	24/29	60	200.00
H	24/29	80	230.00
I	24/29	100	240.00
J	24/29	120	250.00
K	24/29	140	260.00
L	24/29	150	260.00
M	24/29	200	270.00
N	24/29	300	300.00

681 Air / Stem Inlet Tubes, Straight

	Cone size	Stem Length mm	Price each pc
A	10/19	100	280.00
B	14/23	150	290.00
C	19/26	150	300.00
D	24/29	150	310.00




686 Air / Stem Inlet Tubes, Bent

	Cone size	Stem Length mm	Price each pc
A	10/19	100	280.00
B	14/23	150	290.00
C	19/26	150	300.00
D	24/29	150	310.00



Adapters



691 Drying Tubes, Angled

	Cone size	Pack per box	Price each pc
A	10/19	10 Pc	240.00
B	14/23	10 Pc	280.00
C	19/26	10 Pc	300.00
D	24/29	6 Pc	310.00

696 Stop Cock Adapter for **Vaccum Release**


	Cone size	Pack per box	Price each pc
A	19/26	6 Pc	1000.00
B	24/29	6 Pc	1100.00



701 Stop Cock Adapter for **Vaccum Control**

	Cone size	Pack per box	Price each pc
A	19/26	6 Pc	1000.00
B	24/29	6 Pc	1100.00

Micro Pipette




Refer Cat No. JVR





Condensers




711 Condensers, Leibig Sealed Inner Tube

	Length mm.	Pack per box	Price each pc
B	200	6 Pc	410.00
C	300	6 Pc	450.00
E	400	6 Pc	500.00

726 Condensers, Davies (Double Surface) with I/C joint

- Comes with Interchangeable cone & socket joints
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Highly efficient condenser with large cooling surface
- Inner cooling surface gives 'better' effect

	Length mm.	Socket/Cone size	Pack per box	Price each pc
A	200	B-19 B19	1 Pc	1800.00
B	300	B-19 B19	1 Pc	1900.00
C	400	B-19 B19	1 Pc	2000.00
D	200	B-24 B24	1 Pc	1850.00
E	300	B-24 B24	1 Pc	1950.00
F	400	B-24 B24	1 Pc	2100.00



751 Condensers, Air with I/C joint

- A simple condenser comprising of a plain column with ground socket and cone joint
- Simple condensers for condensation of materials with boiling points above 150°C

	Joints size	Length mm.	Pack per box	Price each pc
A	10/19	100	6 Pc	250.00
B	14/23	200	6 Pc	280.00
C	19/26	200	6 Pc	330.00
D	24/29	300	6 Pc	550.00
H	24/29	400	6 Pc	650.00
I	24/29	600	6 Pc	760.00





756 Condensers, Air Vertical



- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Comes with Interchangeable cone & plain tube

	Cone size	Length mm.	Pack per box	Price each pc
A	24/29	400	10 Pc	550.00
B	24/29	500	10 Pc	580.00
C	24/29	750	10 Pc	620.00
D	24/29	1000	10 Pc	670.00

761 Condensers, Leibig with I/C joint

- Water-jacketed condenser suitable for distillation and reflux applications
- Comes with Interchangeable cone & socket joints
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Joints size	Length mm.	Pack per box	Price each pc
A	14/23	100	6 Pc	400.00
B	14/23	150	6 Pc	450.00
C	19/26	150	6 Pc	480.00
D	19/26	200	4 Pc	500.00
E	19/26	300	6 Pc	515.00
F	19/26	350	6 Pc	620.00
G	19/26	400	6 Pc	650.00
M	24/29	200	4 Pc	500.00
N	24/29	300	6 Pc	530.00
O	24/29	350	6 Pc	640.00
P	24/29	400	6 Pc	675.00
Q	24/29	450	6 Pc	750.00
R	24/29	500	4 Pc	810.00



771 Condensers, Bulb (Allihn) with I/C joint

- Comes with I/C cone & socket joints
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Inner bulb sections provide increased cooling surface

	Joints size	Length mm.	Pack per box	Price each pc
A	19/26	200	1 Pc	680.00
B	24/29	300	1 Pc	720.00
C	24/29	400	1 Pc	850.00
D	24/29	500	1 Pc	1100.00





Condensers



781 Condensers, **Coil** (Graham) with I/C joint


- Comes with Interchangeable cone & socket joints
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Coil provides large cooling surface

	Cone size	Length mm.	Pack per box	Price each pc
A	19/26	150	1 Pc	780.00
B	19/26	300	1 Pc	850.00
C	24/29	150	1 Pc	800.00
D	24/29	300	1 Pc	860.00

796 Condensers, Jacketed Coil with I/C joint

- Comes with Interchangeable cone & socket joints
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Coil condenser has an outer water jacket
- Coil provides large cooling surface

	Joints size	Length mm.	Pack per box	Price each pc
A	19/26	200	1 Pc	1550.00
B	24/29	200	1 Pc	1600.00
C	29/32	200	1 Pc	1900.00
D	19/26	300	1 Pc	1700.00
E	24/29	300	1 Pc	1800.00
F	29/32	300	1 Pc	2100.00




801 Condensers, Coiled Distillage with I/C joint

- Comes with Interchangeable cone & socket joints
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Coil provides large cooling surface

	Joints size	Length mm.	Pack per box	Price each pc
A	19/26	200	1 Pc	750.00
B	19/26	300	1 Pc	880.00
C	19/26	400	1 Pc	1025.00
D	24/29	300	1 Pc	900.00
E	24/29	400	1 Pc	1050.00





811 Condensers, Combined Double Coil


- Comes with Interchangeable cone & socket joints
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- double coil design with a common water inlet
- Suitable for pilot plant work

	Joints size	Length mm.	Pack per box	Price each pc
A	24/29	300	1 Pc	6100.00
B	24/29	400	1 Pc	6700.00
C	34/35	300	1 Pc	6200.00
D	34/35	400	1 Pc	6850.00

816 Condensers, Inland Revenue (Thorpe)

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Joint size	Length mm	Pack per box	Price each pc
A	19/26	200	1 Pc	4000.00
B	24/29	200	1 Pc	4200.00



831 Condensers, Immersion (Cold Finger)

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- For immersion into distillation assemblies etc
- Can be used to regulate refluxing by adjusting the flow of air or water through it


	Cone size	Length mm	Price each pc
A	19/26	100	650.00

Micro Pipette



Refer Cat No. JVE on Page No. 6





846 Columns, Chromatography, plain with Teflon Rota flow Stop-cock



- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Comes with **Rota Flow** PTFE key, lubricate free stop cock

	Bore mm	Length mm.	Pack per box	Price each pc
A	10	150	2 Pc	450.00
B	10	300	2 Pc	490.00
C	18	300	2 Pc	565.00
D	18	450	2 Pc	710.00
E	30	450	1 Pc	830.00
F	30	600	1 Pc	970.00

861 Columns, Chromatography, plain with **sintered disc** and Rotaflow Teflon stop-cock

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Comes with Rotaflow PTFE key, lubricate free stop cock
- Comes with G-1 Integral Sintered disc for filtration

	Bore mm	Length mm.	Pack per box	Price each pc
A	10	150	2 Pc	620.00
B	10	300	2 Pc	640.00
C	18	300	2 Pc	710.00
D	18	450	2 Pc	960.00
E	30	450	1 Pc	1080.00
F	30	600	1 Pc	1225.00



866 Chromatography Columns, with **Integral Sintered**

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Comes with G-1 Integral Sintered disc for filtration
- Comes with socket at top & cone at bottom

	Bore mm	Length mm	Socket size	Cone size	Pack per box	Price each pc
A	10	200	14/23	19/26	2 Pc	550.00
B	10	300	14/23	19/26	2 Pc	570.00
C	10	400	14/23	19/26	2 Pc	650.00
D	18	200	24/29	19/26	2 Pc	670.00
E	18	300	24/29	19/26	2 Pc	720.00
F	18	400	24/29	19/26	2 Pc	850.00
G	30	200	34/35	24/29	1 Pc	900.00
H	30	300	34/35	24/29	1 Pc	1050.00
I	30	400	34/35	24/29	1 Pc	1200.00



Columns & Sprayer



876 Fractionating columns, **Vigrex**


- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Vigreux type with combination of horizontal and angled indentations
- Horizontal indentations ensure intimate vapour and liquid contact
- Angled indentations redistribute liquid from walls to centre of column

	Socket size	Cone mm	Length each	Pack per box	Price each pc
A	24/29	24/29	300	6 Pc	900.00
B	24/29	24/29	400	6 Pc	1020.00
C	24/29	24/29	600	4 Pc	1150.00

881 Fractionating columns, **Plain**

- Borosilicate glass fractionating column is perfect for Raschig rings or Fenske helices.
- Comes with one vigrex at bottom to hold rings

	Socket size	Cone mm	Length each	Pack per box	Price each pc
A	24/29	24/29	250	6 Pc	500.00
B	24/29	24/29	300	6 Pc	550.00
C	24/29	24/29	400	6 Pc	640.00
D	24/29	24/29	500	4 Pc	690.00



886 Sprayer for Chromatography with Rubber Bulb and Clip

- Used as an accessory in thin layer chromatography applications.
- Glass Laboratory Sprayer includes a clear glass conical flask with a sprayer-stopper assembly and rubber bulb.


	Cap ml.	Pack per box	Price each pc
A	25	1 Pc	720.00
B	50	1 Pc	785.00
C	100	1 Pc	875.00
D	Spare Head Only		700.00
E	Spare Bulb for Sprayer		350.00

889 Tong, Flask & Beaker Euro Design

- Available in fibre Glass yarn to resist direct flame
- Hold Flask upto 60 mm OD
- Made of Stainless Steel

	Use	Max Temp	Pack per box	Price each pc
A	For Flask	500°C	4 pc	425.00
B	For Beaker	500°C	4 pc	425.00

For Crucible Tong, Kindly refer Cat. No. 1370



891 Flasks, Round Bottom, Short/Medium neck with I/C joint

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Uniform wall thickness & round shape for ideal heating comes with Interchangeable joint



	Cap. ml.	Socket size	Pack per box	Price each pc
A	5	10/19	4 PC	125.00
B	5	14/23	4 PC	125.00
C	10	10/19	4 PC	130.00
D	10	14/23	4 PC	130.00
E	25	14/23	4 PC	140.00
F	25	19/26	4 PC	140.00
G	50	14/23	4 PC	150.00
H	50	19/26	4 PC	150.00
I	50	24/29	4 PC	150.00
J	100	14/23	4 PC	160.00
K	100	19/26	4 PC	160.00
L	100	24/29	4 PC	160.00
M	150	19/26	4 PC	200.00
N	150	24/29	4 PC	200.00
Q	250	14/23	4 PC	210.00
R	250	19/26	4 PC	210.00
S	250	24/29	4 PC	210.00
T	250	29/32	4 PC	250.00
U	250	34/35	4 PC	300.00
V	500	14/23	2 PC	310.00
W	500	19/26	2 PC	310.00
X	500	24/29	2 PC	310.00
Y	500	29/32	2 PC	360.00
AA	1000	19/26	2 PC	470.00
BB	1000	24/29	2 PC	470.00
CC	1000	29/32	2 PC	495.00
GG	2000	24/29	2 PC	1200.00
HH	2000	29/32	2 PC	1250.00
II	2000	34/35	2 PC	1300.00
KK	3000	24/29	1 PC	2300.00
LL	3000	29/32	1 PC	2350.00
MM	3000	34/35	1 PC	2400.00
PP	5000	24/29	1 PC	4000.00
QQ	5000	29/32	1 PC	4050.00
RR	5000	34/35	1 PC	4050.00





896 Flasks, Flat Bottom, Short/Medium neck with I/C joint

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Uniform wall thickness & flat shape for ideal heating on hot plate comes with Interchangeable joint



	Cap. ml.	Socket size	Pack per box	Price each pc
A	5	10/19	4 PC	135.00
B	5	14/23	4 PC	135.00
C	10	10/19	4 PC	140.00
D	10	14/23	4 PC	140.00
E	25	14/23	4 PC	150.00
F	25	19/26	4 PC	150.00
G	50	14/23	4 PC	165.00
H	50	19/26	4 PC	165.00
I	50	24/29	4 PC	165.00
J	100	14/23	4 PC	175.00
K	100	19/26	4 PC	175.00
L	100	24/29	4 PC	175.00
M	150	19/26	4 PC	220.00
N	150	24/29	4 PC	220.00
Q	250	14/23	4 PC	225.00
R	250	19/26	4 PC	225.00
S	250	24/29	4 PC	220.00
T	250	29/32	4 PC	250.00
U	250	34/35	4 PC	300.00
V	500	14/23	2 PC	325.00
W	500	19/26	2 PC	325.00
X	500	24/29	2 PC	325.00
Y	500	29/32	2 PC	350.00
AA	1000	19/26	2 PC	480.00
BB	1000	24/29	2 PC	480.00
CC	1000	29/32	2 PC	510.00
GG	2000	24/29	2 PC	1250.00
HH	2000	29/32	2 PC	1250.00
II	2000	34/35	2 PC	1300.00
KK	3000	24/29	1 PC	2400.00
LL	3000	29/32	1 PC	2450.00
MM	3000	34/35	1 PC	2500.00
PP	5000	24/29	1 PC	4000.00
QQ	5000	29/32	1 PC	4200.00
RR	5000	34/35	1 PC	4200.00





901 Round Bottom Flasks, Plain Neck with Beaded

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- uniform wall thickness & round shape for ideal heating

	Cap. ml.	Neck Dia	x	Bulb Dia	x	Height mm	Pack per box	Price each pc
E	100	26	x	65	x	110	4 Pc	100.00
F	250	34	x	85	x	145	4 Pc	120.00
G	500	34	x	105	x	170	2 Pc	215.00
H	1000	42	x	131	x	205	2 Pc	365.00
I	2000	50	x	166	x	243	1 Pc	700.00
J	3000	50	x	190	x	310	1 Pc	1300.00

902 Flat Bottom Flasks, Plain Neck with Beaded

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- uniform wall thickness & flat shape for ideal heating

	Cap. ml.	Neck Dia	x	Bulb Dia	x	Height mm	Pack per box	Price each pc
A	100	26	x	65	x	105	4 Pc	120.00
B	250	34	x	85	x	140	4 Pc	160.00
C	500	34	x	105	x	165	2 Pc	250.00
D	1000	42	x	131	x	190	2 Pc	400.00
E	2000	50	x	166	x	237	1 Pc	760.00
F	3000	50	x	186	x	300	1 Pc	1500.00



906 Round Bottom Flasks, Two Neck, Centre Neck and One Parallel side neck

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Fitted with two I/C glass sockets - One vertical (centre neck) - One parallel side neck

	Cap. ml.	Socket size	Side Socket	Pack per box	Price each pc
A	100	24/29	14/23	2 Pc	600.00
B	100	24/29	19/26	2 Pc	600.00
D	250	24/29	19/26	2 Pc	650.00
F	500	24/29	19/26	2 Pc	800.00
H	1000	24/29	19/26	2 Pc	1050.00
J	2000	24/29	19/26	1 Pc	1600.00
K	2000	34/35	24/29	1 Pc	1600.00
L	3000	34/35	19/26	1 Pc	2750.00
M	5000	34/35	24/29	1 Pc	4200.00





911 Round Bottom Flasks, **Three Necks**, Centre Neck and Two Parallel Side Necks

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Fitted with two I/C glass sockets - One vertical (centre neck) - two parallel side neck

	Cap. ml.	Centre Socket	Side Sockets	Pack per box	Price each pc
C	100	24/29	14/23	2 Pc	700.00
D	100	24/29	19/26	2 Pc	700.00
G	250	24/29	14/23	2 Pc	750.00
H	250	24/29	19/26	2 Pc	750.00
M	500	24/29	14/23	2 Pc	880.00
N	500	24/29	19/26	2 Pc	880.00
S	1000	24/29	19/26	2 Pc	1200.00
T	1000	34/35	24/29	2 Pc	1200.00
U	2000	24/29	19/26	1 Pc	1750.00
W	2000	34/35	24/29	1 Pc	1750.00
X	3000	34/35	24/29	1 Pc	3100.00
Y	5000	34/35	24/29	1 Pc	4550.00
Z	10000	34/35	24/29	1 Pc	8350.00

916 Round Bottom Flasks **Four Neck** Centre Neck and Three Parallel Side Necks

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Fitted with two I/C glass sockets - One vertical (centre neck) - three parallel side neck
- one side free for clamping the flask

	Cap. ml.	Centre Socket	Side Sockets	Pack per box	Price each pc
A	100	B-34	B-24	2 Pc	850.00
B	250	B-34	B-24	2 Pc	900.00
C	500	B-34	B-24	2 Pc	930.00
D	1000	B-34	B-24	2 Pc	1300.00
E	2000	B-34	B-24	1 Pc	2000.00
F	3000	B-34	B-24	1 Pc	3700.00
G	5000	B-34	B-24	1 Pc	5430.00
H	10000	B-34	B-24	1 Pc	10800.00
I	20000	B-34	B-24	1 Pc	19000.00





918 Round Bottom Flasks, 4 neck with **Thermometer Pocket** & one stopper for side neck

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Fitted with two I/C glass sockets - One vertical (centre neck) - three parallel side neck
- one side free for clamping the flask

	Cap. ml.	Centre Socket	Side Sockets	Pack per box	Price each pc
A	100	B-34	B-24	2 Pc	1050.00
B	250	B-34	B-24	2 Pc	1100.00
C	500	B-34	B-24	2 Pc	1150.00
D	1000	B-34	B-24	2 Pc	1500.00
E	2000	B-34	B-24	1 Pc	2200.00
F	3000	B-34	B-24	1 Pc	3900.00
G	5000	B-34	B-24	1 Pc	5680.00
H	10000	B-34	B-24	1 Pc	11100.00
I	20000	B-34	B-24	1 Pc	19300.00

919 Flask Stand Autoclavable. PP

- Made of polypropylene.
- Gives stairs like top view.
- Every step provides excellent support to the round bottom flask.

	Capacity	Pack per box	Price per pc
A	Upto 10 ltr	6 pcs.	300.00



920 Teflon Bush / Gland with O Ring

	Size	Bore mm	Price each pc
A	B-19	6	Rate on Request
B	B-24	8	Rate on Request
C	B-34	8	Rate on Request
D	B-34	10	Rate on Request
E	B-40	10	Rate on Request
F	B-45	12	Rate on Request
G	B-50	12	Rate on Request
H	B-55	12	Rate on Request





921 Kjeldhal Flasks, Long Neck

- For the determination of nitrogen by the Kjeldahl method
- comes with I/C joint for adapter fitting

	Cap. ml.	Socket size	Pack per box	Price each pc
A	50	24/29	2 Pc	350.00
B	100	24/29	2 Pc	400.00
C	300	24/29	2 Pc	550.00
D	500	24/29	2 Pc	680.00
E	800	24/29	2 Pc	750.00

926 Kjeldhal Flasks, Long Neck Plain

- Suitable for the determination of nitrogen etc
- comes with plain neck, without joint

	Cap. ml.	Pack per box	Price each pc
A	10	2 Pc	110.00
B	30	2 Pc	130.00
C	100	2 Pc	230.00




Refer Cat. No. 1161



Refer Cat. No. 1166







931 Conical Flasks with I/C joint

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- comes with I/C joint for fitting

	Cap. ml.	Socket size	Pack per box	Price each pc
A	5	10/19	4 Pc	120.00
B	5	14/23	4 Pc	120.00
C	10	10/19	4 Pc	130.00
D	10	14/23	4 Pc	130.00
E	25	10/19	4 Pc	140.00
F	25	14/23	4 Pc	140.00
G	25	19/26	4 Pc	140.00
H	50	14/23	4 Pc	160.00
I	50	19/26	4 Pc	160.00
J	50	24/29	4 Pc	160.00
K	100	14/23	4 Pc	190.00
L	100	19/26	4 Pc	190.00
M	100	24/29	4 Pc	190.00
N	150	19/26	4 Pc	210.00
O	150	24/29	4 Pc	210.00
Q	250	19/26	4 Pc	230.00
R	250	24/29	4 Pc	230.00
S	250	34/35	4 Pc	310.00
T	500	19/26	2 Pc	330.00
U	500	24/29	2 Pc	330.00
X	1000	24/29	2 Pc	700.00
Z	1000	34/35	2 Pc	750.00
AA	2000	24/29	2 Pc	1150.00
CC	2000	34/35	2 Pc	1200.00




Refer Cat. No.476



Refer Cat. No.481






936 Stopped Conical Flasks

- comes with I/C joint hollow stopper
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Joint size	Pack per box	Price each pc
A	10	10/13	4 Pc	180.00
B	25	14/15	4 Pc	190.00
C	50	19/26	4 Pc	200.00
D	100	24/29	4 Pc	250.00
E	150	24/29	4 Pc	280.00
F	250	24/29	4 Pc	290.00
G	250	29/29	4 Pc	350.00
H	500	24/29	2 Pc	430.00
I	500	29/32	2 Pc	490.00
J	1000	24/29	2 Pc	750.00

941 Pear Shaped Flasks, Single Neck

	Cap. ml.	Socket size	Pack per box	Price each pc
A	5	10/19	4 Pc	190.00
B	5	14/23	4 Pc	190.00
C	10	10/19	4 Pc	210.00
D	10	14/23	4 Pc	210.00
E	25	10/19	4 Pc	250.00
F	25	14/23	4 Pc	250.00
G	50	14/23	4 Pc	270.00
H	50	19/26	4 Pc	270.00
I	100	14/23	4 Pc	310.00
J	100	19/26	4 Pc	310.00
K	100	24/29	4 Pc	365.00



946 Filter Flasks with I/C joint

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- comes with I/C joint for fitting

	Cap. ml.	Socket size	Pack per box	Price each pc
A	100	19/26	4 Pc	330.00
B	100	24/29	4 Pc	340.00
C	250	24/29	4 Pc	390.00
D	250	29/32	4 Pc	450.00
E	500	24/29	2 Pc	580.00
F	500	29/32	2 Pc	640.00
G	1000	24/29	2 Pc	875.00
H	1000	29/32	2 Pc	950.00
I	1000	34/35	2 Pc	1020.00






951 Flask Filter-Conical flask with side tube, Heavy Wall

	Cap. ml.	Pack per box	Price each pc
A	100	4 Pc	420.00
B	250	4 Pc	530.00
C	500	2 Pc	700.00
D	1000	2 Pc	980.00

956 Iodine Flasks, Cup Shaped

- For the determination of iodine
- Flask has ground socket and cup shaped top
- Stopper has extended handle which projects well above the liquid seal
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Socket size	Pack per box	Price each pc
A	50	24/29	2 Pc	250.00
B	100	24/29	2 Pc	250.00
C	250	24/29	2 Pc	270.00
D	500	24/29	2 Pc	340.00
E	1000	24/29	1 Pc	800.00




957 Oxygen Flasks with **Platinum** basket

- For determination of Halogens, Sulphur, Phosphorous metal etc.
- Comes with Platinum wire gauge basket for sample
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Pack per box	Price each pc
A	250	1 Pc	on Request
B	500	1 Pc	on Request
C	1000	1 Pc	on Request





960 Pear Shape Flasks for Rotary Evaporators App With I/C Joint

- Especially for **rotary evaporators**, with center neck and conical joint, pear-shaped
- Robust heavy wall construction for added strength
- Shape allows easy removal of solid residues
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Socket size	Pack per box	Price each pc
A	50	29/32	4 Pc	460.00
B	100	29/32	4 Pc	675.00
C	250	29/32	4 Pc	900.00
D	500	29/32	2 Pc	1000.00
E	1000	29/32	2 Pc	1400.00
F	2000	29/32	1 Pc	2350.00

976 Distillation Flasks

- All side arm tubes are sealed at an angle of 75° from the neck

	Cap. ml.	Pack per box	Price each pc
A	10	6 Pc	165.00
B	25	6 Pc	180.00
C	50	6 Pc	200.00
D	100	6 Pc	235.00



981 Distillation Flasks, as per I. P. ASTM Specification


- ASTM standards and test methods for petroleum products
- Uniform wall thickness distribution makes these flasks ideal for heating applications and distillations.
- All sidearm tubes are sealed at an angle of 75° from the neck
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Pack per box	Price each pc
A	100	6 Pc	250.00
B	125	6 Pc	300.00
C	150	6 Pc	375.00
D	200	3 Pc	600.00
E	250	3 Pc	700.00





Funnels



986 Funnel with Thistle Head

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Stem Dia x Length mm.	Pack per box	Price each pc
A	7 x 300	6 Pc	200.00

991 Filter Funnels with Long Stem

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- funnel is constructed with heavy uniform walls, molded rims, and fire polished stems to give a long service life

	Dia mm	Pack per box	Price each pc
A	15	10 Pc	65.00
B	25	10 Pc	70.00
C	35	10 Pc	90.00
D	50	4 Pc	95.00
E	65	4 Pc	110.00
F	75	4 Pc	120.00
G	100	4 Pc	175.00
H	125	2 Pc	400.00



996 Funnels, Short Stem (For Powder)

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- with a heavy uniform wall and a strong, fire-polished rim and stem to provide a long service life
- available with standered joint also

	Dia mm	Pack per box	Price each pc
A	75	4 Pc	200.00
B	100	4 Pc	290.00
C	125	2 Pc	475.00
E	75 B-24	4 Pc	410.00
F	100 B-24	4 Pc	500.00
G	125 B-24	2 Pc	625.00
EE	75 B-19	4 Pc	400.00
H	150 B-24	1 Pc	800.00





997 Analytical funnel. (PP)

- Made of polypropylene.

	Size mm	Pack per box	Price per pc
A	50 mm (inside 50, stem length 50, pp)	72 Pc	15.00
B	62 mm (inside 62, stem length 67, pp)	36 Pc	27.00
C	75 mm (inside 75, stem length 77, pp)	36 Pc	28.50
D	100 mm (inside 100, stem length 95, pp)	36 Pc	57.00

998 Funnel Holder, Polypropylene

- Corrosion free Holder can hold funnels with dia 3' to 6'
- Two tapered wedges ensure firm grip on rods having dia 1/2" or 3/8"

	Type	Pack per box	Price per pc
A	Funnel Holder Single	12 pc	156.00
B	Funnel Holder Double	6 pc	306.00



Refer Cat. No.1911



Refer Cat. No.1906






Filter Pump



1001 Filter Pump

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass


	Description	Price each pc
A	University Pattern	360.00



1006 Filter Tube

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Height mm.	Pack per box	Price each pc
A	12	100	12 Pc	80.00
B	15	125	12 Pc	90.00
C	18	150	12 Pc	95.00
D	25	150	12 Pc	100.00



1011 Digestion Tube


- Digestion tubes designed for uniform wall thickness and impact strength. Made of borosilicate glass with heavy duty retention lip and contoured bottom for optimum heat transfer and minimal thermal shock


	Description	Price each pc
A	24 x 250 mm	250.00

1016 Socket, Plain Shank Double

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Socket size	Shank Length mm	Pack per box	Price each pc
A	10/19	80	25 Pc	80.00
B	14/23	80	25 Pc	85.00
C	19/26	80	25 Pc	90.00
D	24/29	80	25 Pc	105.00
E	29/32	90	10 Pc	170.00
F	34/35	100	6 Pc	260.00





1031 Cone, Plain End Double



- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Socket size	Shank Length mm	Pack per box	Price each pc
A	10/19	80	25 Pc	80.00
B	14/23	80	25 Pc	85.00
C	19/26	80	25 Pc	90.00
D	24/29	80	25 Pc	105.00
E	29/32	90	10 Pc	170.00
F	34/35	100	6 Pc	260.00

1041 Cones, long Tip, Double

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cone size	Shank Length mm	Pack per box	Price each pc
A	14/23	80	25 Pc	150.00
B	19/26	80	25 Pc	160.00
C	24/29	80	25 Pc	175.00
D	29/32	90	10 Pc	230.00
E	34/35	100	6 Pc	360.00



1046 Spherical Joints, Cup


- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Size	Pack per box	Price each pc
A	S 13	25 Pc	130.00
B	S 19	25 Pc	150.00
C	S 29	25 Pc	200.00
D	S 35	25 Pc	250.00

1051 Spherical Joints, Ball.

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Size	Pack per box	Price each pc
A	S 13	25 Pc	130.00
B	S 19	25 Pc	150.00
C	S 29	25 Pc	200.00
D	S 35	25 Pc	250.00





Impinger & Stalgnometer



1056 Standard Impinger for air Sampling as per 'NEERI' Specification

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Pack per box	Price each pc
A	35	1 Pc	380.00
B	50	1 Pc	590.00
C	80	1 Pc	950.00
D	150	1 Pc	1200.00
E	275	1 Pc	1700.00

1061 Stalgnometer

- Tubes specifically designed for checking the surface tension of liquids
- Each tube is accurately calibrated to help quickly determine surface tension by the drop count method.
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

		Pack per box	Price each pc
A	Straight Type with Glass Button	25 Pc	275.00
B	Bend Type with Glass Button	25 Pc	325.00
C	With LDPE Button	25 Pc	260.00
D	With Rotaflow Stopcock	10 Pc	360.00




Syrings Filters



Refer Cat No. 71, 76 & 81



J-SIL
www.j-sil.in



1066 PTFE Bladed Stirrers

- Complete with glass shaft and PTFE Blade
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Description	Pack per box	Price each pc
A	Shaft dia 6 mm Shaft Length 450 mm to Pass through 19/26 Socket Blade width 50 mm Curved to suit Flask 500 ml	10 Pc	190.00
B	Shaft dia 6 mm Shaft Length 450 mm to Pass through 24/29 Socket Blade width 75 mm Curved to suit Flask 1000 ml	10 Pc	200.00
C	Shaft dia 8 mm Shaft Length 600 mm to Pass through 34/35 Socket Blade width 100 mm Curved to suit Flask 2000 ml	10 Pc	295.00
D	Shaft dia 8 mm Shaft Length 600 mm To Pass through 34/35 Socket Blade width 125 mm Curved to suit Flask 5000 ml	10 Pc	375.00
E	Shaft dia 10 mm Shaft Length 450 mm to Pass through 34/35 Socket Blade width 100 mm Curved to suit Flask 5000 ml	5 Pc	400.00
F	Shaft dia 10 mm. Shaft Length 600 mm to Pass through 55/44 Socket Blade width 150 mm Curved to suit Flask 10000 ml	5 Pc	800.00
G	Shaft dia 10 mm. Shaft Length 750 mm to Pass through 55/44 Socket Blade width 200 mm Curved to suit Flask 20000 ml	5 Pc	1400.00





1067 PTFE Bladed Stirrers Glass Rod only

- Glass Shaft without PTFE Blade
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Dia.	Length	Pack per box	Price each pc
G	6 mm	450 mm	10 Pc	115.00
H	8 mm	450 mm	10 Pc	130.00
I	8 mm	600 mm	10 Pc	140.00
J	10 mm	450 mm	10 Pc	160.00
K	10 mm	600 mm	5 Pc	180.00
L	10 mm	750 mm	5 Pc	260.00
M	12 mm	750 mm	5 Pc	330.00
N	50 mm	Teflon Blade	6 mm	80.00
O	75 mm	Teflon Blade	6 mm	100.00
P	100 mm	Teflon Blade	8/10 mm	150.00
Q	125 mm	Teflon Blade	10 mm	300.00
R	150 mm	Teflon Blade	12 mm	650.00
S	200 mm	Teflon Blade	12 mm	1250.00



Refer Cat. No. 920



Refer Cat. No.1388




Refer Cat. No.1389



Refer Cat. No.1492



Refer Cat. No.1494



Refer Cat. No.1495



Refer Cat. No.1497



Refer Cat. No.1499



Stop Cocks



1086 Stop-cock, Straight Bore solid Stopper Plug




- Comes with solid glass plug and screw thread washer, one straight bore
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- comes with both side plain tube

	Bore mm	Price each pc
A	2	170.00
C	4	210.00
D	6	300.00

1091 Stop-cock, Straight Bore with Teflon Plunger

- Comes with solid Teflon plug and screwthread washer,
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Bore mm	Pack per box	Price each pc
A	2	10 Pc	310.00
B	3	10 Pc	350.00
C	4	10 Pc	400.00



1096 Stop-cock, for Burettes


- Comes with solid glass plug or Teflon Plunger and screwthread washer.
- Comes with one side capillary & other end with glass tube
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Bore mm	Plunger	Pack per box	Price each pc
A	2	Glass	25 Pc	160.00
B	2	Teflon	25 Pc	310.00

1101 Stop-cock Teflon Rotaflow for Burette, Sep funnel or column


- Comes with Smooth operating spindle, Made of PTFE
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Bore mm	Use for	Price each pc
A	2	Burette	90.00
B	2	Sep funnel	90.00
C	4	Sep funnel	195.00
D	6	Sep funnel	295.00





Stop Cocks



1106 Stop-cock, Double Oblique Bore Solid Stopper Plug


- Stopcock has a solid glass plug having parallel diagonal bores
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Bore mm	Price each pc
A	2	1000.00
B	3	1100.00

1111 Stop-cock, 'T' Bore Solid Stopper plug Bore

- T-Bore 90° Stopcock with solid Glass Plug
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Bore mm	Pack per box	Price each pc
A	2	10 Pc	1000.00
B	3	10 Pc	1100.00



Volumetric Flask




Refer Cat No. 256 & 257



J-SIL
www.j-sil.in

Stoppers



1140 Stoppers, Solid Penny Head/Flat Head

- Penny-Head Stopper, Glass, Solid, Various Joint Sizes Available.
 - Available in clear & amber colour
- Available for Vol Flask & Reagent Bottle
Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Size New	Plain/ Amber	Suitable for	Pack per box	Price each pc
A	10/19	Plain	Vol Flask	100 Pc	55.00
B	14/23	Plain	Vol Flask	100 Pc	60.00
C	19/26	Plain	Vol Flask	50 Pc	75.00
D	24/29	Plain	Vol Flask	50 Pc	125.00
E	10/19	Amber	Vol Flask	100 Pc	75.00
F	14/23	Amber	Vol Flask	100 Pc	95.00
G	19/26	Amber	Vol Flask	50 Pc	120.00
H	24/29	Amber	Vol Flask	50 Pc	190.00
I	14/23	Plain	Reagent Bottle	50 Pc	95.00
J	19/26	Plain	Reagent Bottle	50 Pc	110.00
K	24/29	Plain	Reagent Bottle	25 Pc	160.00
L	29/32	Plain	Reagent Bottle	10 Pc	210.00
M	34/35	Plain	Reagent Bottle	10 Pc	325.00
N	14/23	Amber	Reagent Bottle	50 Pc	145.00
O	19/26	Amber	Reagent Bottle	50 Pc	160.00
P	24/29	Amber	Reagent Bottle	25 Pc	210.00
Q	29/32	Amber	Reagent Bottle	10 Pc	270.00
R	34/35	Amber	Reagent Bottle	10 Pc	400.00
S	19/20	Plain	B.O.D. Bottle	25 Pc	100.00
AA	12/21	Plain	Vol Flask 50ml	100 Pc	60.00



1141 Stoppers, Hexagonal Hollow

- Flat Bottom Hollow Stoppers, Various Joint Sizes available
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cone size	Pack per box	Price each pc
A	10/19	100 Pc	22.00
B	14/23	100 Pc	25.00
C	19/26	100 Pc	30.00
D	24/29	50 Pc	42.00
E	29/32	50 Pc	80.00
F	34/35	25 Pc	110.00
H	24/29 For Iodine Flask	50 Pc	100.00





1146 Extractor for Soxhlet

- For the extraction of solids with a suitable solvent
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- The vapour tube also protects the siphon tube from accidental damage

	Cap. ml.	Socket size	Cone size	Pack per box	Price each pc
A	60	34/35	24/29	1 Pc	780.00
B	100	40/38	24/29	1 Pc	1060.00
C	200	50/44	24/29	1 Pc	1250.00
D	400	50/44	24/29	1 Pc	3600.00
E	600	50/44	34/35	1 Pc	6000.00
F	1000	50/44	34/35	1 Pc	6400.00

👉 For Complete set, Kindly refer Cat. No. 1706

1151 Condenser only for Soxhlet

- Comes with I/C cone joint
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Inner bulb sections provide increased cooling surface

	Length mm	Cone size	Pack per box	Price each pc
A	150	34/35	2 Pc	680.00
B	200	40/38	2 Pc	850.00
C	250	50/44	1 Pc	1150.00
D	300	50/44	1 Pc	2000.00



1156 Extraction Flasks

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- uniform wall thickness & flat shape for ideal heating on hot plate
- comes with Interchangeable joint

	Cap. ml.	Socket size	Pack per box	Price each pc
A	100	24/29	4 Pc	175.00
B	150	24/29	4 Pc	220.00
C	250	24/29	4 Pc	220.00
D	500	24/29	2 Pc	325.00
E	1000	24/29	2 Pc	480.00
F	2000	34/35	1 Pc	1250.00
G	3000	34/35	1 Pc	2400.00





1161 Reagent Bottle, Narrow mouth with I/C Stopper Flat Head,

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Comes with solid glass flat head stopper

	Cap. ml.	Joint size	Approx O.D. x Height mm	Pack per box	Price each pc
A	30	14/19	36x82	6 Pc	250.00
B	60	14/23	45x97	6 Pc	275.00
C	125	19/26	55x120	6 Pc	310.00
D	250	19/26	70x140	6 Pc	365.00
E	500	24/29	86x177	2 Pc	415.00
F	1000	29/32	101x220	1 Pc	630.00
G	2000	34/35	136x260	1 Pc	1180.00


1162 Reagent Bottle Amber Colour Narrow mouth with I/C Stopper Flat Head Graduates

- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
- Comes with solid glass flat head stopper
- Comes with Amber colour for light sensitive media.

	Cap. ml.	Joint size	Approx O.D. x Height mm	Pack per box	Price each pc
A	30	14/19	36x82	6 Pc	440.00
B	60	14/23	45x97	6 Pc	465.00
C	125	19/26	55x120	6 Pc	500.00
D	250	19/26	70x140	6 Pc	610.00
E	500	24/29	86x177	2 Pc	725.00
F	1000	29/32	101x220	1 Pc	1000.00
G	2000	34/35	136x260	1 Pc	1850.00



Micropipette



Refer Cat No. JVR





1163 Bottles,B.O.D.,with interchangeable stopper

- Comes with solid glass penny head stopper
- Individual number printed on the bottle
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass

	Cap. ml.	Joint size	Approx O.D. x Height mm	Pack per box	Price each pc
A	125	19/20	55x120	2 Pc	425.00
B	300	19/20	70x143	2 Pc	550.00
C	60	19/20	42x100	2Pc	380.00



1164 Bottles,B.O.D.,with I/C Stopper **Amber Colour**

- Comes with solid glass penny head stopper
- Individual number printed on the bottle
- Made from Heat Resistant, Low Expansion 3.3 Borosilicate Glass
Comes with uniform Amber colour

	Cap. ml.	Joint size	Approx O.D. x Height mm	Pack per box	Price each pc
A	125	19/20	55x120	2 Pc	1000.00
B	300	19/20	70x143	2 Pc	1150.00
C	60	19/20	42x100	2Pc	900.00

Vials



Refer Cat No. 61 & 66



J-SIL
www.j-sil.in



1166 Reagent Bottle, Graduated with P.P. Screw Cap (Blue) Wide mouth and Pouring ring

- Autoclavable, Made from 3.3 borosilicate glass
- Permanent amber graduations and marking spots
- Approx. 30mm I.D. **GL 45 screw-thread opening**
- Linerless blue polypropylene cap autoclavable to 140°C and replaceable, drip-free, clear polypropylene pour ring included

	Cap. ml.	Neck Specification	Approx O.D. x Height mm	Pack per box	Price each pc
A	10	GL 25	36x50	10 Pc	350.00
B	25	GL 25	36x70	10 Pc	260.00
C	50	GL 45	46x87	6 Pc	310.00
D	100	GL 45	56x100	6 Pc	330.00
E	150	GL 45	60x120	6 Pc	345.00
F	250	GL 45	70x138	6 Pc	350.00
G	500	GL 45	86x176	2 Pc	370.00
H	1000	GL 45	101x225	1 Pc	590.00
I	2000	GL 45	136x260	1 Pc	1040.00
J	5000	GL 45	183x323	1 Pc	5400.00
K	3000	GL 45	158x295	1 Pc	4700.00



1167 Reagent Bottle Amber Colour Graduated with P.P. Screw Cap (Blue) Wide mouth and Pouring ring

- Autoclavable, Made from 3.3 borosilicate amber glass
- White graduations and marking spots
- Approx. 30mm I.D. **GL 45 screw-thread opening**
- Linerless blue polypropylene cap autoclavable to 140°C and replaceable, drip-free, clear polypropylene pour ring included

	Cap. ml.	Neck Specification	Approx O.D. x Height mm	Pack per box	Price each pc
A	10	GL 25	36x50	10 Pc	450.00
B	25	GL 25	36x70	10 Pc	410.00
C	50	GL 45	46x87	6 Pc	450.00
D	100	GL 45	56x100	6 Pc	500.00
E	250	GL 45	70x138	6 Pc	540.00
F	500	GL 45	86x176	2 Pc	650.00
G	1000	GL 45	101x225	1 Pc	1100.00



1169 Screw Cap For Bottle (with Pouring Ring)

- Made from polypropylene.
- Autoclavable up to 140°C.



	Colour	Size	Pack per box	Price each pc
A	Blue	GL 45	10 pc	80.00
B	Yellow	GL 45	10 pc	90.00
C	Green	GL 45	10 pc	90.00
D	Red	GL 45	10 pc	90.00
E	Orange	GL 45	10 pc	90.00
F	Grey	GL 45	10 pc	90.00

1171 Reagent Bottle wide mouth with I/C Hollow Stopper

- Autoclavable, Made from 3.3 borosilicate glass
- comes with wide mouth & Hollow stopper

	Cap. ml.	Joint size	Approx O.D. x Height mm	Pack per box	Price each pc
A	30	19/26	36x82	6 Pc	300.00
B	60	24/29	45x105	6 Pc	350.00
C	125	29/32	56x145	6 Pc	395.00
D	250	34/35	70x150	6 Pc	560.00



Refer Cat. No. 1891



Refer Cat. No. 1892





Tall



Squat



1181 Weighing Bottle with I/C Stopper

- Autoclavable, Made from 3.3 borosilicate glass


	Cap. ml.	Dia x Ht. mm	Shape	Pack per box	Price each pc
A	5	20 x 40	Tall	20 Pc	250.00
B	10	40 x 20	Tal	10 Pc	350.00
D	15	25 x 60	Tall	10 Pc	275.00
E	15	40 x 30	Squat	10 Pc	380.00
H	25	40 x 40	Squat	10 Pc	500.00
J	25	30 x 65	Tall	10 Pc	380.00
K	20	50 x 35	Squat	10 Pc	485.00
M	40	60 x 40	Squat	6 Pc	725.00
N	60	40 x 90	Tall	10 Pc	650.00
O	80	70 x 40	Squat	10 Pc	1050.00
P	20	50 x 35 *	Squat	10 Pc	800.00
Q	40	60 x 40 *	Squat	6 Pc	990.00

* Amber Colour

1182 Glass Weighing Boat/Scoop

- Autoclavable, Made from 3.3 borosilicate glass for weighing small quantities of liquid or powder

	Cap. ml.	Pack per box	Price each pc
A	0.5	50 Pc	150.00
B	1	50 Pc	150.00
C	2	50 Pc	150.00
D	3	50 Pc	160.00
E	6	50 Pc	180.00
F	10	25 Pc	190.00



1196 Liquid Weighing Bottle with Rubber Teat

- Autoclavable, Made from 3.3 borosilicate glass

	Cap. ml.	Joint size	Pack per box	Price each pc
A	25	19/17	2 Pc	750.00
B	50	19/17	2 Pc	850.00





1201 Dropping Bottle with Pipette and Rubber Teat

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with interchangeable ground joint transparent
- dropper with rubber teat

	Cap. ml.	Joint size	Pack per box	Price each pc
A	30	14/19	6 Pc	275.00
B	60	14/19	6 Pc	335.00
C	125	19/26	6 Pc	365.00
D	250	19/26	6 Pc	550.00
E	Dropper	B-14	25 Pc	90.00
F	Dropper	B-19	25 Pc	110.00

1202 Dropping Bottle Complete Amber Colour with Pipette and Rubber Teat

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with interchangeable ground joint amber
- Comes with uniform Amber colour

	Cap. ml.	Joint size	Pack per box	Price each pc
A	30	14/19	6 Pc	430.00
B	60	14/19	6 Pc	540.00
C	125	19/26	6 Pc	595.00
D	250	19/26	6 Pc	765.00



Refer Cat. No. 71



Refer Cat. No. 76



Refer Cat. No. 81



Wash Bottles



1205 Bottles, Wash, PE

- Made of low Density **polyethylene**.
- These Bottles contain built in delivery tube making it as one piece wash Bottle & it is translucent & unbreakable.
- These Bottles are Graduated on the outer surface to give an approximate volume of solution filled or dispersed.
- The cap is fitted with a flexible polyethylene delivery tube which can be aimed wherever required.

	Capacity ml	Pack per box	Price per pc
A	125	12 Pc	72.00
B	250	12 Pc	90.00
C	500	6 Pc	120.00
D	1000	6 Pc	210.00

1206 Wash Bottle Complete

- Autoclavable, Made from 3.3 borosilicate glass

	Cap. ml.	Joint size	Pack per box	Price each pc
A	250	24/29	2 Pc	600.00
B	500	24/29	2 Pc	750.00
C	1000	24/29	2 Pc	950.00



1211 Wash Bottle Fitting only

- Autoclavable, Made from 3.3 borosilicate glass

	Cap. ml.	Pack per box	Price each pc
A	250	6 Pc	450.00
B	500	6 Pc	470.00
C	1000	6 Pc	490.00





Wash Bottles



1216 Gas Wash Bottle with Interchangeable fitting

- Autoclavable, Made from 3.3 borosilicate glass

	Cap. ml.	Joint size	Pack per box	Price each pc
A	125	29/32	1 Pc	900.00
B	250	29/32	1 Pc	1000.00
C	500	29/32	1 Pc	1300.00

1221 Gas Wash Bottle only

- Autoclavable, Made from 3.3 borosilicate glass

	Cap. ml.	Joint size	Pack per box	Price each pc
A	125	29/32	1 Pc	700.00
B	250	29/32	1 Pc	750.00
C	500	29/32	1 Pc	990.00



1226 Gas Wash Bottle, fitting only

- Autoclavable, Made from 3.3 borosilicate glass

	Cap. ml.	Joint size	Pack per box	Price each pc
A	125	29/32	6 Pc	450.00
B	250	29/32	6 Pc	460.00
C	500	29/32	6 Pc	490.00

1231 Gas Wash Bottle, with sintered end for downward gas distribution, complete

- Autoclavable, Made from 3.3 borosilicate glass
- comes with sintered dias at bottom in head

	Cap. ml.	Joint size	Pack per box	Price each pc
A	125	29/32	1 Pc	1800.00
B	250	29/32	1 Pc	2000.00
C	500	29/32	1 Pc	2300.00



Separating Funnels



1251 Separating Funnels with Glass Stop Cock, Pear Shape

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with standered glass stop cock & solid glass stoper

	Cap. ml.	Stopper size	Pack per box	Price each pc
A	25	14/19	1 Pc	600.00
B	50/60	14/19	1 Pc	700.00
C	100/125	19/26	1 Pc	775.00
D	250	24/29	1 Pc	910.00
E	500	24/29	1 Pc	1150.00
F	1000	29/32	1 Pc	1900.00
G	2000	29/32	1 Pc	3500.00

1256 Separating Funnels, with Teflon Rotaflow Plug, Pear Shape

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with boroflow stop cock with teflon key & solid glass stopper

	Cap. ml.	Stopper size	Pack per box	Price each pc
A	25	14/19	1 Pc	550.00
B	50/60	14/19	1 Pc	600.00
C	100/125	19/26	1 Pc	700.00
D	250	24/29	1 Pc	800.00
E	500	24/29	1 Pc	1050.00
F	1000	29/32	1 Pc	1700.00
G	2000	29/32	1 Pc	3400.00



1257 Separating Funnels Amber Colour with Teflon Rotaflow/ Stop-cock

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with boroflow stop cock with teflon key & solid glass stopper
- With uniform Amber colour

	Cap. ml.	Stopper size	Pack per box	Price each pc
A	25	14/19	1 Pc	1450.00
B	50/60	14/19	1 Pc	1600.00
C	100/125	19/26	1 Pc	1800.00
D	250	24/29	1 Pc	1900.00
E	500	24/29	1 Pc	2500.00
F	1000	29/32	1 Pc	4000.00





Separating Funnels



1261 Separating Funnels with PTFE Plunger Stop-cock and I/C stopper

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with PTFE/ Teflon key stop cock & solid glass stopper

	Cap. ml.	Stopper size	Pack per box	Price each pc
A	100/125	19/26	1 Pc	900.00
B	250	24/29	1 Pc	950.00
C	500	24/29	1 Pc	1300.00
D	1000	29/32	1 Pc	2000.00
E	2000	29/32	1 Pc	3600.00

1266 Separating Funnels with Glass Stop Cock, Globe Shape

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with standered glass stop cock & solid glass stopper

	Cap. ml.	Pack per box	Price each pc
A	25	2 Pc	600.00
B	50	1 Pc	700.00
C	100	1 Pc	775.00
D	250	1 Pc	910.00
E	500	1 Pc	1150.00
F	1000	1 Pc	1900.00
G	2000	1 Pc	3500.00
H	3000	1 Pc	5800.00
I	5000	1 Pc	7600.00



1271 Separating Funnels, with Teflon Rotaflow Plug, Globe Shape


- Autoclavable, Made from 3.3 borosilicate glass
- Comes with boroflow stop cock with teflon key & solid glass stopper

	Cap. ml.	Pack per box	Price each pc
A	25	2 Pc	550.00
B	50	1 Pc	600.00
C	100	1 Pc	700.00
D	250	1 Pc	800.00
E	500	1 Pc	1050.00
F	1000	1 Pc	1700.00
G	2000	1 Pc	3400.00
H	3000	1 Pc	5700.00
I	5000	1 Pc	7500.00



J-SIL
www.j-sil.in

Separating Funnels



1272 Separating Funnel Holder, Polypropylene

- Non-corrosive front opening provides clear view of the solution, Allows unobstructed placement of separating Funnels
- Two tapered wedges ensures firm grip on rods having dia ½ or 3/5

	Capacity ml.	Pack per box	Price per pc
A	Separating Funnel holder	12 pc	190.00

1273 Separating Funnel, Polypropylene

- Strenghth
- Pear per shaped
- Has leakproof screw cap threaded
- PP stopcock fitted with self lubricating PTEE plug, is push fitted to the stem of funnels


	Capacity ml.	Pack per box	Price per pc
A	100	2 pc	600.00
B	250	2 pc	690.00
C	500	2 pc	840.00



1282 Separating Funnels, with Teflon plug, Cylindrical Graduated

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with boroflow stop cock with teflon key & solid glass stopper
- comes with Mohr type graduation in amber colour

	Cap. ml.	Pack per box	Price each pc
A	100	1 Pc	800.00
B	250	1 Pc	1000.00
C	500	1 Pc	1250.00



Refer Cat. No. JVR on Page on 8





Separating Funnels



1291 Separating Funnels, Cylindrical with socket & cone with glass stop cock


- Autoclavable, Made from 3.3 borosilicate glass
- Comes with standered glass stop cock
- comes with I/C joints, cone at bottam & socket at top

	Cap. ml.	Socket size	Cone size	Pack per box	Price each pc
A	20	10/19	10/19	1 Pc	840.00
B	25	14/23	14/23	1 Pc	880.00
C	50	19/26	19/26	1 Pc	880.00
E	100	19/26	19/26	1 Pc	990.00
G	250	24/29	24/29	1 Pc	1150.00
I	500	24/29	24/29	1 Pc	1600.00
K	1000	24/29	24/29	1 Pc	2800.00

1296 Separating Funnels, Open top

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with standered glass stop cock
comes with I/C cone at bottom with stem

	Cap. ml.	Cone size	Pack per box	Price each pc
A	50	19/26	2 Pc	600.00
B	100	19/26	2 Pc	625.00
C	100	24/29	2 Pc	750.00
D	250	19/26	1 Pc	1000.00



Refer Cat. No. 1388



Refer Cat. No. 1389



J-SIL
www.j-sil.in

Separating Funnels



1301 Separating Funnels Cylindrical with Socket, Cone & Teflon Plunger Stop Cock

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with PTFE/ Teflon key stop cock
- comes with I/C joints, cone at bottom & socket at top

	Cap. ml.	Socket size	Cone size	Pack per box	Price each pc
A	100	24/29	24/29	1 Pc	1020.00
B	250	24/29	24/29	1 Pc	1150.00
C	500	24/29	24/29	1 Pc	1550.00

1306 Pressure Equalising Funnels, Cylindrical with Teflon Plunger Stop Cock

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with PTFE/ Teflon key stop cock & side
- comes with I/C joints, cone at bottom & socket at top

	Cap. ml.	Socket size	Cone size	Pack per box	Price each pc
A	50	19/26	19/26	1 Pc	1350.00
B	100	24/29	24/29	1 Pc	1400.00
C	250	24/29	24/29	1 Pc	1550.00
D	500	24/29	24/29	1 Pc	1950.00



1311 Pressure Equalising Funnels, Cylindrical with Glass Stop Cock

- Autoclavable, Made from 3.3 borosilicate glass
- Comes with standard glass stop cock
- comes with I/C joints, cone at bottom & socket at top

	Cap. ml.	Socket size	Cone size	Pack per box	Price each pc
A	100	24/29	24/29	1 Pc	1350.00
B	250	24/29	24/29	1 Pc	1500.00
D	500	24/29	24/29	1 Pc	1900.00





Tilt Measure



1316 Test Jar (Pour Point Tube)


- Autoclavable, Made from 3.3 borosilicate glass

	Flat bottom 30 to 34 mm inside dia 115 to 125 mm height marked with line 54 ± 3 mm	Pack per box	Price each pc
	above in side bottom	12 Pc	265.00

1321 Automatic Tilt Measures

- Autoclavable, Made from 3.3 borosilicate glass
- For accurate dosing of acid and alcohol.
- comes with the conical flask

	Cap. ml.	Flask Cap. ml.	Pack per box	Price each pc
A	1	100	1 Pc	800.00
B	5	100	1 Pc	800.00
C	10	250	1 Pc	850.00
D	20	250	1 Pc	950.00
E	25	250	1 Pc	1060.00
F	50	500	1 Pc	1300.00



Syringe filters



Refer Cat. No. 71, 76, 81



J-SIL
www.j-sil.in



1326 Petri Dishes Plain

- Autoclavable, Made from 3.3 borosilicate glass
- withstand repeated sterilization (wet or dry)
- The edges are beaded to provide greater mechanical strength

	O.D. mm.	Height mm.	Pack per box	Price each pc
A	50	17	5 Pc	175.00
B	80	17	10 Pc	200.00
C	100	17	10 Pc	225.00
D	150	20	5 Pc	650.00
E	90	15	10 Pc	180.00
F	100	20	10 Pc	240.00

1327 Dishes, Evaporating, Flat Bottom, with Pour Out

- Autoclavable, Made from 3.3 borosilicate glass
- withstand repeated sterilization (wet or dry)
- The edges are beaded to provide greater mechanical strength

	O.D. mm.	Height mm.	Pack per box	Price each pc
E	80	45	10 Pc	345.00
F	105	55	5 Pc	500.00
G	150	80	1 Pc	925.00



1328 Dishes, Crystallizing

- Autoclavable, Made from 3.3 borosilicate glass
- withstand repeated sterilization (wet or dry)
- The edges are beaded to provide greater mechanical strength

	O.D. mm.	Height mm.	Pack per box	Price each pc
A	100	50	5 Pc	425.00
B	150	75	1 Pc	670.00

Micropipette



Refer Cat. No. JVR





1329 Dish, Evaporating, Euro Design

- Made of Porcelain
- Goes upto 1000°C

	Cap. mm	OD x H mm	Pack per box	Price each pc
A	35	50 X 24	6 PC	55.00
B	70	75 X 30	6 PC	70.00
E	150	100 X 42	6 PC	110.00

1330 Petri Dish (Disposable) Gama Sterile

- These transparent, ready to use.
- These are made of optically clear, non-toxic, polystyrene.
- These are packed individually in polyethylene bags.

	Size mm	Pack per box	Price per pack
A	90	460 pcs.	9900.00



Refer Cat. No. 1891



Refer Cat. No. 1892



Tubes & Connectors



1336 Melting Point Tubes, Thiel Tube

- Autoclavable, Made from 3.3 borosilicate glass
- are used for melting point determination tests

	Length x Dia mm.	Side tube Dia mm	Pack per box	Price each pc
A	150 x 18	10	10 Pc	135.00
B	150 x 25	10	10 Pc	155.00

1346 Calcium Chloride Tubes Straight with one Bulb

- Autoclavable, Made from 3.3 borosilicate glass

	Length x Dia mm.	Pack per box	Price each pc
A	100 x 12	10 Pc	95.00
B	125 x 15	10 Pc	110.00
C	150 x 18	10 Pc	130.00



1356 Calcium Chloride Towers

- Autoclavable, Made from 3.3 borosilicate glass

	Length x Dia mm.	Joint size	Pack per box	Price each pc
A	250 x 40	34/35	1 Pc	1850.00
B	300 x 40	34/35	1 Pc	2000.00
C	300 x 60	34/35	1 Pc	2300.00

1361 Connection Tubes, T-Shape

- Autoclavable, Made from 3.3 borosilicate glass

	Length x Dia mm.	Pack per box	Price each pc
A	50 x 8	12 Pc	95.00
B	60 x 10	12 Pc	105.00
C	80 x 12	12 Pc	160.00



1366 Connection Tubes, Y-Shape

- Autoclavable, Made from 3.3 borosilicate glass

	Length x Dia mm.	Pack per box	Price each pc
A	50 x 8	12 Pc	95.00
B	60 x 10	12 Pc	105.00
C	80 x 12	12 Pc	160.00





Gooch Crucible

1370 Tong, Crucible

- Made of SS, Grade 202




	Size mm	Pack per box	Price each pc
A	150	4 pc	125.00
B	200	4 pc	140.00
C	250	4 pc	250.00
D	300	4 pc	300.00
E	375	4 pc	495.00
F	450	4 pc	600.00
G	600	4 pc	720.00

👉 For Beaker & Flask Tong, Kindly refer Cat. No. 889

1371 Gooch Crucible with Sintered Disc

- Autoclavable, Made from 3.3 borosilicate glass
- comes with acid resistant sintered glass disc for filtration
- an easy and fast filtering without the need of any kind of filter paper

	Cap. ml.	Disc. Dia mm.	Porosity Grade	Pack per box	Price each pc
A	15	20	1	25 Pc	180.00
B	15	20	2	25 Pc	180.00
C	15	20	3	25 Pc	180.00
D	15	20	4	25 Pc	180.00
E	30	30	1	25 Pc	250.00
F	30	30	2	25 Pc	250.00
G	30	30	3	25 Pc	250.00
H	30	30	4	25 Pc	250.00
I	50	40	1	25 Pc	360.00
J	50	40	2	25 Pc	360.00
K	50	40	3	25 Pc	360.00
L	50	40	4	25 Pc	360.00



1376 Adapter for Gooch Crucible

- Autoclavable, Made from 3.3 borosilicate glass
- use with crucibles with fritted discs. Also for use with Gooch crucibles or directly with porcelain plates

	To Fit	Pack per box	Price each pc
A	15 ml Crucible	2 Pc	200.00
B	30 ml Crucible	2 Pc	250.00
C	50 ml Crucible	2 Pc	400.00



J-SIL
www.j-sil.in

Buckner Funnel



1380 Buckner funnel, with cone on stem.

- Autoclavable made from 3.3 borosilicate glass.
- With an integral ground glass cone joint on the funnel stem for connecting to filter flask with standard socket joint.
- Useful where conventional filter paper would be destroyed by acid and alkalis.

	Cap. ml.	Disc. Dia mm.	Joint size	Porosity Grade	Pack per box	Price each pc
A	35	30	19/26	1	2	750.00
B	35	30	19/26	2	2	750.00
C	35	30	19/26	3	2	750.00
D	35	30	19/26	4	2	750.00
E	80	40	24/29	1	2	1200.00
F	80	40	24/29	2	2	1200.00
G	80	40	24/29	3	2	1200.00
H	80	40	24/29	4	2	1200.00
I	200	65	24/29	1	2	1400.00
J	200	65	24/29	2	2	1400.00
K	200	65	24/29	3	2	1400.00
L	200	65	24/29	4	2	1400.00
M	500	90	24/29	1	2	3000.00
N	500	90	24/29	2	2	3000.00
O	500	90	24/29	3	2	3000.00
P	500	90	24/29	4	2	3000.00



Refer Cat. No. 1391



Refer Cat. No. 1911



Refer Cat. No. 1906





Silica & Quartz



1387 Crucible Silica Lid only


- Only for Crucible Lid
- Temperature upto 1000°C

	Cap. ml.	Type	Pack per box	Price each pc
A	15	Lid only	6 pc	330.00
B	25	Lid only	6 pc	370.00
C	50	Lid only	6 pc	550.00
D	80	Lid only	4 pc	920.00
E	100	Lid only	4 pc	1060.00
F	150	Lid only	4 pc	1350.00

1388 Crucible Silica

- Available without Lid
- Temperature upto 1000°C

	Cap. ml.	Dia mm.	Height mm.	Pack per box	Price each pc
A	15	41	25	6 pc	330.00
B	25	47	28	6 pc	350.00
C	50	57	37	6 pc	650.00
D	80	67	45	4 pc	1100.00
E	100	71	50	4 pc	1380.00
F	150	82	51	4 pc	1950.00



1389 Quartz Crucible Only


- Made of Quartz Gel & Powder
- Temperature upto 2000°C

	Size mm	Type	Pack per box	Price each pc
A	15	Without Lid	4 pc	770.00
B	25	Without Lid	4 pc	1040.00
C	50	Without Lid	2 pc	1320.00
D	80	Without Lid	2 pc	2300.00
E	100	Without Lid	2 pc	3500.00

1390 Quartz Crucible Lid only

- Made of Quartz Gel & Powder
- Temperature upto 2000°C

	Size mm	Type	Pack per box	Price each pc
A	15	Lid only	4 pc	1265.00
B	25	Lid only	4 pc	1500.00
C	50	Lid only	2 pc	2100.00
D	80	Lid only	2 pc	2900.00
E	100	Lid only	2 pc	3795.00





1391 Buckner Filter with **Sintered Dish** and Plain Stem

- Autoclavable, Made from 3.3 borosilicate glass
- comes with acid resistant sintered glass disc for filtration
- an easy and fast filtering without the need of any kind of filter paper are available in various disc diameters and porosities.

	Cap. ml.	Disc. Dia mm.	Porosity Grade	Pack per box	Price each pc
A	35	30	1	2 Pc	400.00
B	35	30	2	2 Pc	400.00
C	35	30	3	2 Pc	400.00
D	35	30	4	2 Pc	400.00
E	80	40	1	2 Pc	800.00
F	80	40	2	2 Pc	800.00
G	80	40	3	2 Pc	800.00
H	80	40	4	2 Pc	800.00
I	200	65	1	2 Pc	950.00
J	200	65	2	2 Pc	950.00
K	200	65	3	2 Pc	950.00
L	200	65	4	2 Pc	950.00
M	500	90	1	2 Pc	2200.00
N	500	90	2	2 Pc	2200.00
O	500	90	3	2 Pc	2200.00
P	500	90	4	2 Pc	2200.00
Q	1000	100	1	1 Pc	5500.00
R	1000	100	2	1 Pc	5500.00
S	1000	100	3	1 Pc	5500.00
T	1000	100	4	1 Pc	5500.00



Refer Cat. No. 1380



Refer Cat. No. 1891M



Refer Cat. No. 1892J





Buckner Funnel



1393 Funnel, Buchner, **Porcelain**

- Made of porcelain
- Temperature upto 1000°C

	Dia (mm)	Pack per box	Price each pc
A	35	10 pc	150.00
B	50	10 pc	200.00
C	75	6 pc	420.00
D	100	4 pc	540.00
E	125	2 pc	660.00
F	150	2 pc	960.00
G	200	1 pc	2880.00
H	250	1 pc	3960.00
I	300	1 pc	7140.00

1394 Funnel, Buchner, **Polypropylene**

- Light weight.
- Easy to clean and autoclavable.
- Top portion snap fitted with the lower one forms a vacuum seal which requires no greasing.
- Base of the top portion has perforated filter plate.

	Size mm	Pack per box	Price per pc
A	70	12 pcs.	240.00
B	90	12 pcs.	315.00
C	110	6 pcs.	360.00



Refer Cat. No. 1891



Refer Cat. No. 1892



J-SIL
www.j-sil.in



1421 Viscometers College Pattern, Ostwald

- Autoclavable, Made from 3.3 borosilicate glass
- useful for student & training purpose

	Empties in	Pack per box	Price each pc
A	60-90 Sec.	1 Pc	650.00


1426 Viscometers B.S.S.

- Autoclavable, Made from 3.3 borosilicate glass

	B.S.S. No.	Pack per box	Price each pc
A	1	1 Pc	4000.00
B	2	1 Pc	4000.00
C	3	1 Pc	4000.00
D	4	1 Pc	4000.00



Viscometer Holder



Price on Request





1431 Viscometers, U - tube, BS / U (w/o Constant Certificate)

- Autoclavable, Made from 3.3 borosilicate glass

	Size	Viscosity Range mm ² /sec (cSt)*	Pack per box	Price each pc
A	A	0.9-3	1 Pc	4500.00
B	B	2-10	1 Pc	4500.00
C	C	6-30	1 Pc	4500.00
D	D	20-100	1 Pc	4500.00
E	E	60-300	1 Pc	4500.00
F	F	200-1000	1 Pc	4500.00
G	G	600-3000	1 Pc	4500.00
H	H	2000-10000	1 Pc	4500.00



1432 Viscometers BS/U with **Constant Certificate** & Calibrated @40°C

- Autoclavable, Made from 3.3 borosilicate glass
- For measurement of kinematic viscosity of transparent Newtonian liquids according to ASTM D 445 and ISO 3104. Specifications conform to ASTM D 446 and ISO 3105.
- Minimum sample volume 7 mL in sizes A to C, 12 mL in sizes D to F, and 23 mL in sizes G and H.
- comes with certificate of constant

	Cap. ml.	Viscosity Range mm ² /sec (cSt)*	Approx. Constant mm ² /s ² (cSt/s)	Pack per box	Price each pc
A	A	0.9-3	0.003	1 Pc	11500.00
B	B	2-10	0.01	1 Pc	11500.00
C	C	6-30	0.03	1 Pc	11500.00
D	D	20-100	0.1	1 Pc	11500.00
E	E	60-300	0.3	1 Pc	11500.00
F	F	200-1000	1	1 Pc	11500.00
G	G	600-3000	3	1 Pc	11500.00
H	H	2000-10000	10	1 Pc	11500.00

☞ Constant on different Temp. are available against special order on additional cost.





1436 Viscometers, Suspended Level BS / IP / SL (w/o Constant Certificate)

- Autoclavable, Made from 3.3 borosilicate glass

	Size	Viscosity Range mm ² /sec (cSt)*	Pack per box	Price each pc
A	1	2-10	1 Pc	5500.00
B	1 A	6-30	1 Pc	5500.00
C	2	20-100	1 Pc	5500.00
D	2 A	60-300	1 Pc	5500.00
E	3	200-1000	1 Pc	5500.00
F	3 A	600-3000	1 Pc	5500.00
G	4	2000-10000	1 Pc	5500.00
H	4 A	6000-30000	1 Pc	5500.00
I	5	20000-100000	1 Pc	5500.00

1437 Viscometers Suspended Level with **Constant Certificate** & Calibrated @ 40°C

- Autoclavable, Made from 3.3 borosilicate glass
- For measurement of kinematic viscosity of transparent Newtonian liquids according to ASTM D 445 and ISO 3104. Specifications conform to ASTM D 446 and ISO 3105
- Minimum sample volume 11 mL.
- comes with certificate of constant

	Size	Viscosity Range mm ² /sec (cSt)*	Approx. Constant mm ² /s ² (cSt/s)	Pack per box	Price each pc
A	1	2-10	0.01	1 Pc	12250.00
B	1 A	6-30	0.03	1 Pc	12250.00
C	2	20-100	0.1	1 Pc	12250.00
D	2 A	60-300	0.3	1 Pc	12250.00
E	3	200-1000	1.00	1 Pc	12250.00
F	3 A	600-3000	3.0	1 Pc	12250.00
G	4	2000-10000	10	1 Pc	12250.00
H	4 A	6000-30000	30	1 Pc	12250.00
I	5	20000-100000	100	1 Pc	12250.00

☞ Constant on different Temp. are available against special order on additional cost.





Viscometers



1446 Viscometers, Cannon-Fenske type, For Transparent Liquids, BS/ IP/CF (w/o Constant Certificate)


- Autoclavable, Made from 3.3 borosilicate glass

	Size	Viscosity Range mm ² /sec (cSt)*	Pack per box	Price each pc
A	25	0.5-2	1 Pc	4500.00
B	50	0.8-4	1 Pc	4500.00
C	100	3-15	1 Pc	4500.00
D	150	7-35	1 Pc	4500.00
E	200	20-100	1 Pc	4500.00
F	300	50-250	1 Pc	4500.00
G	350	100-500	1 Pc	4500.00
H	400	240-1200	1 Pc	4500.00
I	450	500-2500	1 Pc	4500.00
J	500	1600-8000	1 Pc	4500.00
K	600	4000-20000	1 Pc	4500.00

1451 Viscometers Cannon-Fenske U-Tube for Transparent Liquid with **Constant Certificate** & Calibrated @ 40°C

- Autoclavable, Made from 3.3 borosilicate glass
- For measuring kinematic viscosity of transparent Newtonian liquids, particularly petroleum products or lubricants, according to ASTM D 445 and ISO 3104. Specifications conform to ASTM D 446 and ISO 3105
- Minimum sample volume 7 mL
- comes with certificate of constant

	Size	Viscosity Range mm ² /sec (cSt)*	Approx. Constant mm ² /s ² (cSt/s)	Pack per box	Price each pc
A	25	0.5-2	0.002	1 Pc	11500.00
B	50	0.8-4	0.004	1 Pc	11500.00
C	100	3-15	0.015	1 Pc	11500.00
D	150	7-35	0.035	1 Pc	11500.00
E	200	20-100	0.1	1 Pc	11500.00
F	300	50-250	0.25	1 Pc	11500.00
G	350	100-500	0.5	1 Pc	11500.00
H	400	240-1200	1.2	1 Pc	11500.00
I	450	500-2500	2.5	1 Pc	11500.00
J	500	1600-8000	8	1 Pc	11500.00
K	600	4000-20000	20	1 Pc	11500.00
L	75	1.6-8	0.008	1 Pc	11500.00



☞ Constant on different Temp. are available against special order on additional cost.





1456 Viscometers, Cannon-Fenske, for Opaque Liquids, ASTM-D 445 (w/o Constant Certificate)

- Autoclavable, Made from 3.3 borosilicate glass


	Size	Viscosity Range mm ² /sec (cSt)*	Pack per box	Price each pc
A	100	3-15	1 Pc	4500.00
B	150	7-35	1 Pc	4500.00
C	200	20-100	1 Pc	4500.00
D	300	50-250	1 Pc	4500.00
E	350	100-500	1 Pc	4500.00
F	400	240-1200	1 Pc	4500.00
G	450	500-2500	1 Pc	4500.00
H	500	1600-8000	1 Pc	4500.00
I	600	4000-20000	1 Pc	4500.00

1457 Viscometers C/F for Opaque Liquid with **Constant Certificate** & Calibrated @ 40°C

- Autoclavable, Made from 3.3 borosilicate glass
- Opaque Viscometer measures kinematic viscosity of dark Newtonian liquids according to ASTM D445, ISO 3104 and ASTM D2170. Minimum sample volume is 12 mL
- Minimum sample volume 12 mL
- comes with certificate of constant

	Size	Viscosity Range mm ² /sec (cSt)*	Approx. Constant mm ² /s ² (cSt/s)	Pack per box	Price each pc
A	100	3-15	0.015	1 Pc	12000.00
B	150	7-35	0.035	1 Pc	12000.00
C	200	20-100	0.1	1 Pc	12000.00
D	300	50-250	0.25	1 Pc	12000.00
E	350	100-500	0.5	1 Pc	12000.00
F	400	240-1200	1.2	1 Pc	12000.00
G	450	500-2500	2.5	1 Pc	12000.00
H	500	1600-8000	8	1 Pc	12000.00
I	600	4000-20000	20	1 Pc	12000.00

☞ Constant on different Temp. are available against special order on additional cost.





1461 Viscometers, U-tube, Reverse Flow, BS/IP/Rf (w/o Constant Certificates)

- Autoclavable, Made from 3.3 borosilicate glass


	Size	Viscosity Range mm ² /sec (cSt)*	Pack per box	Price each pc
A	2	2-10	1 Pc	4500.00
B	3	6-30	1 Pc	4500.00
C	4	20-100	1 Pc	4500.00
D	5	60-300	1 Pc	4500.00
E	6	200-1000	1 Pc	4500.00
F	7	600-3000	1 Pc	4500.00
G	8	2000-10000	1 Pc	4500.00
H	9	6000-30000	1 Pc	4500.00
I	10	20000-100000	1 Pc	4500.00
J	11	60000-300000	1 Pc	4500.00

1462 U-Tube Reverse flow Viscometers with **Constant Certificate** & Calibrated @ 40°C

- Autoclavable, Made from 3.3 borosilicate glass
- For measurement of kinematic viscosity of opaque Newtonian liquids according to ASTM D 445 and ISO 3104. Specifications conform to ASTM D 446 and ISO 3105.
- Minimum sample volume 7 ml
- comes with certificate of constant

	Size	Viscosity Range mm ² /sec (cSt)*	Approx. Constant mm ² /s ² (cSt/s)	Pack per box	Price each pc
A	2	2-10	0.01	1 Pc	11500.00
B	3	6-30	0.03	1 Pc	11500.00
C	4	20-100	0.1	1 Pc	11500.00
D	5	60-300	0.3	1 Pc	11500.00
E	6	200-1000	1	1 Pc	11500.00
F	7	600-3000	3	1 Pc	11500.00
G	8	2000-10000	10	1 Pc	11500.00
H	9	6000-30000	30	1 Pc	11500.00
I	10	20000-100000	100	1 Pc	11500.00
J	11	60000-300000	300	1 Pc	11500.00

☞ Constant on different Temp. are available against special order on additional cost.



Pipeline Filter



1481 Gas Sampling Tubes, Straight Bore Stop-cock

- Autoclavable, Made from 3.3 borosilicate glass
- have a stopcock at each end to allow for a flow of gas into the body.
- used for collecting and transporting a known Gas Sampling Tube, Glass Key, different volume of gas for testing or analysis

	Cap. ml.	Pack per box	Price each pc
A	125	1 Pc	1100.00
B	250	1 Pc	1200.00
C	500	1 Pc	1550.00
D	1000	1 Pc	2200.00




JNG

NITRILE GLOVES

BOX OF 100 PC.

RATE ON REQUEST





1491 Micro Kjeldhal Distillation Apparatus

- comprising of steam generation flask, steam trap, vacuum jacketed flask with funnel, condenser tube of silver and receiving flask, complete with heavy duty stand, brass clamps with burner

		Price each pc
A	Pregl's for determination of Nitrogen Automatic discharge according to parnes Wanger, complete with silver tube Clamp, Stand and Burner,	23000.00
B	Spare Vacuum Mental	9500.00
C	Spare Silver Condenser	13000.00



1492 Spatula


- Material stainless steel, Grade 202

	Size mm	Type	Pack per box	Price each pc
A	125	One end spoon and one end flat	10 pc	40.00
B	150	One end spoon and one end flat	10 pc	47.00
C	200	One end spoon and one end flat	10 pc	60.00
D	250	One end spoon and one end flat	10 pc	140.00
E	300	One end spoon and one end flat	10 pc	160.00

1494 Spatula Double (Tapered/Pointed)

- Material stainless steel, Grade 202
- One end tapered, & one end flat, flexible type

	Size mm	Pack per box	Price each pc
A	150	10 pc	90.00
B	200	10 pc	100.00
C	250	10 pc	145.00



1495 Spatula Micro

- Material stainless steel, Grade 202
- One end spoon and one side spatula

	Size mm	Pack per box	Price each pc
A	150	10 Pc	70.00
B	200	10 Pc	85.00





1497 Forceps (Pointed)

- Material SS 410 Grade
- Straight fine points




	Size mm	Pack per box	Price each pc
A	100	10 pc	40.00
B	125	10 pc	55.00
C	150	10 pc	70.00
D	200	10 pc	100.00
E	250	10 pc	160.00
F	300	10 pc	235.00

1499 Forceps (Blunt)

- Material SS 410 Grade

	Size mm	Pack per box	Price each pc
A	100	10 pc	40.00
B	125	10 pc	55.00
C	150	10 pc	70.00
D	200	10 pc	100.00
E	250	10 pc	160.00
F	300	10 pc	235.00




Refer Cat. No. 889




Refer Cat. No. 1370






1501 Distillation Apparatus Vertical with Electrically Heated all Glass assembly complete with stand, Clamp, Switch, Boards with Auto Cut off

	Stage	Flaks Cap Lt.	Price each pc
A	Single Stage	3 Lt	23000.00
B	Single Stage	5 Lt	26500.00
C	Single Stage	10 Lt	30000.00
D	Double Stage	3 Lt	46000.00
E	Double Stage	5 Lt	53000.00
F	Double Stage	10 Lt	60000.00
G	Triple Stage	3 Lt	69000.00
H	Triple Stage	5 Lt	79500.00
I	Triple Stage	10 Lt	90000.00
SPARE PART FOR SINGLE STAGE			
J	Flask with heater	3 Lt.	13000.00
K	Flask with heater	5 Lt	16000.00
L	Flask with heater	10 Lt	20000.00
M	Condenser	3 Lt	7000.00
N	Condenser	5 Lt	8000.00
O	Condenser	10 Lt	9000.00
P	Stand with cut off		7000.00
Q	Stand without cut off		6500.00




Refer Cat. No. 471



Refer Cat. No. 476





1506 Volatile Acid Determination Apparatus



- Autoclavable, Made from 3.3 borosilicate glass
- For the determination of soluble and insoluble volatile acids in fats and oils
- comes with 300 ml flask, stillhead, condenser and sugar flask

	Description	Pack per box	Price each pc
A	Polenske and Reichertmeissl for the determination dilute and indilute volatile acid in fats and oils. For analysis of butter and butter fats Comprising of 300 ml. flasks, stillhead condenser and receiver flask graduated at 100 and 110 ml.		2000.00
	SPARE PARTS		
B	R.M. Condenser	4 Pc	700.00
C	R. M. Still Head	4 Pc	600.00
D	R. M. Flak 300 ml.	4 Pc	450.00
E	R. M. Sugar Flask	2 Pc	550.00

1511 Essential Oil (Clevenger) Apparatus for ascertaining the percentage of essential oil (Lighter than water)

- Autoclavable, Made from 3.3 borosilicate glass

	Cap. of flask ml.	Pack per box	Price each pc
A	1000	1 Pc	6000.00
B	2000	1 Pc	7000.00
C	5000	1 Pc	12200.00




1516 Coking bulb for carbon

- for the absorption of CO₂ in carbon determination, Weight, when filled, is approximately 135 grams

		Pack per box	Price each pc
	Residue apparatus in Ramshotton methods	25 Pc	160.00






1521 Water Estimation Apparatus (Dean & Stark) Light Entrainer without Stop-cock

	Components	Price each pc
A	a. Receiver 2 ml 24/29, 19/26 b. R. B. Flask, 250 ml 24/29 c. Leibig condenser 19/26	2700.00
B	a. Receiver 5 ml 24/29, 19/26 b. R. B. Flask, 500 ml 24/29 c. Leibig condenser 19/26	2800.00
C	a. Receiver 10 ml 24/29, 19/26 b. R. B. Flask, 1000 ml 24/29 c. Leibig condenser 19/26	2700.00
D	a. Receiver 25 ml 24/29, 19/26 b. R. B. Flask 1000 ml 24/29 c. Leibig condenser 19/26	4300.00
E	a. Receiver 100 ml 24/29, 19/26 b. R. B. Flask 1000 ml 24/29 c. Leibig condenser 19/26	5000.00
	Spare Receiver for Dean & Stark	
	Cap. ml.	
F	2	2100.00
G	5	2200.00
H	10	2100.00
I	25	3800.00
J	100	4300.00


Micropipette



Refer Cat. No. JVR



J-SIL
www.j-sil.in




1526 Water Estimation Apparatus (Dean & Stark) Light Entrainer with Stop-cock

	Components	Price each pc
A	a. Receiver 2 ml 24/29, 19/26 b. R. B. Flask, 250 ml 24/29 c. Leibig condenser 19/26	2900.00
B	a. Receiver 5 ml 24/29, 19/26 b. R. B. Flask, 500 ml 24/29 c. Leibig condenser 19/26	3000.00
C	a. Receiver 10 ml 24/29, 19/26 b. R. B. Flask, 1000 ml 24/29 c. Leibig condenser 19/26	2800.00
D	a. Receiver 25 ml 24/29, 19/26 b. R. B. Flask 1000 ml 24/29 c. Leibig condenser 19/26	4500.00
E	a. Receiver 100 ml 24/29, 19/26 b. R. B. Flask 1000 ml 24/29 c. Leibig condenser 19/26	5100.00
Spare Receiver for Dean & Stark with stop cock		
	Cap ml.	
F	2	2200.00
G	5	2300.00
H	10	2200.00
I	25	4000.00
J	100	4500.00
K	10ml (for oil heavier than water)	5500.00

1531 Arsenic Determination Apparatus as per IP

- Autoclavable, Made from 3.3 borosilicate glass comes with rubber cock fitting & brass clip

	Description	Pack per box	Price each set
A	Comprising of 100 ml Borosilicate bottle arsenic tube in two parts with rubber cork & Brass Clip	1 Set	1300.00
B	As per I.P. (Rubber Cock)	1 Set	1300.00





1532 Arsenic Determination Apparatus as per BP

- Autoclavable, Made from 3.3 borosilicate glass
- comes with 100ml Conical flask & spring

	Description	Pack per box	Price each set
A	As per BP	1 Set	2100.00

1533 Arsenic Determination Apparatus as per USP

- Autoclavable, Made from 3.3 borosilicate glass
- comes with 250ml conical flask with ball joint U tube & ss clamp

	Description	Pack per box	Price each set
A	As per USP	1 Set	2100.00




1536 Arsenic Determination Apparatus Gutzeit

- Autoclavable, Made from 3.3 borosilicate glass
- comes with 100ml wide mouth bottle & brass clip

	Description	Pack per box	Price each set
	Comprising of 100 ml bottle and inner fitting. The arsenic tube is in two parts with ground flange ends and can be connected together by means of a spring or rubber bend	1 Set	3700.00



Orsat Apparatus





1546 Landsburger Boiling Apparatus (Walker Lumsden)

	Description	Price each set
	Comprising of graduated inner tube with outer jacket and conical flask with connecting tubes- Cork fitted	2800.00

1551 Orsat Gas Analysis Apparatus 3 TEST


	Description	Price each set
	Comprising of Levelling bottle, gas burette with outer jacket, three absorption pipettes, 3-test manifold, calcium chloride tube and rubber below, in wooden-case with sliding doors	28000.00



1556 Orsat Gas Analysis Apparatus 4 TEST

	Description	Price each set
	Comprising of Levelling bottle, gas burette with outer jacket, four absorption pipettes, 4-test manifold, calcium chloride tube and rubber bellow, in wooden case with sliding doors	31100.00







1561 Spares Parts for Orsat Apparatus

	Description	Pack per box	Price each pc
A	100 ml. Gas Burette with jacket	1 Pc	3000.00
B	Jacket for 100 ml gas burette	1 Pc	1200.00
C	3-Test Manifold (Orsat bridge)	1 Pc	2500.00
D	4-Test manifold (Orsat bridge)	1 Pc	2900.00
E	Absorption Pipette fitted with glass tubes	1 Pc	2100.00
F	Levelling bottle 250 ml.	2 Pc	650.00

1566 Distillation Apparatus


	Components	Pack per box	Price each pc
A	R. B. Flask with side tube 500 ml 24/29 Coil condenser 19/26 Stopper 24/29	1 Set	2700.00
B	Spare condenser	1 Pc	1000.00
C	Spare flask	1 Pc	1600.00



1571 Ammonia Distillation Assembly

	Components	Pack per box	Price each pc
A	Kjeldhal flask 500 ml B24 Splash head 24/29, 19/26 Coil condenser 19/26		2700.00
B	Spare splash head	1 Pc	1100.00
C	Spare condenser	1 Pc	1000.00






1576 Distillation Apparatus, Compact

	Components	Price each pc
A	Conical flask 1000 ml Friedrich Condenser 34/35	4800.00

1581 Distillation Apparatus with Friedrich Condenser

	Components	Price each pc
A	R. B. Flask with side tube 2000 ml 24/29 Friedrich Condenser Stopper 24/29	5300.00



Weighing Bottle



Weighing Scoop







Refer Cat. No. 1181 & 1182






1586 Distillation Assemblies Similar to Quickfit



	Components	Price each pc
A	R. B. Flask 2000 ml 24/29 Splash head 24/29, 24/29 Double surface condenser, 200 mm 24/29, 24/29 Receiver adapter Bend 24/29	4950.00
B	R. B. Flask 500 ml 24/29 Still head 24/29, 19/26, 19/26 Leibig condenser 19/26, 19/26 Receiver adapter Plain Bend 19/26 Separating funnel 100 ml 19/26 Stopper 19/26	3600.00
C	R. B. Flask 1000 ml 24/29 Recovery Bend, Vertical 24/29, 24/29 Thrope condenser, 200 mm 24/29 Receiver adapter Straight 24/29 (Alcohol distillation unit)	5700.00
D	R. B. Flask, 250 ml 24/29 Recovery Bend 24/29, 19/26 Leibig condenser, 200 mm 19/26 Receiver adapter Plain Bend 19/26	1700.00

1596 Preparation Assembly Similar to Quickfit


	Components	Price each pc
A	R. B. Flask 500 ml 24/29 Still head 24/29, 24/29, 14/23 Thermometer 14/23 Leibig condenser 24/29 Receiver adapter Bend 24/29	3200.00





1611 Vaccum Distillation Assemblies, Similar to Quickfit

	Components	Price each pc
	R. B. Flask 1000 ml Recovery bend Jacketed coil condenser, 150 mm Receiver adapter, bend with vent R. B. Flask, 250 ml	5500.00



Refer Cat. No. 71




Refer Cat. No. 76





Assemblies






1616 Kjeldhal Distillation Assemblies

	Components	Price each pc
A	R. B. Flask 1000 ml Multiple adapter two necks Splash head vertical Leibig condenser Receiver adapter Tap funnel	3750.00

1621 Kjeldhal Distillation Assemblies with 800 ml Flask



	Components	Price each pc
A	Kjeldhal flask 800ml. B24 B19 B24 Multiple adapter two necks B19 B19 Splash head vertical 300mm B19 B19 Leibig condenser Receiver adapter St. B19 Top funnel 100ml B19	3900.00



1626 Steam Distillation Assemblies

	Components	Price each pc
A	R. B. Flask two necks 500 ml Splash head Leibig condenser Receiver adapter Air leak tube	3100.00
B	R. B. Flask, 500 ml Still Head Air leak tube Leibig condenser Receiver adapter	2950.00







1631 Steam Distillation Assemblies

	Components	Price each pc
A	R. B. Flask, 500 ml Steam distillation head Leibig condenser Receiver adapter	3500.00
B	R. B. Flask 1000 ml Multiple adapter, two necks Air leak tube, bent Recovery bend Jacketed coil condenser, 150 mm Receiver adapter	3950.00

1636 Sulphur Dioxide Apparatus as per USP

	Cap. ml.	Price each pc
A	500	6950.00




1646 Reflux assemblies, Similar to Quickfit

	Components	Price each pc
A	R. B. Flask, 100 ml. Air condenser 500 mm	800.00
B	R. B. Flask, 250 ml. Air condenser, 750 mm	990.00





Assemblies





1671 Gas Absorption assembly

	Components	Price each pc
A	2 -U-tubes with Socket 1-u-tube with cone Drying tube Cone adapter, bent	1650.00

1676 Vaccum Sublimation assembly

	Components	Price each pc
A	Test tube Adapter with 'T' connection Immersion condenser	1050.00



1691 Kulderna Danish Evaporator

	Components	Price each pc
A	R.B. Flask, 150 ml. Upper chamber 250 ml Kulderna column	4900.00





Refer Cat. No. 81



J-SIL
www.j-sil.in

Soxhlet Apparatus




1706 Soxhlet Extraction apparatus, Complete

	Components	Price each pc
A	Extractor 60 ml F.B. Flask 100 ml Allihn condenser 150 mm	1750.00
B	Extractor 100 ml F.B. Flask 250 ml Allihn condenser 200 mm	2000.00
C	Extractor 200 ml F. B. Flask 500 ml Allihn condenser 250 mm	3450.00
D	Extractor 400 ml F.B. Flask 1000 ml Allihn condenser 250 mm	5700.00
E	Extractor 600 ml F.B. Flask 2000 ml Allihn condenser 300 mm	8850.00
F	Extractor 1000 ml F. B. Flask 3000 ml Allihn condenser 300 mm	11750.00

☞ For Spare Parts, Kindly Refer Cat No. 1146, 1151, 1156

Micropipette




Refer Cat. No. JVR






Assemblies

1731 Recovery Assemblies Similar to Quickfit





	Components	Price each pc
A	R. B. Flask 250 ml Multiple Adapter two necks Recovery Bend Leibig condenser Separating funnel 100 ml	3350.00
B	R. B. Flask 2000 ml Still head Jacketed coil condenser, 150 mm Receiver adapter Separating funnel 100 ml Stopper	7100.00



1736 Solvent Recovery assemblies Similar to Quickfit


	Components	Price each pc
A	R. B. Flask 150 ml Recovery bend Jacketed coil condenser, 150 mm Receiver adapter	4100.00



1746 Chromatography Apparatus, SMALL SIZE Similar to Quickfit

	Components	Price each pc
A	Buchner Flask 100 ml Column 200 X 10 mm with Sintered disc Reservoir 50 ml Cone adapter	4500.00
B	300 mm X 10 mm columns	5000.00






1751 Chromatography Apparatus, MEDIUM SIZE Similar to Quickfit

	Components	Price each pc
A	Buchner Flask 250 ml Tap adapter Column 200 mm X 18 mm with Sintered disc Reduction adapter Reservoir, 250 ml	5000.00
B	300 mm X 18 mm columns	5350.00
C	400 mm X 18 mm columns	5550.00

1752 Andreasen apparatus for determining particle size by the gravity sedimentation method

		Price each pc
A	A much quicker and less tedious Method of sedimentation analysis is the Andreasen Pipette Technique	4700.00



Vials



Refer Cat No. 61 & 66 on Page No. 1







Filter Holder Assembly



1891 Millipore Type all Glass Filter Holder Assembly with I/C Joints

- Manufactured from borosilicate glass particularly suitable for analysis of microscopic residues, particle controls and multipurpose filtration application. Thick wall heavy base and reinforced rim prevent breakage during washing and transportation.
- Suitable to be used with 47/50 mm diameter membrane filters. Supplied complete with filter funnel of 300 ml volume capacity, support base with sintered disc, vacuum connection side arm, connecting clamp and vacuum flask of 1000 ml volume capacity.

	Filter mm.	Flask ml.	Price each pc
A	13	250 B29	5000.00
B	25	500 B29	5300.00
C	47	1000 B40	7500.00
D	47	2000 B40	8900.00
E	90	5000 B55	49500.00
F	Spare Flask	1000 B40	2000.00
G	Spare Flask	2000 B40	3000.00
H	Spare Flask	5000 B55	15000.00
I	Spare Cup	300ml 47mm	2900.00
J	Spare Cup	500ml 47mm	5200.00
K	Spare Cup	1000ml 90mm	15500.00
L	Spare Sintered Part	47mm B40	3000.00
M	Spare Sintered Part	90mm B55	25000.00
N	Spare Clamp	47mm	2000.00
O	Spare Clamp	90mm	10000.00




Filter Holder Assembly



1892 All Glass Filter Holder Assembly with Silicone Stopper.

- Made for borosilicate glass with glass filter.
- The filtration base has a silicone stopper.

	Filter mm.	Filter Flask	Price each pc
A	13	250 ml	5000.00
B	25	500 ml	5300.00
C	47	1000 ml	7500.00
D	47	2000 ml	9500.00
E	Spare cup	35 ml	1400.00
F	Spare cup	50 ml	1500.00
G	Spare cup	300 ml	2900.00
H	Spare sintered part	13 mm	2700.00
I	Spare sintered part	25 mm	2800.00
J	Spare sintered part	47 mm	3000.00
K	Spare clamp	25 mm	2000.00
L	Spare clamp	47 mm	2000.00
M	SS filter Holder only	47 mm	6000.00



Refer Cat. No. 1892M



Refer Cat. No. 1892J



Refer Cat. No. 951





Dessicators



1896 Kipp's Apparatus Complete

	Cap. ml.	Pack per box	Price each pc
A	500	1 Set	6000.00
B	1000	1 Set	7000.00

1906 Glass Dessicators **plain with Lid** without Mark (Imported)

- Desiccators allow you to cool and store objects in a low-humidity atmosphere.
- its cone with round glass knob

	Inner Dia.	Pack per box	Price without Mark
A	100/120 mm	1 Pc	1700.00
B	150 mm	1 Pc	3300.00
C	210 mm	1 Pc	5100.00
D	240 mm	1 Pc	7000.00
E	300 mm	1 Pc	12500.00



1911 **Vaccum** Dessicators with Lid Imported without Mark

- can be used for precision production test
- it is specially useful for volatile liquid.

	Inner Dia.	Pack per box	Price each pc
A	150 mm.	1 Pc	5200.00
B	210 mm.	1 Pc	8200.00
C	240 mm.	1 Pc	9400.00
D	300 mm.	1 Pc	16500.00





1913 Desiccators plain PP

- Made from polypropylene and Polycarbonate
- Provided with a thick perforate polypropylene disc

	Inner Dia.	Pack per box	Price each pc
A	150 mm	1 Pc	2220.00
B	200 mm	1 Pc	4050.00
C	250 mm	1 Pc	5400.00
D	300 mm	1 Pc	9600.00


1915 Desiccators Vaccum PP

- Made from polypropylene and Polycarbonate
- Fitted with Polypropylene stopcock with PTFE Plug
- Autoclavable vaccum desiccator can hold vaccum upto 740 mm of HG for 24 Hrs without any greasing

	Inner Dia.	Pack per box	Price each pc
A	150 mm	1 Pc	3750.00
B	200 mm	1 Pc	5850.00
C	250 mm	1 Pc	7500.00
D	300 mm	1 Pc	11700.00



Refer Cat. No. 51



Refer Cat. No. 56



Refer Cat. No. 61





INDEX

ITEM NAME	PAGE NO	CAT. NO
VIALS & SEPTA	1	
SYRINGE FILTER	2-4	
NYLON SYRINGE FILTER	2	71
PVDF SYRINGE FILTER	3	76
PTFE SYRINGE FILTER	4	81
MICROPIPETTES	5-9	
FIXED MICROPIPETTE FOR SCHOOL, COLLEGES & EDUCATION LABS	5	JFE
VARIABLE MICROPIPETTE FOR SCHOOL, COLLEGES & EDUCATION LABS	6	JVE
FIXED MICROPIPETTE FOR QC & RND LABS.	7	JFR
VARIABLE MICROPIPETTE FOR QC & RND LABS.	8	JVR
NANO MICROPIPETTE FOR DIAGNOSTIC KITS.	9	JNP
SAFETY RINGS FOR VOLUMETRIC FLASK	9	RVF
BURETTES	10 - 18	
BURETTES WITH GLASS STOP-COCK CLASS B	10	101
BURETTES WITH ROTA FLOW CLASS B	10	106
BURETTES WITH AMBER COLOUR ROTA FLOW CLASS B	11	107
BURETTES WITH TEFLON KEY STOP-COCK STRAIGHT BORE CLASS B	11	111
BURETTES WITH TEFLON KEY STOP-COCK STRAIGHT BORE CLASS A	12	116
USP BURETTES WITH TEFLON KEY STOP-COCK STRAIGHT BORE CLASS A	12	117
BURETTES WITH GLASS KEY STOP-COCK STRAIGHT BORE CLASS A	13	121
USP BURETTES WITH TEFLON KEY (ROTA FLOW)SINGLE BORE CLASS A	13	125
BURETTE SINGLE BORE ROTA FLOW CLASS A	14	126
BURETTE AMBER COLOUR SINGLE BORE ROTA FLOW CLASS A	14	127
USP BURETTE AMBER COLOUR SINGLE BORE ROTA FLOW CLASS A	15	128
BURETTE FOR PINCK COCK CLASS B	15	131
FISHER CLAMP	15	132
BURETTE AUTOMATIC WITH GLASS STOP COCK CLASS B	16	176
BURETTE AUTOMATIC WITH TEFLON STOP COCK CLASS B	16	181
BURETTE AMBER COLOUR AUTOMATIC ZERO TEFLON STOP COCK CLASS B	17	182
BURETTE AUTOMATIC WITH GLASS STOP COCK CLASS A	17	186
BURETTE AUTOMATIC WITH TEFLON STOP COCK CLASS A	18	191
BURETTE AMBER COLOUR AUTOMATIC WITH TEFLON STOP COCK CLASS A	18	192
BURETTES CLAMP	11 & 15	112, 113, 132
MEASURING CYLINDERS	19 - 23	
MEASURING CYLINDERS WITH ROUND BASE CLASS B	19	196
MEASURING CYLINDERS WITH HEXAGONAL BASE CLASS B	19	197
MEASURING CYLINDERS WITH HEXAGONAL BASE CLASS A	20	202
USP MEASURING CYLINDERS WITH WORKS CERTIFICATE CLASS A	21	204
AMBER MEASURING CYLINDERS WITH HEXAGONAL BASE CLASS A	21	205
MEASURING CYLINDERS WITH ROUND BASE & STOPPER CLASS B	22	206
MEASURING CYLINDERS WITH HEXAGONAL BASE & STOPPER CLASS B	22	207
MEASURING CYLINDERS WITH HEXAGONAL BASE & STOPPER CLASS A	23	212
USP MEASURING CYLINDERS WITH STOPPER CLASS A	23	213
USP TAPPED DENSITY CYLINDERS CLASS A	20	203
CYLINDERS FOR SUSPENSIBILITY TEST CLASS B	24	216
CYLINDERS FOR SUSPENSIBILITY TEST CLASS A	24	217
POLYPROPYLENE CYLINDERS WITH HEXAGONAL BASE	24	218
NESSLER CYLINDERS	25	221
NESSLER CYLINDERS CLASS A	25	222
CROW RECEIVER CLASS B	25	226
CROW RECEIVER WITH I/C STOPPER CLASS B	26	231
CROW RECEIVER CLASS A	26	236
CROW RECEIVER WITH I/C STOPPER CLASS A	26	241
VOLUMETRIC FLASK	27-33	
VOLUMETRIC FLASK CLASS B	27	246
AMBER VOLUMETRIC FLASK CLASS B	27	247
VOLUMETRIC FLASK WITH LDPE STOPPER CLASS B	28	251



INDEX



ITEM NAME	PAGE NO	CAT. NO
VOLUMETRIC FLASK CLASS A	29	256
AMBER VOLUMETRIC FLASK CLASS A	30	257
USP VOLUMETRIC FLASK CLASS A	30	261
USP AMBER VOLUMETRIC FLASK CLASS A	31	262
AMBER VOLUMETRIC FLASK WIDE MOUTH CLASS A	31	266
WIDE MOUTH VOLUMETRIC FLASK CLASS A	32	267
USP VOLUMETRIC FLASK WIDE MOUTH CLASS A	32	268
AMBER USP VOLUMETRIC FLASK WIDE MOUTH CLASS A	33	269
SUGAR FLASK WITH TWO MARKS CLASS B	33	271
SUGAR FLASK WITH TWO MARKS CLASS A	33	276
MUD FLASK CLASS B	34	281
MUD FLASK CLASS A	34	286
VISCOSITY FLASK	34	291
CASSIA FLASK GRADUATED LONG NECK CLASS B	35	296
CASSIA FLASK GRADUATED LONG NECK CLASS A	35	297
LE-CHATETIER FLASK	35	301
LE-CHATETIER FLASK WITH WORK CERTIFICATE CLASS A	36	306
PHENOL FLASK	36	311
PIPETTES	37-42	
VOLUMETRIC PIPETTES CLASS B	37	316
VOLUMETRIC PIPETTES CLASS A	37	321
AMBER VOLUMETRIC PIPETTES CLASS A	38	322
USP VOLUMETRIC PIPETTES CLASS A	38	323
SEROLOGICAL PIPETTES CLASS B	39	326
SEROLOGICAL PIPETTES CLASS A	39	331
USP GRADUATED PIPETTES (SEROLOGICAL) CLASS A	40	332
AMBER COLOUR GRADUATED PIPETTE CLASS A	40	336
MOHR TYPE GRADUATED PIPETTES CLASS B	41	341
MOHR TYPE GRADUATED PIPETTES CLASS A	41	346
USP MOHR TYPE GRADUATED PIPETTES CLASS A	42	347
LAMBDA PIPETTES CLASS A	42	356
MILK TEST PIPETTES , BACTERIOLOGICAL	43	366
MILK TEST PIPETTES, GERBER METHOD	43	371
PIPETTES STAND HORIZONTAL POLYPROPYLENE	43	372
PIPETTES STAND VERTICAL POLYPROPYLENE	43	373
TUBES	44-52	
CENTRIFUGE TUBE PLAIN WITH RIM	44	376
CENTRIFUGE TUBE PLAIN WITH I/C STOPPER	44	381
CENTRIFUGE TUBE GRADUATED WITH RIM	44	386
CENTRIFUGE TUBE GRADUATED & I/C STOPPER	44	391
CENTRIFUGE TUBE ROUND BOTTOM GRADUATED	45	396
CENTRIFUGE TUBE POLYPROPYLENE & AUTOCLAVABLE	45	397
CENTRIFUGE TUBE ROUND BOTTOM POLYPROPYLENE	45	398
OAK RIDGE CENTRIFUGE TUBES POLYPROPYLENE	45	399
RACK FOR MICRO CENTRIFUGE TUBE	46	400
CENTRIFUGE VESSEL (SPHERICAL)	46	401
MICRO CENTRIFUGE TUBE (PP)	46	402
CENTRIFUGE VESSEL (CYLINDRICAL)	46	406
CENTRIFUGE TUBE VESSEL CLASS A	47	407
TEST TUBE WITH RIM	47	408
TEST TUBE WITHOUT RIM	47	408A
TUBE CULTURE	48-49	
CULTURE TUBES ROUND BOTTOM WITH SCREW CAP	48	409
CULTURE TUBES FLAT BOTTOM WITH SCREW CAP	48	410
TEST TUBE WITH RIM GRADUATED	48	411
GLASS ROD FOR STIRRING	49	412
AMBER CULTURE TUBES FLAT BOTTOM WITH SCREW CAP	49	413
AMBER CULTURE TUBES ROUND BOTTOM WITH SCREW CAP	49	414
TEST TUBE STOPPERED WITH HOLLOW STOPPER	50	416
TEST TUBE STOPPERED WITH PENNY HEAD SOLID STOPPER	50	417
TEST TUBE GRADUATED, STOPPERED WITH HOLLOW STOPPER	51	421





INDEX

ITEM NAME	PAGE NO	CAT. NO
TEST TUBE GRADUATED, STOPPERED WITH SOLID STOPPER	51	422
STORAGE VIAL (POLYPROPYLENE)	51	423
SV5 VIAL (POLYPROPYLENE)	52	424
TEST TUBE WITH SCREW CAP, (POLYPROPYLENE)	52	425
BLOOD SUGAR TUBE	52	426
TEST TUBE STAND POLYPROPYLENE & AUTOCLAVABLE	53	427
TEST TUBE BASKETS POLYPROPYLENE	53	428
DRYING TUBES	70	691
SPECIFIC GRAVITY BOTTLES, PYKNOMETERS	54-55	
SP. GR. BOTTLES WITH GLASS STOPPER	54	431
SP. GR. BOTTLES WITH TEFLON STOPPER	54	436
SP. GR. BOTTLES WITH WORK CERTIFICATE CLASS A	54	441
SP. GR. BOTTLES WITH THERMOMETER (PYKNOMETER)	55	446
SP. GR. BOTTLES WITH CERTIFIED THERMOMETER (PYKNOMETER) CLASS A	55	447
BEAKERS	56- 58	
BEAKERS WITH SPOUT DOUBLE GRADUATED	57	471
BEAKERS AMBER COLOUR WITH SPOUT	57	472
POLYPROPYLENE BEAKER	56	451
MEASURING JUGS WITH HANDLE POLYPROPYLENE	56	456
BEAKERS LOW FORMS WITH SPOUT	56	466
BEAKERS TALL FORMS WITH GRADUATION	58	474
BEAKER & FLASK TONGS	76	889
CONICAL FLASK	59, 83-85	
CONICAL FLASK ERLLENMEYER,PLAIN	59	476
CONICAL FLASK WITH SCREW CAP	59	481
BAFFLED ERLLENMEYER FLASK	60	486
CONICAL FLASK WITH I/C JOINT	83	931
STOPPERED CONICAL FLASK	84	936
CONICAL FLASK FOR FILTER ASSEMBLY	145	1891
BEAKER & FLASK TONGS	76	889
ADAPTERS	61-70	
ADAPTERS REDUCTION	61	496
ADAPTERS EXPANSION	61	501
ADAPTERS MULTIPLE , TWO NECKS PARALLEL	61	506
ADAPTERS TWO NECK ONE VERTICAL & ONE AT ANGLE	62	511
ADAPTERS THREE NECK, TWO PARALLEL & ONE AT ANGLE	62	516
ADAPTERS SWAN NECK FOR USE WITH AIR LEAK TUBE	62	521
ADAPTER WITH 'T' CONNECTION	62	526
CONE ADAPTERS, STRAIGHT CONNECTION PLAIN	63	531
CONE ADAPTERS, STRAIGHT CONNECTION WITH STOP-COCK	63	536
CONE ADAPTERS ANGLE CONNECTION	63	541
CONE ADAPTERS, RIGHT ANGLE CONNECTION WITH STOP-COCK	63	546
SOCKET ADAPTERS STRAIGHT CONNECTION WITH STOP-COCK	64	551
ADAPTERS RECEIVER PLAIN BEND SHORT STEM	64	556
ADAPTERS RECEIVER PLAIN BEND LONG STEM	64	561
ADAPTERS RECEIVER STRAIGHT DELIVERY	64	566
ADAPTER RECEIVER, PLAIN BEND, SOCKET TO CONE	65	571
ADAPTER RECEIVER, BEND WITH VENT	65	576
ADAPTER RECEIVER, BEND WITH VACCUM CONNECTION	65	581
ADAPTER RECEIVER, STRAIGHT WITH VACCUM CONNECTION	65	586
ADAPTER RECEIVER, BEND WITH SOCKET CONNECTION	66	591
RECOVERY BEND SLOPING END	66	621
RECOVERY BEND VERTICAL	66	626
PLAIN STILL HEAD, SLOPING	67	631
CLAISEN HEAD, SLOPING	67	636
SPLASH HEAD ADAPTER	67	641
SPLASH HEAD SLOPING	68	646
SPLASH HEAD VERTICAL	68	651
STEAM DISTILLATION HEAD, SLOPING	68	656



INDEX



ITEM NAME	PAGE NO	CAT. NO
THERMOMETERS WITH STEM & TEFLON CONE	68	671
THERMOMETERS POCKET	69	676
AIR/STEM INLET TUBES, STRAIGHT	69	681
AIR/STEM INLET TUBES, BENT	69	686
DRYING TUBES	70	691
STOP COCK ADAPTER FOR VACCUM RELEASE	70	696
STOP COCK ADAPTER FOR VACCUM CONTROL	70	701
CONDENSERS	71 -74	
CONDENSERS LEIBIG	71	711
CONDENSERS DOUBLE SURFACE	71	726
CONDENSERS AIR WITH I/C JOINT	71	751
CONDENSERS AIR VERTICAL	72	756
CONDENSERS LEIBIG WITH I/C JOINT	72	761
CONDENSERS BULB WITH I/C JOINT	72	771
CONDENSERS COIL (GRAHAM) WITH I/C JOINT	73	781
CONDENSERS, JACKETED COIL WITH I/C JOINT	73	796
CONDENSERS COILED DISTILLAGE WITH I/C JOINT	73	801
CONDENSERS COMBINED DOUBLE COIL	74	811
CONDENSERS INLAND REVENUE (THORPE)	74	816
CONDENSERS IMMERSION (COLD FINGER)	74	831
COLUMNS	75-76	
COLUMNS, CHROMATOGRAPHY WITH TEFLON ROTA FLOW STOP-COCK	75	846
COLUMNS, CHROMATOGRAPHY WITH SINTERED DISC	75	861
COLUMNS, CHROMATOGRAPHY WITH INTEGRAL SINTERED	75	866
FRACTIONATING COLUMNS. VIGREX	76	876
FRACTIONATING COLUMNS, PLAIN	76	881
SPRAYER	76	
SPRAYER FOR CHROMATOGRAPHY WITH RUBBER BULB & CLIP	76	886
FLASK	77-84	
ROUND BOTTOM FLASKS WITH I/C JOINT	77	891
FLAT BOTTOM FLASKS WITH I/C JOINT	78	896
ROUND BOTTOM FLASKS, PLAIN NECK WITH BEADED	79	901
FLAT BOTTOM FLASKS, PLAIN NECK WITH BEADED	79	902
ROUND BOTTOM WITH TWO NECK	79	906
ROUND BOTTOM WITH THREE NECK	80	911
ROUND BOTTOM WITH FOUR NECK	80	916
ROUND BOTTOM FLASKS FOUR NECK WITH THERMOMETER POCKET	81	918
BEAKER & FLASK TONGS	76	889
FLASK STAND AUTOCLAVABLE PP	81	919
TEFLON BUSH/ GLAND WITH O RING	81	920
KJELDHAL FLASKS WITH I/C JOINT	82	921
KJELDHAL FLASKS, LONG NECK PLAIN	82	926
CONICAL FLASKS WITH I/C JOINT	83	931
STOPPERED CONICAL FLASKS	84	936
PEAR SHAPPED FLASKS, SINGLE NECK	84	941
FILTER FLASK WITH I/C JOINT	84	946
FLASK FILTER- CONICAL FLASK WITH SIDE TUBE, HEAVY WALL	85	951
IODINE FLASKS, CUP SHAPED	85	956
OXYGEN FLASKS WITH PLATINUM BASKET	85	957
PEAR SHAPE FLASKS FOR ROTARY EVAPORATORS WITH I/C JOINT	86	960
DISTILLATION FLASK	86	976
DISTILLATION FLASK AS PER ASTM SPECIFICATION	86	981
FUNNELS	87- 88	
FUNNEL WITH THISTLE HEAD	87	986
FILTER FUNNELS WITH LONG STEM	87	991
FILTER FUNNELS WITH SHORT STEM	87	996
ANALYTICAL FUNNEL POLYPROPYLENE	88	997
FUNNEL HOLDER POLYPROPYLENE	88	998





INDEX

ITEM NAME	PAGE NO	CAT. NO
FILTER PUMP	89	1001
FILTER TUBE	89	1006
DIGESTION TUBE	89	1011
CONES & SOCKET	89- 90	
SOCKET, PLAIN SHANK DOUBLE	89	1016
CONE, PLAIN END DOUBLE	90	1031
CONES, LONG TIP DOUBLE	90	1041
SPHERICAL JOINTS, CUP	90	1046
SPHERICAL JOINTS, BALL	90	1051
IMPINGER & STALGOMETER	91	
STANDARD IMPINGER	91	1056
STALGNOMETER	91	1061
STIRRERS	92-93	
PTFE BLADED STIRRERS	92	1066
PTFE BLADED STIRRERS GLASS ROD ONLY	93	1067
GLASS ROD FOR STIRRING	48	412
STOP COCKS	94-95	
STOP COCKS, STRAIGHT BORE WITH SOLID STOPPER PLUG	94	1086
STOP COCKS STRAIGHT BORE WITH TEFLON PLUNGER	94	1091
STOP COCKS FOR BURETTES	94	1096
STOP-COCK TEFLON ROTAFLOW FOR BURETTE, SEP . FUNNEL OR COLUMN	94	1101
STOP-COCK, DOUBLE OBLIQUE BORE SOLID STOPPER PLUG	95	1106
STOP-COCK 'T' BORE SOLID STOPPER PULG BORE	95	1111
STOPPERS	96	
STOPPERS SOLID PENNY HEAD	96	1140
STOPPERS HEXAGONAL HOLLOW	96	1141
SOXHLET	97- 142	
EXTRACTOR FOR SOXHLET	97	1146
CONDENSER ONLY FOR SOXHLET	97	1151
EXTRACTION FLASKS	97	1156
SOXHLET EXTRACTION APPARATUS, COMPLETE	142	1706
BOTTLES	98-105	
B.O.D. BOTTLES WITH INTERCHANGEABLE STOPPER	99	1163
B.O.D. BOTTLES WITH I/C STOPPER AMBER COLOUR	99	1164
REAGENT BOTTLE, NARROW MOUTH WITH I/C STOPPER FLAT HEAD	98	1161
REAGENT BOTTLE, AMBER COLOUR NARROW MOUTH WITH I/C STOPPER FLAT HEAD	98	1162
REAGENT BOTTLE WITH P.P SCREW CAP WIDE MOUTH	100	1166
REAGENT BOTTLE AMBER COLOUR WITH P.P. SCREW CAP WIDE MOUTH	100	1167
SPARE SCREW CAP FOR MOBILE PHASE BOTTLE	101	1169
REAGENT BOTTLE WIDE MOUTH WITH I/C HOLLOW STOPPER	101	1171
DROPPING BOTTLE WITH PIPETTE & RUBBER TEAT	103	1201
DROPPING BOTTLE AMBER COLOUR WITH PIPETTE & RUBBER TEAT	103	1202
WEIGHING BOTTLES WITH I/C STOPPER	102	1181
LIQUID WEIGHING BOTTLE WITH RUBBER TEAT	102	1196
BOTTLES, WASH, PE	104	1205
WASH BOTTLE COMPLETE	104	1206
WASH BOTTLE FITTING ONLY	104	1211
GAS WASH BOTTLE WITH INTERCHANGEABLE FITTING	105	1216
GAS WASH BOTTLE ONLY	105	1221
GAS WASH BOTTLE, FITTING ONLY	105	1226
GAS WASH BOTTLE WITH SINTERED END FOR DOWNWARD GAS DISTRIBUTION	105	1231
GLASS WEIGHING SCOOP/BOAT	102	1182





ITEM NAME	PAGE NO	CAT. NO
SEPARATING FUNNELS	106-110	
SEPARATING FUNNELS WITH GLASS STOP COCK (PEAR SHAPE)	106	1251
SEPARATING FUNNELS WITH TEFLON ROTAFLOW PLUG (PEAR SHAPE)	106	1256
SEPARATING FUNNELS AMBER WITH TEFLON STOP COCK	106	1257
SEPARATING FUNNELS WITH PTFE PLUNGER STOP-COCK	107	1261
SEPARATING FUNNELS WITH GLASS STOP-COCK, GLOBE SHAPE	107	1266
SEPARATING FUNNELS WITH TEFLON STOP-COCK, GLOBE SHAPE	107	1271
SEPARATING FUNNELS HOLDER, POLYPROPYLENE	108	1272
SEPARATING FUNNELS POLYPROPYLENE	108	1273
SEPARATING FUNNELS WITH TEFLON PLUG	108	1282
SEPARATING FUNNELS WITH SOCKET & CONE WITH GLASS STOP-COCK	109	1291
SEPARATING FUNNELS OPEN TOP	109	1296
SEPARATING FUNNEL CYLINDRICAL WITH SOCKET, CONE & TEFLON PLUNGER STOP COCK	110	1301
PRESSURE EQUALISING FUNNELS WITH TEFLON STOP COCK	110	1306
PRESSURE EQUALISING FUNNELS WITH GLASS STOP COCK	110	1311
TILT MEASURE	111	
POUR POINT TUBE	111	1316
AUTOMATIC TILT MEASURES	111	1321
DISHES	112-113	
PETRI DISHES PLAIN	112	1326
EVAPORATING DISHES WITH FLAT BOTTOM WITH POUR OUT	112	1327
CRYSTALLIZING DISHES	112	1328
EVAPORATING DISHES PORCELAIN	113	1329
DISPOSABLE PETRI DISH (GAMA STERILE)	113	1330
MELTING POINT TUBE, THIEL TUBE	114	1336
CALCIUM CHLORIDE TUBES STRAIGHT WITH ONE BULB	114	1346
CALCIUM CHLORIDE TOWERS	114	1356
CONNECTION TUBES, T-SHAPE	114	1361
CONNECTION TUBES, Y-SHAPE	114	1366
GOOCH CRUCIBLE	115	
TONG CRUCIBLE	115	1370
GOOCH CRUCIBLE WITH SINTERED DISC	115	1371
ADAPTER FOR GOOCH CRUCIBLE	115	1376
BUCKNER FUNNEL	116, 118 - 119	
BUCKNER FUNNEL WITH CONE ON STEM	116	1380
BUCKNER FILTER WITH SINTERED DISH & PLAIN STEM	118	1391
FUNNEL BUCHNER PORCELAIN	119	1393
FUNNEL BUCHNER POLYPROPYLENE	119	1394
SILICA & QUARTZ CRUCIBLE	117	
CRUCIBLE SILICA LID ONLY	117	1387
CRUCIBLE SILICA	117	1388
QUARTZ CRUCIBLE ONLY	117	1389
QUARTZ CRUCIBLE LID ONLY	117	1390
VISCOMETERS	120-125	
VISCOMETERS COLLEGE PATTERN, OSTWALD	120	1421
VISCOMETERS B.S.S.	120	1426
VISCOMETERS U-TUBES BS/U (W/O CONSTANT CERTIFICATE)	121	1431
VISCOMETERS U-TUBES BS/U WITH CONSTANT CERTIFICATE	121	1432
VISCOMETERS SUSPENDED LEVEL BS/IP/SL (W/O CONSTANT CERTIFICATE)	122	1436
VISCOMETERS SUSPENDED LEVEL BS/IP/SL WITH CONSTANT CERTIFICATE	122	1437
VISCOMETERS CANNON-FENSKE TYPE BS/IP/CF (W/O CONSTANT CERTIFICATE)	123	1446
VISCOMETERS CANNON-FANSKE TYPE BS/IP/CF WITH CONSTANT CERTIFICATE	123	1451
VISCOMETERS CANNON-FANSKE FOR OPAQUE LIQUIDS (W/O CONSTANT CERTIFICATE)	124	1456
VISCOMETERS C/F FOR OPAQUE LIQUID WITH CONSTANT CERTIFICATE	124	1457
VISCOMETERS U-TUBES, REVERSE FLOW, BS/IP/RF (W/O CONSTANT CERTIFICATE)	125	1461
VISCOMETERS U-TUBES REVERSE FLOW BS/IP/RF WITH CONSTANT CERTIFICATE	125	1462
GLASS SAMPLING TUBES, STRAIGHT BORE GLASS STOP-COCK	126	1481





INDEX

ITEM NAME	PAGE NO	CAT. NO
SPATULA & FORCEPS	127-128	
MICRO KJELDHAL DISTILLATION APPARATUS	127	1491
SPATULA S.S	127	1492
SPATULA DOUBLE (TAPERED/POINTED)	127	1494
MICRO SPATULA	127	1495
FORCEP (POINTED)	128	1497
FORCEP (BLUNT)	128	1499
APPARATUS	129-136	
DISTILLATION APPARATUS VERTICAL	129	1501
VOLATILE ACID DETERMINATION APPARATUS	130	1506
ESSENTIAL OIL APPARATUS	130	1511
COKING BULB FOR CARBON	130	1516
WATER ESTIMATION APPARATUS (DEAN & STARK) WITHOUT STOP COCK	131	1521
WATER ESTIMATION APPARATUS (DEAN & STARK) WITH STOP COCK	132	1526
ARSENIC DETERMINATION APPARATUS AS PER IP	132	1531
ARSENIC DETERMINATION APPARATUS AS PER BP	133	1532
ARSENIC DETERMINATION APPARATUS AS PER USP	133	1533
ARSENIC DETERMINATION APPARATUS GUTZEIT	133	1536
LANDSBURGER BOILING APPARATUS (WALKER LUMSDEN)	134	1546
ORSAT GAS ANALYSIS APPARATUS 3 TEST	134	1551
ORSAT GAS ANALYSIS APPARATUS 4 TEST	134	1556
SPARES PARTS FOR ORSAT APPARATUS	135	1561
DISTILLATION APPARATUS	135	1566
AMMONIA DISTILLATION ASSEMBLY	135	1571
DISTILLATION APPARATUS, COMPACT	136	1576
DISTILLATION APPARATUS WITH FRIEDRICH CONDENSER	136	1581
DISTILLATION ASSEMBLIES SIMILAR TO QUICKFIT	137	1586
PREPARATION ASSEMBLY SIMILAR TO QUICKFIT	137	1596
ASSEMBLIES	137 - 144	
VACCUM DISTILLATION ASSEMBLIES	138	1611
KJELDHAL DISTILLATION ASSEMBLIES	139	1616
KJELDHAL DISTILLATION ASSEMBLIES WITH 800 ML FLASK	139	1621
STEAM DISTILLATION ASSEMBLIES	139	1626
STEAM DISTILLATION ASSEMBLIES	140	1631
SULPHUR DIOXIDE APPARATUS AS PER USP	140	1636
REFLUX ASSEMBLIES	140	1646
GAS ABSORPTION ASSEMBLY	141	1671
VACCUM SUBLIMATION ASSEMBLY	141	1676
KULDERNA DANISH EVAPORATOR	141	1691
SOXHLET EXTRACTION APPARATUS, COMPLETE	142	1706
RECOVERY ASSEMBLIES SIMILAR TO QUICKFIT	143	1731
SOLVENT RECOVERY ASSEMBLIES SIMILAR TO QUICKFIT	143	1736
CHROMATOGRAPHY APPARATUS, SMALL SIZE SIMILAR TO QUICKFIT	143	1746
CHROMATOGRAPHY APPARATUS, MEDIUM SIZE SIMILAR TO QUICKFIT	144	1751
ANDREASEN APPARATUS	144	1752
FILTER HOLDER	145 -146	
MILIPORE TYPE ALL GLASS FILTER HOLDER ASSEMBLY WITH I/C JOINT	145	1891
ALL GLASS FILTER HOLDER ASSEMBLY WITH SILICON FITTING	146	1892
KIPPS APPARATUS COMPLETE	147	1896
DESSICATORS	147 - 148	
GLASS DESSICATORS PLAIN WITH LID	147	1906
VACCUM DESSICATORS WITH LID	147	1911
DESSICATORS PLAIN PP	148	1913
DESSICATORS VACCUM PP	148	1915





1. **Price** : Prices are nett Ex-Works Agra. Subject to market fluctuation and alteration without notice.
2. **Packing** : Packing, forwarding, freight, and any other expenses at cost will be charge extra. The company will not be responsible for any breakage or loss in transit.
3. **GST** : GST or IGST or any other tax will be charge extra as per Govt. Rules at the time of supply.
4. **Delivery** : Delivery from ready stock or within 8-10 days.
5. **Payment** : Full Payment Advance against Proforma Invoice.

SUBJECT TO AGRA JURISDICTION ONLY.

Request

While placing orders please quote our
Price-List serial number & please place your order
quantity should be multiple with packing quantity only.

